

AIDA compact NEO  
HL7 Interface Description



PRODUCT INFO  
OR1™

**KARL STORZ  
GmbH & Co. KG**

Mittelstraße 8  
78532 Tuttlingen/Germany  
Postfach 230  
78503 Tuttlingen/Germany

Phone: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105

E-Mail: [info@karlstorz.de](mailto:info@karlstorz.de)  
Web: [www.karlstorz.com](http://www.karlstorz.com)

### ***Change History***

| <b>Version</b> | <b>Date</b> | <b>Changes</b> | <b>Reason</b> | <b>Editor</b> |
|----------------|-------------|----------------|---------------|---------------|
|                |             |                |               |               |
|                |             |                |               |               |
|                |             |                |               |               |
| BA             | 2010-03-03  | Whole document | New document  | UNG           |

## Table of Content

|   |          |
|---|----------|
| <b>Change History .....</b>                     | <b>2</b> |
| <b>Table of Content .....</b>                   | <b>3</b> |
| <b>1 Introduction.....</b>                      | <b>4</b> |
| 1.1 Purpose of the Document.....                | 4        |
| 1.2 Abbreviations .....                         | 4        |
| 1.3 Definitions .....                           | 5        |
| 1.4 References.....                             | 5        |
| <b>2 HL7 Interface of AIDA compact NEO.....</b> | <b>6</b> |
| 2.1 Introduction .....                          | 6        |
| 2.2 Query Messages .....                        | 6        |
| 2.2.1 Patient Query .....                       | 6        |
| 2.2.2 Account Query .....                       | 6        |
| 2.3 Evaluated Message Segments.....             | 6        |
| 2.3.1 PID Segment.....                          | 6        |
| 2.3.2 PV1 Segment.....                          | 7        |
| 2.4 Outgoing Messages .....                     | 7        |
| 2.4.1 ORU Messages.....                         | 7        |
| 2.4.2 MDM Messages .....                        | 7        |
| 2.5 Created Message Segments.....               | 8        |
| 2.5.1 PID + PV1 Segments .....                  | 8        |
| 2.5.2 ORC Segment.....                          | 8        |
| 2.5.3 OBR Segment.....                          | 8        |
| 2.5.4 OBX Segment .....                         | 8        |
| 2.5.5 TXA Segment.....                          | 8        |

## 1 Introduction

### 1.1 Purpose of the Document

This document gives a short description of the HL7 Interface of AIDA compact NEO. The intended audiences are interested customers and IT departments.

### 1.2 Abbreviations

#### General Abbreviations

|          |  |
|----------|--|
| AIDA     | Advanced Image and Data Acquisition / Archiving System = KARL STORZ AIDA™<br>(KST applications for endoscopic image and video acquisition with various functional options) |
| AVA      | AIDA Voice Activation  |
| CCU      | Camera Control Unit  |
| DICOM    | Digital Imaging and Communication in Medicine  |
| HIMSS    | Healthcare Information and Management Systems Society of North America   |
| HIPAA    | Health Insurance Portability and Accountability Act of the United States of America  |
| HIS      | Hospital Information System  |
| HL7      | Health Level Seven - Communication Standard in Medicine  |
| IHE      | Integrating the Healthcare Enterprise  |
| KIS      | Krankenhaus Informations-System  |
| KST      | KARL STORZ GmbH & Co. KG   |
| NEMA     | National Electrical Manufacturers Association of the United States of America  |
| OR       | Operating Room   |
| OR1      | Operating Room all in One = KARL STORZ OR1™<br>(KST solution for operating room integration of OR table, equipment, light and documentation)                               |
| PACS     | Picture Archiving and Communication System   |
| RSNA     | Radiological Society of North America  |
| R-UI     | Realistic User Interface (KST application to control devices over the SCB)   |
| SAM      | Storz Application Manager (KST application to switch between applications)   |
| SCB      | Storz Communication Bus  |
| SEPS     | Storz Endoskop Produktions GmbH Schaffhausen   |
| VA / DoD | Veteran Administration / Department of Defense of the United States of America   |

#### Video specific Abbreviations

|      |  |
|------|--|
| AV   | Audio / Video  |
| HD   | High Definition Television   |
| NTSC | National Television Standard Committee, standard for TV broadcast in USA and Japan |
| PAL  | Phase-Alternating Line, standard for TV broadcast in Europe                        |
| SD   | Standard Definition Television (PAL or NTSC)                                       |
| TV   | Television   |

#### Computer specific Abbreviations

|              |   |
|--------------|---|
| BIOS         | Basic Input / Output System             |
| DLL          | Dynamic Link Library                    |
| GUI          | Graphic User Interface                  |
| KVM-switch   | Keyboard, Video, Mouse - switch         |
| KVM-S-switch | Keyboard, Video, Mouse and RS232 switch |
| NDS          | National Display Systems                |

#### Network specific Abbreviations

|        |   |
|--------|---|
| IP     | Internet Protocol                                 |
| LLP    | Lower Layer Transport Protocol                    |
| TCP    | Transmission Control Protocol                     |
| TCP/IP | Transmission Control Protocol / Internet Protocol |
| TLS    | Transport Layer Security                          |

#### DICOM specific Abbreviations

|     |   |
|-----|---|
| SCP | DICOM Service Class Provider (= Server) |
| SCU | DICOM Service Class User (= Client)     |

## HL7 specific Abbreviations

### Messages

|         |  |
|---------|--|
| ACK     | Acknowledgement  |
| ADT     | Admission, Discharge, Transfer (of patients)                       |
| MDM     | Medical Document Management ()                                     |
| ORM     | General Order Message (Request)                                    |
| ORU     | Unsolicited Transmission of Observation (Result)                   |
| QRY     | Query  |
| QRY-DEM | Demographic Query (patient is identified by Patient ID)            |
| QRY-APA | Account Number Query (patient is identified by Visit ID / Case ID) |

### Message Segments

|       |   |
|-------|---|
| AL1   | Allergy Information Segment               |
| MSH   | Message Header Segment                    |
| ORC   | Common Order Segment                      |
| OBR   | Observation Request Segment               |
| OBX   | Observation/Results Segment               |
| PID   | Patient Identification Segment            |
| PV1   | Patient Visit Segment                     |
| QRD   | Query Definition Segment (original style) |
| QRF   | Query Filter Segment (original style)     |
| TXA   | Transcription Document Header Segment     |
| Z.... | Private Segments                          |

## 1.3 Definitions

### General Definitions

|                          |   |
|--------------------------|---|
| DICOM                    | Digital <u>I</u> maging and <u>C</u> ommunication in <u>M</u> edicine is a standard of NEMA defining how imaging devices in medicine communicate.   |
| HL7                      | <u>H</u> ealth <u>L</u> evel <u>7</u> is a set of standards developed by the Health Level Seven organization for the management, exchange and integration of electronic healthcare information.   |
| IHE                      | <u>I</u> ntegrating the <u>H</u> ealthcare <u>E</u> nterprise is an initiative by HIMSS and RSNA to improve the way computer systems in healthcare share information based on DICOM and HL7.  |
| IHE Integration Profiles | IHE Integration Profiles define the workflow processes and data contents which must be supported by IHE compliant applications.   |
| Modality                 | Refers in medical imaging to an equipment used to acquire images of the body, such as radiography, ultrasound, magnetic resonance imaging or endoscopy.   |
| Procedure                | A medical procedure is a course of action intended to achieve a result in the care of a patient, normally a diagnosis or a therapy. Procedures are often complex and include a number of different steps over an extended time period.                          |
| Procedure Step           | A procedure step is a part of a procedure performed at a given time and location. In a diagnostic procedure a procedure step is e.g. a CT scan or a MR scan. In surgery, a procedure step includes a single surgical event, e.g. a Laparoscopy or a Cystoscopy. |

### HL7 specific Definitions

|              |  |
|--------------|--|
| HL7 Listener | A permanently running application listening for HL7 messages.                |
| HL7 Server   | A HIS application capable of sending and receiving HL7 messages.             |
| QRY Server   | An HL7 server in the HIS capable of handling and responding to QRY messages. |

## 1.4 References

|           |   |   |        |
|-----------|---|---|--------|
| /REF_001/ | DICOM Standard 2007                             | <a href="http://medical.nema.org/">http://medical.nema.org/</a>   | NEMA   |
| /REF_002/ | IHE Technical Framework                         | <a href="http://www.ihe.net/">http://www.ihe.net/</a>             | IHE    |
| /REF_003/ | VA / DoD DICOM Conf. Req. for Digital Acq. Mod. | <a href="http://www.va.gov/imaging">http://www.va.gov/imaging</a> | VA/DoD |
| /REF_004/ | HL7 Standard Version 2.3.1, 2.4 and 2.5.1       | <a href="http://hl7.org/">http://hl7.org/</a>                     | HL7    |

## 2 HL7 Interface of AIDA compact NEO

### 2.1 Introduction

AIDA compact NEO is intended to capture image data (images, videos and audio clips) in the operating room. It can communicate through Ethernet with a HL7 or DICOM based HIS.

Through its HL7 interface, AIDA compact NEO can query for patient data when a new procedure is started and send result messages when a procedure is finished. These functions can be activated separately, but require purchase of the corresponding license option.

AIDA compact NEO is compatible with HL7 versions 2.3 to 2.5. The HL7 communication uses TCP/IP. The HL7 messages are framed by Hex 0B (start) and Hex 1C 0D (stop).

### 2.2 Query Messages

AIDA compact NEO can query for patient data any time during a procedure. It is recommended to do this when a new procedure is started. This requires that there is an HL7 server present which responds to the query. To get a meaningful query, the user must first enter or read-in the Patient ID or Visit ID.

#### 2.2.1 Patient Query

The HL7 Interface of AIDA compact NEO sends the following patient query:

- Query message type = QRY^A19
- What Subject Filter = DEM

In the query responses selected fields from the PID and PV1 segments are evaluated.

#### Example of sent demographic query

```
MSH|^~\&|AIDA|DEMOKIS||201002161200||QRY^A19|Q140|P|2.3.1||
QRD|201002161200|R|I|Q9||1^RD|PID-00-004|DEM|
```

#### Example of processed response

```
MSH|^~\&|DEMOKIS|AIDA|201002161200|ADR^A19|P|2.3.1|||||
MSA|AA|Q140|||||
QRD|201002161200|R|I|Q140||1^RD|PID-00-004|DEM|
PID|1651||PID-00-004|Patient^Four^R.^CBE^Sir^Ph.D.^L||19300711|M|||||
PV1|I|||||10-JONES^Jones^Miriam^^Dr.|||||||V-ID-004|||||
```

#### 2.2.2 Account Query

The HL7 Interface of AIDA compact NEO sends the following account query:

- Query message type = QRY^A19 or QRY^Q01
- What Subject Filter = APA

#### Example of sent account number query

```
MSH|^~\&|AIDA|DEMOKIS||201002161200||QRY^Q01|Q140|P|2.3.1||
QRD|201002161200|R|I|Q9||1^RD|V-ID-004|APA|
```

#### Example of processed response

```
MSH|^~\&|DEMOKIS|AIDA|201002161200|ADR^Q01|P|2.3.1|||||
MSA|AA|Q140|||||
QRD|201002161200|R|I|Q140||1^RD|V-ID-004|APA|
PID|1651||PID-00-004|Patient^Four^R.^CBE^Sir^Ph.D.^L||19300711|M|||||
PV1|I|||||10-JONES^Jones^Miriam^^Dr.|||||||V-ID-004|||||
```

## 2.3 Evaluated Message Segments

### 2.3.1 PID Segment

| PID Segment in QRY Message |       |                               |       | HL7 Version    |     |     | AIDA compact NEO |                          |
|----------------------------|-------|-------------------------------|-------|----------------|-----|-----|------------------|--------------------------|
| Seg.                       | Field | Element Name                  | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute      |
| PID                        | 2     | Patient ID <sup>1)</sup>      | 1     |                | O   | B   | B                | Patient ID <sup>2)</sup> |
| PID                        | 3     | Patient ID List <sup>1)</sup> | 1     |                | R   | R   | R                |                          |

| PID Segment in QRY Message |       |                          |       | HL7 Version    |     |     | AIDA compact NEO |                      |
|----------------------------|-------|--------------------------|-------|----------------|-----|-----|------------------|----------------------|
| Seg.                       | Field | Element Name             | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute  |
| PID                        | 4     | Alternate Patient ID     | 1     |                | O   | B   | B                |                      |
| PID                        | 5     | Patient Name             | 1     | Family Name    | R   | R   | R                | Patient's Name       |
|                            |       |                          | 2     | First Name     | O   | O   | O                |                      |
|                            |       |                          | 3     | Middle Name    | O   | O   | O                |                      |
|                            |       |                          | 5     | Prefix         | O   | O   | O                |                      |
| PID                        | 7     | Date/Time of Birth       | 1     |                | O   | O   | O                | Patient's Birth Date |
| PID                        | 8     | Admin. Sex <sup>3)</sup> | 1     |                | O   | O   | O                | Patient's Sex        |

<sup>1)</sup> In HL7 V2.3, PID-2 is called "Patient ID (external ID)" and PID-3 "Patient ID (internal ID)".

<sup>2)</sup> The HL7 field used is configurable.

<sup>3)</sup> Defined values F = Female, M = Male, O = Other, U = Unknown, etc.

### 2.3.2 PV1 Segment

| PV1 Segment in QRY Message |       |                  |       | HL7 Version    |     |     | AIDA compact NEO |                            |
|----------------------------|-------|------------------|-------|----------------|-----|-----|------------------|----------------------------|
| Seg.                       | Field | Element Name     | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute        |
| PV1                        | 8     | Referring Doctor | 1     | ID Number      | O   | O   | O                | Referring Physician's Name |
|                            |       |                  | 2     | Family Name    | O   | O   | O                |                            |
|                            |       |                  | 3     | First Name     | O   | O   | O                |                            |
|                            |       |                  | 4     | Middle Name    | O   | O   | O                |                            |
|                            |       |                  | 7     | Degree         | O   | O   | O                |                            |
| PV1                        | 19    | Visit Number     | 1     |                | O   | O   | O                | Admission ID / Case ID     |

## 2.4 Outgoing Messages

AIDA compact NEO can send result messages. This requires that there is an HL7 server present which receives these messages.

### 2.4.1 ORU Messages

The HL7 Interface of AIDA compact NEO creates the following ORU messages:

- ORU message type = R01
- ORU message structure =
  - MSH Message Header segment
  - PID Patient Identification segment
  - PV1 Patient Visit segment
  - ORC Common Order segment
  - OBR Observation Request segment
  - {OBX} Observation Results segment (one OBX segment per image or video).

### Example of created ORU message

```
MSH|^~\&|AIDA|DEMOKIS||201002161400||ORU^R01|10|P|2.3.1||AL|NE|||
PID|||PID-00-004||Patient^Four^R.^^^|19300711|M|||||||||||||||||
PV1|I|||||10-JONES^Jones^Miriam^^^Dr.|||||||V-ID-004|
ORC|RE||||CM|||||||||||||
OBR|1||Y02.56^HYSTEROSCOPY||201002161400|||||||||||||F|||||
ONEIL&O'Neil&Aron&&&&Prof.Dr.||||
OBX|1|RP|ENDOSCOPIC IMAGE^AIDA^L|1|\KST-IHE\HL7Export\....\Images\Image_1.BMP
|||||F|||201002161400|
```

### 2.4.2 MDM Messages

The HL7 Interface of AIDA compact NEO creates the following MDM messages:

- MDM message type = T01
- MDM message structure =
  - MSH Message Header segment
  - PID Patient Identification segment
  - PV1 Patient Visit segment
  - ORC Common Order segment
  - OBR Observation Request segment
  - TXA Transcription Document Header Segment

**Example of created MDM message**

```
MSH|^~\&|AIDA||DEMOKIS||201002161400||MDM^T01|10|P|2.3.1|||AL|NE|||
PID||PID-00-004||Patient^Four^R.^^^|19300711|M|||
PV1|I|||||10-JONES^Jones^Miriam^^^Dr.|||||V-ID-004|
ORC|RE||||CM|||||
OBR|1|||Y02.56^HYSTEROSCOPY|||201002161400|||||F|||||
12-ONEIL&O'Neil&Aron&&&Prof.Dr.|||||
TXA|1|OP|IM|||201002161400|||12-ONEIL^O'Neil^Aron^^^Prof.Dr.||||
^1.2.276.0.67.2.121082829.20091122151133.3^L||||
\\KST-IHE\HL7Export\.....\Images\Image_1.BMP|DO|AV|
```

**2.5 Created Message Segments**

**2.5.1 PID + PV1 Segments**

See above.

**2.5.2 ORC Segment**

| ORC Segment in ORU and MDM Messages |       |               |       | HL7 Version    |     |     | AIDA compact NEO |  |
|-------------------------------------|-------|---------------|-------|----------------|-----|-----|------------------|--|
| Seg.                                | Field | Element Name  | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute                    |
| ORC                                 | 1     | Order Control | 1     |                | R   | R   | R                | Constant value <b>RE</b> <sup>1)</sup> |
| ORC                                 | 5     | Order Status  | 1     |                | O   | O   | O                | Constant value <b>CM</b> <sup>1)</sup> |

<sup>1)</sup> Meaning of used values RE = observation results included, CM = order is completed.

**2.5.3 OBR Segment**

| OBR Segment in ORU and MDM Messages |       |                      |       | HL7 Version    |     |     | AIDA compact NEO |  |
|-------------------------------------|-------|----------------------|-------|----------------|-----|-----|------------------|--|
| Seg.                                | Field | Element Name         | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute                                    |
| OBR                                 | 4     | Universal Service ID | 1     | Identifier     | R   | R   | R                | Procedure Name Part 1<br>Part 2<br>Part 3              |
|                                     |       |                      | 2     | Text           | R   | R   | R                |  |
|                                     |       |                      | 3     | Coding System  | R   | R   | R                |  |
| OBR                                 | 7     | Obs. Date/Time       | 1     |                | C   | C   | C                | Treatment Date   |
| OBR                                 | 25    | Result Status        | 1     |                | C   | C   | C                | Constant Value <b>F</b>                                |
| OBR                                 | 32    | Principal Result     | 1.1   | ID Number      | O   | O   | O                | Performing Physician's<br>Name (Surgeon) <sup>1)</sup> |
|                                     |       |                      | 1.2   | Family Name    | O   | O   | O                |  |
|                                     |       |                      | 1.3   | First Name     | O   | O   | O                |  |
|                                     |       |                      | 1.4   | Middle Name    | O   | O   | O                |  |
|                                     |       |                      | 1.5   | Suffix         | O   | O   | O                |  |
|                                     |       |                      | 1.6   | Prefix         | O   | O   | O                |  |

<sup>1)</sup> Value sent as entered, no re-formatting.

**2.5.4 OBX Segment**

| OBX Segment in ORU and MDM Messages |       |                    |       | HL7 Version    |     |     | AIDA compact NEO |  |
|-------------------------------------|-------|--------------------|-------|----------------|-----|-----|------------------|--|
| Seg.                                | Field | Element Name       | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute                    |
| OBX                                 | 2     | Value Type         | 1     |                | C   | C   | C                | Constant value <b>RP</b> <sup>1)</sup> |
| OBX                                 | 4     | Obs. Sub ID        | 1     |                | C   | C   | C                | Counter                                |
| OBX                                 | 5     | Observation Value  | 1     |                | C   | C   | C                | Reference pointer <sup>2)</sup>        |
| OBX                                 | 11    | Obs. Result Status | 1     |                | R   | R   | R                | Constant value <b>F</b> <sup>1)</sup>  |
| OBX                                 | 14    | Obs. Date/Time     | 1     |                | O   | O   | O                | Treatment Date                         |

<sup>1)</sup> Meaning of used RP = reference pointer, F = final result.

<sup>2)</sup> The observation value is a reference to an image or video file on a file server.

**2.5.5 TXA Segment**

| TXA Segment in MDM Message |       |                    |       | HL7 Version    |     |     | AIDA compact NEO |  |
|----------------------------|-------|--------------------|-------|----------------|-----|-----|------------------|--|
| Seg.                       | Field | Element Name       | Subs. | Component Name | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute                      |
| TXA                        | 1     | Set ID             | 1     |                | R   | R   | R                | Counter                                  |
| TXA                        | 2     | Document Type      | 1     |                | R   | R   | R                | Constant value <b>OP</b> <sup>1)</sup>   |
| TXA                        | 3     | Doc. Content Pres. | 1     |                | C   | C   | C                | Constant value <b>IM</b> <sup>1)</sup>   |
| TXA                        | 6     | Orig. Date/Time    | 1     |                | O   | O   | O                | Acquisition Date/Time                    |
| TXA                        | 9     | Orig. Code/Name    | 1     | ID Number      | O   | O   | O                | Performing Physician's<br>Name (Surgeon) |
|                            |       |                    | 2     | Family Name    | O   | O   | O                |  |



| TXA Segment in MDM Message |       |                    |       | HL7 Version       |     |     | AIDA compact NEO |  |
|----------------------------|-------|--------------------|-------|-------------------|-----|-----|------------------|--|
| Seg.                       | Field | Element Name       | Subs. | Component Name    | 2.3 | 2.4 | 2.5              | Rel. AIDA Attribute                    |
|                            |       |                    | 3     | First Name        | O   | O   | O                |  |
|                            |       |                    | 4     | Middle Name       | O   | O   | O                |  |
|                            |       |                    | 5     | Suffix            | O   | O   | O                |  |
|                            |       |                    | 6     | Prefix            | O   | O   | O                |  |
| TXA                        | 12    | Unique Doc. No.    | 1     | Entity Identifier | O   | O   | O                | Empty                                  |
|                            |       |                    | 2     | Namespace ID      | O   | O   | O                | Empty                                  |
|                            |       |                    | 3     | Universal ID      | R   | R   | R                | Instance UID                           |
|                            |       |                    | 4     | Universal ID Type | O   | O   | O                | Constant value <b>L</b>                |
| TXA                        | 16    | Doc. File Name     | 1     |                   | O   | O   | O                | File Name <sup>2)</sup>                |
| TXA                        | 17    | Doc. Compl. Status | 1     |                   | R   | R   | R                | Constant value <b>DO</b> <sup>1)</sup> |
| TXA                        | 19    | Doc. Availability  | 1     |                   | O   | O   | O                | Constant value <b>AV</b> <sup>1)</sup> |

<sup>1)</sup> Meaning of used values OP = operative report, IM = image data, AU = audio data, TEXT = text, DO = documented, AV = available.

<sup>2)</sup> AIDA compact NEO provides the path and file name on the file server.