STORZ

KARL STORZ—ENDOSKOPE

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**Cleaning, Disinfection, Care and Sterilization of
KARL STORZ Instruments**

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1 Foreword

Hygiene must be treated as an integral responsibility and a fundamental component of the whole process of medical treatment.

This principle begins with the design of instruments: alongside the technical solution, effective cleaning, disinfection and sterilization represent the central factors in the development of all instruments.

The objective is to create an environmentally-friendly and long-lasting reusable instrument, the materials, design and quality of which allow preparation using a suitable process that takes all applicable norms and laws into consideration, while maintaining the full functionality of the instrument.

Optimal handling throughout the instrument cycle goes further still and complements the instrument with individually coordinated accessories for optimized storage during the entire process – from the operating theater to transportation to the preparation department for medical devices, right down to the individual steps in the process.

Intelligent storage systems protect the medical devices and reduce the effort expended in handling. A reduced amount of repacking guarantees maximum protection for the instruments, optimizes operating theater procedures and ensures the desired results of the preparation. KARL STORZ achieves this objective by continuously developing new instruments and storage systems and improving existing ones according to the demands of medical professionals.

As the supplier of these medical devices, KARL STORZ shall only be liable for the safety, reliability and performance of the products if installation, maintenance and repairs have been carried out exclusively by personnel authorized by KARL STORZ and the product has been used in accordance with the instruction manual.

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7 Maintenance	74
7.1 Troubleshooting/fault finding	74
7.2 Warranty	75
7.3 Important information.....	75
7.4 Repairs	75
7.5 Repair and replacement program	76
8 Appendix	77
8.1 Extract from the catalog "HYGIENE, Care, Sterilization and Storage Techniques"	77

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2 Terms used

AAMI-TIR	Association for the Advancement of Medical Instrumentation – Technical Information Report
ANSI	American National Standards Institute
Arbeitsanweisung	
AMP	Preparation of medical devices
BfArM	German Federal Institute for Drugs and Medical Devices
BGV A1	German Guidelines for Health and Safety in the Workplace from the Employers' Liability Insurance Association
BioStoffV	German Ordinance on Biological Substances
CEN	European Committee for Standardization
Charge	Charging unit for the sterilizer or washer and disinfecter
Coagulation	
Completely demineralized water	
CSSD	Central Sterile Supply Department
Device manufacturer	
Disinfection	
DGHM	German Society for Hygiene and Microbiology
DGSV	German Society for Sterile Supply

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MD	Medical device
MDU	Medical device unit (measure of performance)
MPG	German Medical Devices Act
MPBetreibV	German Medical Devices Operator Ordinance
RKI	Robert Koch Institute
MPV	German Medical Devices Ordinance
Oil pen	
OR	Operating room/theatre
Passivation	
PAA	Per(oxy)acetic acid
PI	Procedure instruction
PPE	Personal protective equipment
Process chemical	
Pyrogen-free	
QUAT	Quaternary ammonium compounds
QCK	Quick Connect Kit for the Steris® System 1/1E
QM(S)	Quality management system in accordance with DIN EN ISO 9001 or 13485, for example
RDG	Washer and disinfectant
RDG-E	Washer and disinfectant for endoscopes
Recommendations	Example: RKI (Robert Koch Institute)/BfArM (German Federal Institute for Drugs and Medical Devices) recommendation regarding "Hygiene Requirements for Reprocessing Medical Devices"

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RO	Reverse osmosis
RT	Room temperature
SDS	Safety Data Sheet
SOP	Standard Operating Procedure
Stainless steel	
Sterile barrier system	
StU	Sterilization unit (volume measurement of 600 x 300 x 300 mm)
TRBA	German Technical Rules for Biological Agents
TRGS	German Technical Rules for Hazardous Substances
US	Ultrasound
UVV	German Accident prevention regulation (old name – see BGV A1)
VAH	German Association for Applied Hygiene
Validation	
Wallhäußers „Praxis der Sterilisation, Desinfektion, Antiseptik und Konservierung“ Thieme-Verlag 2008	Wallhäußers “Practice of Sterilization, Disinfection, Antisepsis and Conservation” published by Thieme-Verlag, 2008
Working instructions	

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3 Introduction

Thank you for your expression of confidence in the KARL STORZ brand name. Our medical devices are the result of years of experience and great care in manufacture. With this purchase, you have decided in favor of modern, high-quality instruments.

This manual was created with reference to and under consideration of DIN EN ISO 17664 and is intended to provide guidelines and information for personnel responsible for the preparation of instruments.

KARL STORZ instruments are developed for and with operating surgeons. Therefore, function and design are geared to meet the requirements of the surgeon. Hygiene requirements are given substantial consideration in the design of the instruments. At the same time, finding a balance between function and hygiene represents a compromise in part which must nonetheless ensure that the instrument meets all requirements, with the practicality of the “tool” on the one hand, and efficient preparation (in the context of repair and preservation) for re-use on the other.

The medical devices manufactured by KARL STORZ are, as a rule, sensitive precision instruments which require particular care during preparation. To enable proper, expert preparation, this manual renders assistance with regard to the disassembly, cleaning, disinfection, care, assembly and sterilization of KARL STORZ instruments.

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You will also find the current list of approved chemicals on the Internet at www.karlstorz.com. In addition, an extract from the catalog can be found in the appendix (see catalog “HYGIENE, Care, Sterilization and Storage Techniques”).

WARNING: *Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination indicates that preparation has not been performed, or has been performed incorrectly. Prepare the instruments before initial use and before and after every subsequent use using validated procedures.*

NOTE: *The manual “Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments” can be downloaded or requested at www.karlstorz.com.*

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology may occur for this reason. Therefore, no claims may be derived from the information, illustrations, and descriptions in this manual.

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4 About this manual

4.1 Using this manual

This manual is not an instruction manual but rather provides advice regarding the cleaning, disinfection, care, sterilization and storage/transportation of KARL STORZ instruments.

Instruments of a previous generation, which no longer comply with the current, valid regulatory requirements in respect of their effective hygienic preparation, are not covered by the procedures and processes described in this manual.

To provide an explanation of disassembly/assembly as well as information on the preparation methods to be used in each case, specific instruments have an accompanying instruction leaflet. These are included in the scope of supply of the instrument, but can also be accessed via the Hygienics specialty area at our homepage www.karlstorz.com if necessary.

Instruments which are not supplied with a specific instruction manual can, in principle, be prepared using the following procedures:

Cleaning:	manually at room temperature, by machine at up to 55°C for 10 mins.
Disinfection:	manually at room temperature, by machine thermally at 93°C for 5 mins.
Steam sterilization:	fractionated prevacuum procedure at 134-137°C for 4-5 mins.

i NOTE: Machine cleaning/disinfection should be preferred to manual processes, since machine processes can be standardized, reproduced and therefore validated.

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4.2 Symbols and hazard information

Please read and observe this manual.

The words Warning, Caution, and Note convey special meanings. Wherever they are used in this manual, the following text should be carefully reviewed to ensure the safe and effective use of the instrument. To make the words Warning, Caution and Note stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates a potential hazard. Failure to observe a Warning may result in injury to the patient, user or a third party.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the medical devices.

NOTE: Notes contain special information regarding use or important explanations.

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4.3 Overview of important points

- Observe all Warnings and Notes for your own protection as well as for the protection of the personnel around you.
- Failure to observe this manual poses the risk of damage to the instrument, increased instrument wear and invalidation of the warranty.
- Before brand-new and repaired instruments are sterilized for the first time, they must be subjected to at least one complete preparation cycle.
- Instruments which have been used on patients must be taken for cleaning/disinfection without delay in order to avoid the unnecessary drying of contaminants on the instrument, as this may have a strong negative impact on cleaning in areas which are difficult to access.

i NOTE: Process validations for cleaning/disinfection are carried out by KARL STORZ with a drying time of 1 hour at room temperature.

- Disassemble instruments before cleaning/disinfection as specified in the manual and open joints to approx. 90°.
- To aid the preliminary cleaning of heavily soiled and/or complex instruments, an ultrasound treatment at approx. 35 kHz with a maximal exposure time of 10 minutes can/must be performed. As a final measure, the instruments must be rinsed with cold water to ensure neutralization.
- It is imperative that the specifications of the chemical manufacturer regarding concentration, temperature and exposure time for the chemicals used in cleaning/disinfection are observed.
- Only use validated procedures/processes for machine cleaning/disinfection.
- The operator is responsible for the on-site validation of the preparation processes.
- Manual processes must be described in detailed working instructions and also require validation of their efficacy.

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- Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.
- Only use cleaning accessories which are recommended and suitable (e.g., brushes).

CAUTION: Never use metal brushes or metal sponges for cleaning or repeated cleaning.

NOTE: Aerosols should be removed by means of a suction device.

NOTE: Only use recognized disinfectants for chemical disinfection (In Germany: VAH/RKI-listed).

- To avoid acoustic shadows and unrinsed areas, do not overload washer and disinfectant or ultrasound units.
- KARL STORZ recommends using completely demineralized water for the final rinse following cleaning and chemical disinfection. The instruments must be residue-free upon submission for packaging and sterilization.
- Remove instruments immediately after machine cleaning and disinfection in the washer and disinfectant (danger of corrosion).
- Drying: The level of dryness of the instruments to be sterilized can have a significant influence on the efficacy of the sterilization (depending on the method). It must be ensured that the instruments are absolutely dry when submitted for sterilization. This must be carried out whilst allowing for the various influential factors such as the load, the geometry of the instruments and management of the process on-site.
- Prior to the functional check and specific care measures, each instrument must be inspected for purity/cleanliness and, if necessary, subjected to repeated cleaning.

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- Instruments must be cared for locally (on moving parts, e.g., joints, threads) after each cleaning/disinfection procedure. KARL STORZ recommends the instrument oil with Art. No. 27656 for general care and the oil pen with Art. No. 27656 S for easy-to-reach, movable joints.
- Each instrument must undergo a functional check before packaging and sterilization. Disassembled instruments must be assembled for checking and remain in this state for packaging and sterilization. Exceptions will be specifically indicated.
- Only clean, disinfected and dry instruments which have been cared for and are in full working order may be packed and sterilized.
- In order to check specific instruments, KARL STORZ recommends the use of a suitable inspection device, e.g., magnifying glass/ loupe.
- Defective instruments must be removed from circulation and sent to the manufacturer for repair.

i NOTE: Repairs carried out by third parties may lead to invalidation of the warranty.

- Before returning to the manufacturer, any instruments which require a repair must – as far as practicable – be subjected to a complete preparation cycle. The corresponding documentary proof must be enclosed with the delivery.

i NOTE: Disinfection is not assured without cleaning first! Sterilization is not assured without disinfecting first!

- For sterilization, only validated procedures which are sensible and suitable for the medical device concerned may be used.

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5 General information on preparation

WARNING: The stated cleaning, disinfection and sterilization methods alone do not guarantee the necessary level of hygienic safety in the instruments. This can only be achieved if recognized and validated preparation methods are followed.

NOTE: Due to differences in their designs, properties and applications, the instruments must be assessed prior to initial preparation, taking the manufacturer's specifications given in the instruction manual into consideration.

Bear the following in mind:

- Thermal sensitivity: thermostable or thermolabile
- Composition of the instruments: smooth surfaces – with crevices, folding or purgable – instruments with joints, non-separable – separable instruments with/without hollow spaces – rigid/flexible optical instruments
- Compatibility of process chemicals
- Efficiency test (process validation)

Recommended accessories for cleaning and care as well as tray holders, trays and sterilization containers can be found in the catalog “HYGIENE, Care, Sterilization and Storage Techniques”.

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5.1 Definitions

Instruments and accessories for endoscopic surgery are, as a rule, sensitive precision instruments which require particular care during preparation. Only trained personnel (according to §4 of the Medical Devices Operator Ordinance in Germany) are in a position to carry out disassembly, cleaning, disinfection, inspection/care, assembly and sterilization procedures in such a way that the perfect function of the instruments is maintained.

Some instruments, such as those used in minimally invasive surgery, for example, which are exposed to gas overpressure during operations, become contaminated on surfaces which cannot be seen and are hard to reach. Therefore, where technically feasible or compatible with the application, these instruments have an irrigation channel or are developed for easy disassembly.

Bear the following in mind:

- Thermal sensitivity: thermostable or thermolabile
- Composition of the instruments: smooth surfaces – with crevices, folding or purgable – instruments with joints, non-separable – separable instruments with/without hollow spaces – rigid/flexible optical instruments
- Compatibility of process chemicals
- Efficiency test (process validation)

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Recommended accessories for cleaning and care as well as tray holders, trays and sterilization containers can be found in the catalog “HYGIENE, Care, Sterilization and Storage Techniques”.

5.1.1 Cleaning

Cleaning denotes procedures through which undesirable substances and other contaminations are removed, e.g., from surfaces. Blood, mucus, excretion and residues from pharmaceuticals and disinfectants, amongst others, are all considered undesirable here. Cleaning is normally defined as a reduction in the number of microorganisms by a factor of 10^3 , but can also be defined by the depletion in the level of contamination to a specific limit value.

5.1.2 Disinfection

Disinfection denotes procedures through which microorganisms, e.g., bacteria, spores, viruses and fungi are killed, inactivated and/or removed. Disinfection is a reduction by a factor of at least $10^4 - 10^6$.

High Level Disinfection (HLD)

This is understood to be so-called “cold sterilization” (“sporocidal disinfection with specific test spores” according to Wallhäußer's “Practice of Sterilization, Disinfection, Antisepsis and Conservation”, p. 164 et seq., published by Thieme-Verlag, 2008).

It concerns a process which requires a shorter contact time to kill off all microorganisms than sterilization, but does not kill off large quantities of spores. The reduction stipulated for this process amounts to a minimum of 10^6 (see ANSI AAMI ST58).

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- i** **NOTE:** Sterilization methods for the hygienic preparation of medical devices which are needed for urgent surgery are not always available. In such cases, the following alternatives, which must be documented, can be used to destroy all bacteria, including tuberculosis pathogens, fungi and viruses (except for some endospores):
- Immersion of the medical device in boiling water for a minimum of 20 mins.
 - Application of pressure using saturated steam (without overpressure)
 - Chemical disinfection for an extended period of time

Intermediate Level Disinfection (ILD)

This process ensures the elimination of virus material, mycobacteria, fungi and vegetative bacteria, but not of permanent forms (spores).

Low Level Disinfection (LLD)

This process ensures the elimination of viruses with lipid membranes, some forms of fungi and vegetative bacteria, but not of permanent forms (spores).

5.1.3 Sterilization

Sterilization is characterized as procedures which inactivate viable microorganisms, including spores, and make the transmission of infections impossible.

An object is thus considered sterile if the theoretical value of no more than one living microorganism in 1×10^6 sterilized units (StE) is reached in the end product (European Pharmacopoeia Ph.Eur. 6.0).

This means that from an original figure of 1,000,000 fertile microorganisms, no more than one survives: “Sterility Assurance Level” (SAL). A product may only be described as “sterile” if it was sterilized according to validated procedures in packaging which is resistant to microorganisms (sterile barrier system).

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5.1.4 Validation

Validation is the documentary evidence that a process or system reproducibly fulfills the previously specified requirements (acceptance criteria).

Within the framework of the conformity assessment, the manufacturer of the medical device is responsible for providing the evidence of the efficacy of the cleaning, disinfection and sterilization procedures (validation of the procedures).

The operator is responsible for proving the efficacy of (validating) the procedures on-site with due consideration given to the specific framework conditions and ensuring that the routine tests are carried out periodically.

5.2 Preparation procedures

The regular procedures for cleaning, disinfection and sterilization used in clinical facilities are named here.

i **NOTE:** Where possible, machine cleaning/disinfection should be carried out, followed by fractionated steam sterilization (according to DIN EN ISO 17665-1) at 134-137 °C with a hold time of at least 4-5 mins. These procedures are only suitable for thermostable instruments.

Preparation procedures have numerous variants, which is why care should be taken to ensure that their efficacy has been tested (and validated).

As a fundamental principle, always ensure that cleaning is carried out first, and disinfection afterwards in a separate procedure! Exceptions are governed by the Infection Protection Act (in Germany) or the national regulations regarding the protection of personnel.

Rationale for the separation of cleaning and disinfection processes: Disinfectants may contain fixing components and therefore impede cleaning. Thorough cleaning reduces the number of microorganisms, facilitates the subsequent disinfection and makes it more effective.

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CAUTION: Risk of damaging the products: The use of chemicals which have not been approved by KARL STORZ may cause damage to the instruments. Only use chemicals approved by KARL STORZ for preparation. You will find a complete list on the Internet at www.karlstorz.com.

5.2.1 Cleaning procedures

When selecting the instrument cleaner, the range of action should be considered, depending on the field of application of the medical device. Highly diverse contaminations, such as blood, mucus and bone constituents, make specific demands of the process chemicals. The chemical manufacturer's specifications (e.g., concentration, exposure time, temperature) for

- manual cleaning
- machine cleaning

must be observed (DIN EN ISO 15883).

NOTE: Machine cleaning is preferred.

5.2.2 Disinfection procedures

When selecting the disinfectant, the range of action should be taken into account with regard to microorganisms and disinfection itself should be carried out with a suitable disinfectant for instruments (In Germany: VAH/RKI-listed) in compliance with the specifications of the process chemical manufacturer (e.g., concentration, exposure time, temperature):

- chemical, carried out manually
- thermal, carried out by machine
- chemico-thermal, carried out by machine

Machine thermal disinfection is preferred (80-95°C or upon reaching the A_0 value of 3000 in accordance with DIN EN ISO 15883, see appendix).

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5.2.3 Sterilization methods (physical or chemical)

- Steam sterilization in the fractionated prevacuum procedure (DIN EN ISO 17665)
- Low-temperature steam and formaldehyde sterilization process (LTSF) (DIN EN 14180; DIN 58948-17)
- Hydrogen peroxide (H₂O₂) sterilization (e.g., ASP STERRAD®, Steris® AMSCO® V-PRO™ 1)
- Chemical low-temperature sterilization with peracetic acid (e.g., Steris® System 1/1E)
- Ethylene oxide sterilization (6-100%) (DIN EN ISO 11135)
- Hot air sterilization (not recommended)

i NOTE: Steam sterilization in the fractionated prevacuum procedure is preferred.

5.3 Water and steam quality

5.3.1 Cleaning/disinfection/rinsing

To prolong the useful service life of the instruments, KARL STORZ recommends using completely demineralized water for these three preparation steps. Changes to the surface may impair the functional and hygienic properties of the instrument.

Cases of discoloration due to oxide formation represent purely cosmetic changes which do not affect functioning or hygienic safety. Areas of discoloration as a result of oxide formation are purely cosmetic changes

and must be differentiated from potentially infectious adhesions to the surface, which must be removed but are partially covered by fibrin residues, which pass for water stains.

Mechanically induced changes to the surface such as scratches must be critically assessed in respect of the risks posed to hygiene and function and may lead to corrosion. Corrosion results in an assessable threat to the function of the instrument, perhaps an unexpected break when exerting a force on it during an operation (due to pitting, crevice corrosion or stress corrosion cracking).

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Only completely demineralized water should be used for the final rinse, in accordance with the “Guideline Compiled by the DGKH (German Society for Hospital Hygiene), DGSV (German Society for Sterile Supply) and AKI (Working Group on Instrument Reprocessing) for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors” (which refers to DIN EN ISO 15883-1 Point 6.4.2).

Conductivity:	≤ 15 µs (deviation from DIN EN 285)
pH value:	5-7
Total hardness:	≤ 0.02 mmol CaO/l
Salt content:	≤ 10 mg/l
Phosphate: (as P2O5)	≤ 0.5 mg/l
Silica: (as SiO2)	≤ 1 mg/l
Chloride:	≤ 2 mg/l

Table 1: In accordance with the Guideline, the following values are recommended

These types of water are free from microorganisms and are suitable for the final rinse:

- Sterile-filtered water
- Sterile water

These types of water are not free from microorganisms and are therefore not suitable for the final rinse:

- Distilled water
- Deionized water
- Water from the faucet
- Aqua purificata

i NOTE: Failure to heed these recommendations may lead to corrosion, discoloration and staining (in the Red Brochure “Proper Maintenance of Instruments” from the A-K-I (Working Group on Instrument Reprocessing) www.a-k-i.org).

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Tab. 2: DIN EN 285:2009-08 Anhang B

5.3.2 Sterilization

For sterilization, water quality must comply with DIN EN 285:2009-08 Appendix B for feed water and condensate. This compliance should be checked regularly according to recognized analytical methods.

Feed water:	Conductivity at 25°C ≤ 5 µs/cm pH value (degree of acidity) 5 to 7.5
Condensate:	Conductivity at 25°C ≤ 3 µs/cm pH value (degree of acidity) 5 to 7

Table 2: DIN EN 285:2009-08 Appendix B

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Description	Drinking water	Highly purified water
Bacteria (cfu/mL)	< 200	≤ 10
Endotoxin (EU/mL)	n/a	< 10
Total organic carbon (mg/L)	< 1.0	< 0.05
pH	6.5-8.5	n/a
Hardness of water (ppm CaCO ₃)	< 150	< 1.0
Resistance (M'Ω-cm)	n/a	> 1.0
Total dissolved solid (mg/L CaCO ₃)	< 500	< 0.4
Ionic contamination	< 250	< 0.2
Chlorine (mg/L)	< 0.3	< 0.2
Iron (mg/L)	< 0.1	< 0.1
Copper (mg/L)	< 0.1	< 0.1
Manganese (mg/L)	< 0.1	< 0.1
Color and cloudiness	colorless, clear, no particles	colorless, clear, no particles
Reference: AAMI TIR34:2007 "Water for the reprocessing of medical devices"		

Table 3: Characteristics typically associated with different grades of water quality (Categories)

Descripción	Agua potable	Agua de alta pureza
Bacterias (cfu/mL)	< 200	≤ 10
Endotoxina (EU/mL)	n/a	< 10
Carbono orgánico, total (mg/L)	< 1,0	< 0,05
pH	6,5 – 8,5	n/a
Dureza del agua (ppm CaCO ₃)	< 150	< 1,0
Resistencia (M'Ω-cm)	n/a	> 1,0
Materia sólida disuelta, total (mg/L CaCO ₃)	< 500	< 0,4
Contaminación iónica	< 250	< 0,2
Cloro (mg/L)	< 0,3	< 0,2
Hierro (mg/L)	< 0,1	< 0,1
Cobre (mg/L)	< 0,1	< 0,1
Manganeso (mg/L)	< 0,1	< 0,1
Color y turbidez	incolora, limpia, sin partículas	incolora, limpia, sin partículas
Referencia: AAMI TIR34:2007 "Water for the reprocessing of medical devices"		

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**5.4 Brand-new and repaired
instruments**

WARNING: Risk of infection: These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination indicates that preparation has not been performed, or has been performed incorrectly. Prepare the instruments before initial use and before and after every subsequent use using validated procedures.

NOTE: New instruments or those with finished surfaces exhibit a thinner passivation layer; this only begins to build up fully after several preparation cycles. Alkaline cleaning can support this process. The passivation layer represents the necessary level of protection against corrosion, which even so-called "corrosion-resistant" steels only possess to a certain degree.

Instruments undergo cleaning and a final inspection at KARL STORZ prior to delivery.

Upon receipt of the delivery, the customer checks the instruments to ensure that they are complete, function properly and have not sustained any damage in transit.

Should you wish to make a complaint, a description of the procedure is enclosed with every delivery.

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**5.5 Cleaning, disinfection and
sterilization of thermostable
instruments and accessories****Cleaning**

As a rule, cleaning takes place at temperatures of 40 to 60°C for a minimum of 5 minutes. Bear in mind here that these instruments can be prepared with a pH-neutral or alkaline cleaning agent in accordance with national guidelines and recommendations.

Disinfection

Machine thermal disinfection in the washer and disinfectant reaches a maximum temperature of 95°C with an appropriate exposure time (in accordance with EN ISO 15883) and can therefore be used for thermostable instruments.

Sterilization

Steam sterilization is generally carried out at 134 to 137°C and is preferred for thermostable instruments.

**5.6 Cleaning, disinfection and
sterilization of thermolabile
instruments and accessories****Cleaning**

As a rule, cleaning takes place at temperatures of 40 to 60°C for a minimum of 5 minutes. Bear in mind here that the use of alkaline cleaning agents in the preparation of these instruments is not unrestricted, as recommended by the RKI. Process chemicals which have been approved for use on individual instruments/groups of instruments must also be noted.

Disinfection

Machine chemo-thermal disinfection in the washer and disinfectant reaches a maximum temperature of 65°C (< 60°C in the case of flexible endoscopes) and can therefore be used for thermolabile instruments (in accordance with EN ISO 15883). Cleaning/disinfection of thermolabile instruments must, where necessary, be carried out manually.

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Sterilization:

Bear in mind here that thermolabile instruments are sterilized using a low-temperature sterilization method. These methods work with chemical agents at temperatures of between 37 and 75 °C.

WARNING: When carrying out any work on contaminated instruments, the guidelines of the Employer's Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

5.7 Special procedures

Preparation of medical devices, e.g., when (variant) Creutzfeldt-Jakob disease ((v)CJD) is suspected, must be based on the respective applicable national regulations. In clinically verified cases, the used medical devices are generally to be destroyed.

In suspected cases, the reparable medical devices must be marked and stored in isolation until the diagnosis is clarified. Where the suspicion is unfounded, conventional preparation is performed; where the suspicion is confirmed, the instrument set used must be incinerated.

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As a precautionary measure to avoid the transfer of (v)CJKD pathogens from patients without specific symptoms to subsequently treated patients, we recommend alkaline cleaning with an alkalinity of > pH 10 (in accordance with the recommendation of the RKI in Germany) followed by steam sterilization (also with a sterilization time of 18 minutes depending on the nationally applicable regulation).

For Germany: The recommendation of the Robert Koch Institute (RKI) on the preparation of medical devices contaminated with lymphatic tissue or which have come into direct contact with the eye, spinal cord, and brain (CNS) applies.

5.7.1 Thermostable medical devices

- No fixing preconditioning
- Machine cleaning with alkaline cleaner > pH 10
- Thermal disinfection at 93°C (device-specific)
- Steam sterilization 134°C/18 mins.

5.7.2 Thermolabile medical devices

- No fixing preconditioning
- Machine cleaning with alkaline cleaner > pH 10 preferred

Effective disinfection, preferably by machine, e.g., based on peracetic acid sterilization in a validated low-temperature procedure.

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6 Carrying out the preparation

Proper and effective cleaning is essential to the success of the subsequent disinfection and sterilization. It must therefore be carried out with particular care.

Personnel must have the appropriate training and specialist knowledge required for preparation.

They must have in-depth knowledge of the instruments to enable them to carry out the necessary disassembly and care procedures as well as the test for proper functioning.

In certain cases, in-depth consultation of the instrument-specific manuals/accompanying instruction leaflets may be necessary.

Machine processes must be validated and their use preferred. Manual processes must be described as precisely as possible using procedure instructions and validated to prove their efficacy. This is necessary in order to come close to machine processes, which can be standardized, in terms of validity.

Personnel are to wear personal protective equipment (PPE): protective cap to hold back/cover the hair completely, eye protection, mouth and nose protection, protective gown/apron, gloves, suitable moisture-proof shoes.

The process chemicals are selected to suit the medical devices to be prepared. Neutral cleaning agents are available, the effect of which is based on surfactants or quaternary ammonium compounds (QUAT). Enzyme cleaners enhance cleaning and have a liquefying effect which requires sufficient exposure time at the correct temperature.

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Alkaline cleaning is characterized by high levels of efficacy in respect of dissolving protein and fat residues as well as an antimicrobial effect (In Germany: see recommendation of the RKI). An alkalinity (pH value) of > 10 to 12.5 is recommended for the optimal removal of protein, and also to reduce the potential danger of a transfer of vCJK and CJD.

i NOTE: The applicable national laws, regulations and recommendations must be observed.

6.1 Preparation for cleaning and disinfection

The following measures help to ensure that the value and functionality of the instruments are maintained. They facilitate preparation significantly.

- The coagulation surfaces of coagulation instruments should be wiped off from time to time during use, e.g., with a surgical pad moistened with a 3% H₂O₂ solution (hydrogen peroxide solution).
- Heavy soiling must be removed immediately after use whilst still in the operating/procedure room e.g., by wiping off external contaminations and rinsing out working channels before placing (not throwing) into the disposal container.
- Instruments which are found to be defective during use must be visibly labeled and subjected to preparation all the same before sending them to the manufacturer for repair.

! **WARNING:** In order to ensure personal protection and avoid damage to other instruments (isolation, etc.), sharp instruments must be stored separately or fixed in place using specific storage systems.

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6.2 Transportation to the preparation station

Immediately after use, ensure that the instruments are dry and securely stored and then transport them to the preparation station.

i NOTE: In order to ensure safe transportation and storage, specific holding systems are available for delicate instruments (see catalog "HYGIENE, Care, Sterilization and Storage Techniques").

6.3 Disassembly and steps for manual preliminary cleaning

After delivery of the instruments, unload the transport system immediately.

! CAUTION: Use splash protection for eyes, mouth and nose as well as the body (PPE).

Separable instruments must be opened or separated into their individual components as far as possible. To do so, consult the specific manual for the instrument if necessary (see section 3 "Using this manual").

6.3.1 Rinsing off/out

- Clean off any heavy soiling from the external surfaces with the aid of a brush/sponge under cold running water, until the instrument is visibly clean.
- Instruments with working channels, lumina and similar hollow spaces, as well as tubing, must be dismantled and rinsed out prior to cleaning and disinfection (e.g., using a pressure water gun and a suitable rinsing attachment with at least 3 pulses lasting 5 seconds each and a water pressure of at least 2.5 bar), so that patency is ensured and the process chemicals can take effect.
- Suitable brushes may be used to support the process.
- Jointed instruments must be opened to approx. 90°.

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CAUTION: Do not use fixing agents (e.g., alcohol, aldehyde, peracetic acid) when rinsing off and out and ensure that this is not carried out at temperatures above room temperature.

CAUTION: Where cleaning agents or disinfectants are used other than those which have been approved, KARL STORZ cannot provide a guarantee in case of damage.

NOTE: A drying out time of 60 minutes (between use and cleaning) should not be exceeded (the process validations carried out by KARL STORZ are based on a drying out time of one hour).

6.3.2 Incrustations

Incrustations can be dissolved using a 3% H₂O₂ solution. Incrustations which remain on the surfaces of the instrument may cause irreparable corrosion damage and also impede cleaning (see notes in the Red Brochure “Proper Maintenance of Instruments” from the A-K-I (Working Group on Instrument Reprocessing, www.a-k-i.org).

6.4 Auxiliary pretreatment with ultrasound (US)

US (at approx. 35 kHz) can be used as a pretreatment for heavily contaminated instruments.

NOTE: Ultrasound is not a cleaning method! US dissolves contaminations and serves as an aid to machine or manual cleaning and disinfection within the scope of preliminary treatment.

CAUTION: Instruments with optical components (e.g., telescopes and light carriers) are not intended for ultrasound treatment.

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i NOTE: The effect of US can even reach difficult to access hollow spaces. For this reason, all non-separable scissors and forceps, with or without irrigation channels, rail-sheath instruments, suction devices, circular punches, coagulation instruments, and microinstruments must be pretreated with ultrasound.

US has no effect on soft plastic parts because the oscillations are attenuated.

Recommended procedure:

- When preparing the solution, always add the water first and then the chemicals.
- Fill the ultrasound bath with cold, clear water and prepare the solution according to the chemical manufacturer's specifications.
- Before the first treatment, degassing must take place for approx. 15 mins. (idle).
- The ultrasound unit must be covered and then turned on. After degassing, the ground clearance spacer must be inserted.
- Only one sieve tray of instruments can be loaded for each treatment (acoustic shadows).
- The instruments must be completely covered by the solution and air must be removed from any hollow spaces present in order to fill these with solution.
- The ultrasound unit must be covered during the treatments.
- If the ultrasound unit is switched on, nobody may reach inside with their hands.
- The ultrasound treatments should only continue until the soiling has been dissolved, and must not exceed 10 minutes. Exact timings can be found in the instrument-specific instruction manual.
- Prolonged use of ultrasound raises the temperature of the bath and reduces its effect.
- The level of contamination in the bath must be monitored. Contaminated solutions must be replaced in accordance with the chemical manufacturer's instructions (documentation of time, concentration, personnel, temperature).

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Small parts, such as nuts and bolts, may become loose due to the vibration. Following the ultrasound treatment, check that the instruments are complete. If chemicals are added to ultrasound baths, do not use fixing agents.

i NOTE: The maximum treatment time is 10 minutes at approx. 35 kHz. It is imperative that instruments which have been treated in the ultrasound bath are rinsed off/out sufficiently with cold, clear water prior to further preparation, to ensure that residues (dissolved dirt, residual process chemicals) are removed and not carried over into the subsequent cleaning and disinfection procedures.

! CAUTION: Please observe the section on water and steam quality to avoid damage to the instruments.

6.5 Manual cleaning/manual disinfection

6.5.1 Manual cleaning

! CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, temperature, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

i NOTE: For cleaning accessories which are recommended and have been approved by KARL STORZ, including those especially developed for flexible endoscopes and endoscopic instruments, please see the appendix.

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The chemical manufacturer's specifications regarding the preparation, concentration, temperature, exposure time and service life of the cleaning solution must be observed. Prior to immersion in the cleaning solution, rinse the instruments off and out to effect the removal of residues at an early stage (in Germany: as recommended by the RKI). Open jointed instruments to approx. 90°, dismantle separable instruments in accordance with their disassembly instructions and then immerse. Once the exposure time has elapsed, the instruments, which are fully immersed, are cleaned below the surface of the water using suitable cleaning accessories, e.g., brushes. Open and close the joints several times. When cleaning, be sure to reach all external and internal surfaces as far as possible, including those that are difficult to access. Afterwards, the instruments must be rinsed off and through sufficiently under cold, clear running water to ensure the neutralization of the cleaning agent. The cleaned instruments must undergo intermediate drying before disinfection to avoid diluting the disinfectant solution.

i NOTE: Prompt cleaning (within an hour) is necessary in order to prevent contaminations coagulating and drying out.

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6.5.2 Manual disinfection

The chemical manufacturer's specifications regarding the preparation, concentration, temperature, exposure time and service life of the disinfectant solution must be observed.

 **NOTE:** Disinfection can only be carried out after cleaning!

Having been cleaned, rinsed sufficiently with cold, clear water and been subjected to intermediate drying, the instruments are then immersed in the disinfectant solution. Ensure that all external, internal and hard to reach surfaces are covered. The instruments must be completely covered by the solution and any hollow spaces present must be filled with the solution making sure that all air bubbles have been removed.

Non-separable instruments with an irrigation port must be sufficiently rinsed out with disinfectant solution. Be mindful of a flow towards the distal end.

Following the exposure time, the instruments must be sufficiently rinsed off and out several times with completely demineralized water or microbiologically pure/sterile water in order to remove all chemical residues.

 **WARNING:** It is particularly important to use completely demineralized water for the final rinse!

6.5.3 Drying

Once the instruments have been cleaned and disinfected, all of the surfaces, joints, openings, channels, and lumina are then dried completely, preferably with sterile compressed air. The compressed air gun with accessories (Art. No. 27660) is ideal for this purpose.

Instruments with easy access to all surfaces can be dried with an absorbent low-lint cloth.

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6.6 Machine cleaning and disinfection

In order to prepare for machine cleaning, the following measures must be carried out, depending on the degree and type of contamination:

- Rinse the instruments off and out with cold, clear water before placing them in the washer and disinfectant in order to achieve an initial reduction in the volume of contamination being brought into the unit and prevent process chemicals being carried over from upstream cleaning processes.
- Dismantle separable instruments in accordance with their disassembly instructions and position them on a suitable insertion trolley.

As a rule, the most important processes (standardized programs) are preprogrammed in the washer and disinfectant by the manufacturer of the device. Machine cleaning/disinfection and drying in the washer and disinfectant must be preferred to the manual equivalent. The advantage of machine cleaning/disinfection lies in the reduction of risk to personnel, the standardization of the processes and the creation of a basis for a validated preparation process.

i NOTE: *Process validations for cleaning/disinfection are carried out by KARL STORZ with a drying time of 1 hour (at room temperature).*

The properties of the instruments in terms of materials and geometry are taken into account in the selection of machine cleaning and disinfection processes. KARL STORZ recommends the use of suitable storage and holding systems which provide optimal support to the cleaning process. Use of these systems protects the instruments from mechanical damage on the one hand, and facilitates targeted rinsing of the instruments' specific geometries on the other.

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i NOTE: The temperature of the water flowing into the unit should not exceed 45°C, higher temperatures lead to the coagulation of protein and may impair the cleaning performance as a result.

! CAUTION: Ensure that no extreme changes in temperature occur during machine preparation. Flexible and semiflexible telescopes as well as some miniature telescopes must not be heated above 65°C. This also applies to sterilization, transportation and storage. Further details can be found in the specific manuals of the respective instruments.

We recommend using completely demineralized water in all stages of the program. The selection of the process chemicals and their ratios must also be observed precisely. The pH values should preferably be in the neutral range, the possible requirements placed on alkaline cleaning, taking the applicable regulations into account, are to be observed. The use of care and rinse agents in the machine is discouraged, as deposits may form on the instruments or plastics may be damaged by the agents. Local, manual care and lubrication of joints etc. following machine cleaning and disinfection is more effective.

As there is sometimes strong movement of the water, it is necessary to secure all parts in holding devices to avoid damage. To ensure that the instruments are adequately rinsed out, it is essential to connect up the hollow spaces appropriately.

i NOTE: Use the insertion trolley specially provided for use with the washer and disinfectant. Avoid overloading the trays.

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For telescopes, tray holders from the 39501xx series are most suitable (see appendix). In addition, manufacturers of washers and disinfectors offer various retaining mechanisms for different types of instruments. For small parts (stopcocks, seals etc.), small tray holders are available.

When loading the washer and disinfectant, bear the following in mind:

- Avoid unrinsed areas caused by instruments with large surfaces
- When rinsing the inside of instruments which are connected via their LUER lock irrigation port, the jaws are also cleaned well. For certain instruments, specific adapters are available for rinsing off/out. Information concerning this can be found in the accompanying documentation for the products or in the catalog.
- Jointed instruments must be opened wide enough to enable effective rinsing of all accessible surfaces to take place. The use of an opening spring is recommended for this purpose (Art. No. 39100 L).

i NOTE: For the final rinse, use completely demineralized water (to prevent staining, discoloration and corrosion).

i NOTE: Inspect seals regularly for material fatigue (tears, porosity, embrittlement).

i NOTE: Instruments which are still contaminated following machine cleaning/disinfection must be manually cleaned and then subjected to a full preparation procedure once more.

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6.6.1 Machine cleaning/thermal disinfection

Depending on the medical devices to be prepared, machine thermal disinfection should be chosen over chemical or chemico-thermal disinfection as far as possible for ecological and economic reasons.

The A₀ concept (see DIN EN ISO 15883) must be taken into consideration in the on-site process validation carried out by the user.

Depending on the type of machine, thermal disinfection occurs at a temperature of 80-95°C and with an exposure time which is appropriate for both reaching the A₀ value, which is defined as essential, and achieving a level of bacteria reduction of at least 4 log levels.

1.	Prerinse for 1 min. with cold water from the faucet
2.	Empty the bottle
3.	Prerinse for 3 mins. with tempered water from the faucet
4.	Empty the bottle
5.	Clean for 10 mins. with 0.8% Neodisher Septo Clean in completely demineralized water (begin to introduce the chemical from 40°C, then heat up to 55°C)
6.	Empty the bottle
7.	Neutralize for 1 min. using completely demineralized water at a temperature of 20°C
8.	Empty the bottle
9.	Rinse for 2 mins. using completely demineralized water at a temperature of 20°C
10.	Empty the bottle
11.	End of the cleaning process (for validation of cleaning without disinfection)

Table 4: The process used by KARL STORZ for the validation of cleaning

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The subsequent process steps regarding thermal disinfection, neutralization and drying are not carried out within the scope of a validation of the cleaning process.

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Example of a process flow – machine cleaning with thermal disinfection

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Fig. 1: Example of a process flow – machine cleaning with thermal disinfection

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6.6.2 Machine cleaning/chemico-thermal disinfection

The cleaning phase corresponds in principle to that of the thermal procedure. This procedure is used for thermolabile medical devices and avoids temperatures above 65°C. The process step of disinfection, with the aim of a reduction in bacteria of at least 4 log levels, uses chemical disinfectant agents. The information provided by the manufacturer of the washer and disinfectant (for endoscopes) should be observed particularly closely in the case of chemico-thermal disinfection. Specifically:

- Connection of endoscopes with channels/channel systems to the washer and disinfectant for endoscopes via device-specific adapters
- Leakage test/correct connection
- Water quality (see section “General information on preparation”)
- Parameters of preliminary cleaning
- Selection and ratios of the process chemicals (in the cleaning process as well as in the separate disinfection process)
- Process control (irrigation pressure, temperature, exposure times)
- Exposure time and exposure temperature of the process chemicals
- Neutralization (sufficient rinsing afterwards)
- Drying by machine/manually/in a compartment dryer
- Regular inspection of the process quality by studying the prepared product and by studying the device during routine checks (water samples, maintenance etc.)
- Process validation

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Example of a process flow – machine cleaning with chemo-thermal disinfection

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Fig. 2: Example of a process flow – machine cleaning with chemo-thermal disinfection

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6.7 Assembly, inspection and care

Following machine cleaning, disinfection and drying, the instruments must cool down to room temperature.

The cleaned/disinfected instruments must be inspected by visual/tactile means to ensure that they are clean, complete, dry and have not been damaged. For the inspection, the use of a magnifying lamp is recommended.

If residues or contamination are still present, the instruments must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

Damaged or corroded instruments must be withdrawn from use and either sent for repair or replaced. In any guarantee claims, we ask that you send the instrument concerned in to our Service department with a detailed description of the defect (see section on “Brand-new and repaired instruments”).

Deposits on metal surfaces can be removed with polishing wool (Art. No. 27658) if necessary. The instrument must then undergo a complete preparation cycle once more.

i NOTE: *The formation of deposits warrants an analysis of the causes, taking into account the process flows, process chemicals and materials used as well as water quality (see notes in the Red Brochure “Proper maintenance of instruments” from the A-K-I (Working Group on Instrument Reprocessing) www.a-k-i.org).*

Plastic parts and elastomers should be inspected regularly for discoloration, porosity and flexibility. If the parts show differences from their “as-new” condition which are relevant to function, then they must be replaced.

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Joints, threads and glide surfaces should be lubricated locally with instrument oil/special grease for stopcocks prior to packaging and sterilization.

i NOTE: The oil used for this purpose must be suitable for the subsequent sterilization process (silicone-free and paraffin- or white oil-based).

Opening and closing several times ensures the uniform distribution of the lubricant in the crevices between the moving components and thereby prevents damage such as friction corrosion.

Jaws and joints must be checked particularly carefully for damage, as must the insulation on instruments for HF applications. The jaws of scissors should enable smooth cuts and those of forceps should exhibit non-slip edges, be correctly aligned with one another and allow complete opening and closing.

Apply a thin coat of special grease (Art. No. 27657) to sealing and running surfaces, such as stopcocks and plungers. This thin layer of grease not only serves to lubricate the metal surfaces which are in contact with one another, but also makes the stopcocks leakproof. The suitability of oil for use as a sealant is limited.

If the instruments are clean, complete, dry and have been cared for, then they must be subjected to a test for proper functioning.

For this purpose, disassembled instruments must be assembled in accordance with their assembly manual.

In order to test HF instruments and their cables for proper operation, a suitable testing device should be used.

The optical system of rigid endoscopes should be inspected for damage and perfect image transmission. In flexible fiber endoscopes, the fiber bundle responsible for image transmission should also be inspected. In the case of videoscopes, an inspection is only possible if a camera control unit is used.

The inspection of the cold light cable is described in detail in the specific manual.

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6.7.1 Care products

i NOTE: During care procedures, use items from the catalog "HYGIENE, Care, Sterilization and Storage Techniques".

i NOTE: Greased components must be dismantled for sterilization so that steam penetration can be ensured.

Not all care products are suitable for low-temperature processes.

6.8 Packaging and labeling

i NOTE: The instruments must be as free from residue as possible upon submission for packaging and sterilization.

As a basic principle, only clean, dry instruments which are in full working order may be packed for sterilization.

Requirements for correct packaging can be found in the following standards:

- ISO 11607-1:2006 (English version)
Part 1: Requirements for materials, sterile barrier systems and packaging systems
- EN 868 is also applicable: Packaging materials and systems for medical devices which are to be sterilized
- DIN 58953 1-10 for sterilization as a whole (in Germany)
- When using fleece sheets, the folding techniques conforming to DIN 58953-7 (in Germany) must be observed (parallel or diagonal).

6.8.1 Packaging systems

Only standardized packaging materials and systems may be used. Containers, for example, are considered rigid packaging and fleece sheets or film/paper, for example, are considered soft packaging.

i NOTE: National regulations must be taken into account.

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6.8.2 Types of packaging

Single packaging = primary packaging
Double packaging = primary + secondary
packaging

**6.8.3 Requirements of packaging
materials**

- No change during sterilization
- Consistent for the respective sterilization medium
- Suitable for labeling
- Suitable for storage with a permanent barrier function against microorganisms
- Sterile opening possible
- Sufficient protection during transportation

6.8.4 Aspects of packaging to consider

- Suitable packaging material for the sterilization method selected
- Associated process indicators with the specifications of the packaging manufacturer in mind
- Adapt the packaging size to the volume of the goods to be sterilized (maximum fill of the packaging 2/3)
- Packaging undamaged before and after sterilization
- Dry packaging following sterilization
- Sufficient width of the sealing seams, no trapped air or folds
- Labels, imprints, tags for film/paper packaging not on the paper side
- Packaging easy to handle
- Choose an economical and ecological type of packaging

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6.8.5 Securing the instruments in the packaging

KARL STORZ offers holding and storage trays and systems for specific instruments, e.g., rigid and flexible endoscopes, instruments for minimally invasive surgery and freely configurable modular systems (see catalog “HYGIENE, Care, Sterilization and Storage Techniques”).

Systems of this kind facilitate protected storage for the instruments throughout the entire preparation process and are generally suitable for soft packaging and container systems.

6.9 Sterilization

CAUTION: When selecting the sterilization method, it is imperative to consider the compatibility of the medical device with the process selected.

6.9.1 Selection of a suitable sterilization method

As a rule, steam sterilization at 134°C to 137°C (132°C in the USA) is the most gentle and economical method of sterilization. This also applies to autoclavable rigid endoscopes.

The low-temperature steam and formaldehyde sterilization process (LTSF), EO sterilization (6-100%), H₂O₂ sterilization (e.g., ASP STERRAD®, Steris® AMSCO® V-PRO™ 1) and chemical low-temperature sterilization with peracetic acid (Steris® System 1/1E) are all possible for thermolabile medical devices.

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Key:

- ✓ Recommended method for sterilization
- A Telescopes, prismatic light deflectors and cameras with the marking “AUTOCLAV” can only be steam sterilized using the prevacuum method – do not use flash sterilization.

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Table 5: Assignment of various KARL STORZ products to suitable sterilization methods

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KAMERAKÖPFE / CAMERA HEADS / CABEZALES DE LA CÁMARA	STERILISATION / ESTERILIZACIÓN							Desinfektion / Disinfection / Desinfección
	EO Gas / EO gas / Gas OE	STERRAD® (Dampf H ₂ O ₂ , Plasma) / (Steam H ₂ O ₂ , plasma) / (vapor H ₂ O ₂ , plasma)			Steris®			
		100S (US) / 100S (US) / 100S (US)	100S "Kurz"-Zyklus (Int'l) / 100S "Short" cycle (Int'l) / 100S Ciclo "breve" (Int'l)	NX "Standard"-Zyklus / NX "Standard" cycle / NX Ciclo "Estándar"	100NX "Standard"-Zyklus / 100NX "Standard" cycle / 100NX Ciclo "estándar"	V-PRO™ 1 (Dampf H ₂ O ₂) / V-PRO™ 1 (steam H ₂ O ₂) / V-PRO™ 1 (vapor H ₂ O ₂)	System 1 (Peressigsäure-Lösung) / System 1 (peracetic acid solution) / System 1 (solución de ácido peracético)	
H3-Z (mit blauer Kabelbuchse) / H3-Z (with blue cable socket) / H3-Z (con toma azul para cable)	✓	✓	✓	✓	✓	✓	✓	✓
H3 (mit blauer Kabelbuchse) / H3 (with blue cable socket) / H3 (con toma azul para cable)	✓	✓	✓	✓	✓	✓	✓	✓
H3 (schwarzer Kabelbuchse) / H3 (black cable socket) / H3 (con toma negra para cable)	✓	✓	✓				✓	✓
S3/S1 (mit blauem zylindrischen Kabelschuh) / S3/S1 (with blue cylindrical cable lug) / S3/S1 (terminal de cable cilíndrico, azul)	✓	✓	✓	✓	✓	✓	✓	✓
S3/S1	✓	✓	✓				✓	✓
P3/P1	✓	✓	✓	✓	✓		✓	✓
F3	✓	✓	✓	✓	✓		✓	✓
D1	✓	✓	✓	✓	✓		✓	✓
TELECAM	✓	✓	✓				✓	✓
DCI®-1/DCI®-2	✓	✓	✓				✓	✓
TRICAM®	✓	✓	✓				✓	✓
A3/A1	✓	✓	✓	✓	✓	✓	✓	✓
TELECAM A/C	✓	✓	✓				✓	✓
TRICAM® A/C	✓	✓	✓				✓	✓

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Table 6: Assignment of various KARL STORZ products to suitable low-temperature sterilization methods

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Information regarding efficacy and material compatibility is available from the manufacturer of the instrument or the distributor of the method. All of the models listed in Table 7 have been validated in respect of their efficacy in (high-level) disinfection, with exemplary use of the products CI-DEX (glutaraldehyde, 2.4%) and Resert XL HLD (2.0% H₂O₂).

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**6.9.2 Overview of important points
regarding sterilization**

The instruments must be freed from all organic and inorganic materials as well as residues from process chemicals prior to sterilization. The sterilization of all internal and external surfaces is only possible if the instruments have been dried, cleaned and disinfected.

i NOTE: In case this is not expressly indicated elsewhere, all thermostable medical devices can be steam sterilized in their fully assembled state using the fractionated prevacuum procedure described in this manual.

Pull only lightly on screws and fasteners. Only close handles with catches using the first notch, so that no tensile cracks are caused during sterilization. Heating and cooling during sterilization may cause tension in instruments with closed catch locks which may lead to a reduction in tensile forces or even cracks at the joints. The methods described are valid where maintenance is carried out as prescribed and the sterilization device is validated.

During sterilization, ensure that all parts are stored safely. To do so, items from the catalog “HYGIENE, Care, Sterilization and Storage Techniques” may be used.

The quality of the water used for steam generation should be taken into account in order to avoid damage to the instruments and sterilizers (see section on “Water and steam quality”).

The steam sterilization method recommended in some countries, i.e., 18 minutes at 134-137°C, can be used where necessary (e.g., for prophylactic prion inactivation).

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6.9.3 Steam sterilization

Pure, saturated steam is used for sterilization. Steam is economical, non-toxic, environmentally friendly, easily controllable and, under the correct conditions, leaves no residue.

An increased temperature and the increased pressure resulting from it facilitate a shorter sterilization time whilst maintaining the same elimination rate of the microorganisms which are potentially present, including spores (121°C/15 mins., 134°C, 4-5 mins.).

The energy transfer which occurs when the steam condenses on the cooler surfaces of the instruments, combined with the moisture penetrating the cell structure, has a biocidal effect.

Requirements for correctly carrying out sterilization using saturated steam can be found in the following standards:

- EN 285 Steam sterilizers – Large sterilizers (from 1 StU upwards) – requirements and test procedures
- EN 13060 Small steam sterilizers
- EN ISO 17665-1 for routine monitoring

The standard applies to sterilization methods for medical devices whereby moist heat is used as a sterilizing medium. The area of application is not limited to sterilization methods using saturated steam, but rather also encompasses the steam-air mixture method, water sprinkling method and methods involving immersion in water. There are specific requirements for

- the elements of the quality management system
- the characterization of the sterilizing medium, method and equipment
- product and procedure definitions
- validation including the principal elements – installation qualification, operation qualification, performance qualification, monitoring and confirmation of validation
- routine monitoring, checks and maintenance of the efficacy of the method
- approval of the product following sterilization

A detailed, informative appendix includes explanations of the essential requirements.

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6.9.4 Steam sterilization methods

- Fractionated prevacuum (repeated evacuation of the chamber alternating with injections of steam)
- Prevacuum (single evacuation of the chamber)
- Gravity displacement (use of gravity – heavy air is displaced downwards and forced out of the sterilization chamber by lighter steam)
- Fractionated flow process (repeated filling of the chamber with steam until a prescribed operating pressure is reached)
- Steam injection process (simultaneous evacuation of the chamber and injection of steam)

6.9.5 Fractionated prevacuum

The fractionated prevacuum procedure must be used for instruments with lumina, since it offers the highest degree of penetration of the lumina and thereby safety, due to the way in which the air is repeatedly removed (fractionation).

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It consists of three basic phases:

- Conditioning phase
- Sterilization phase
- Drying phase

The conditioning phase consists of the repeated removal of air (evacuation) from the sterilization chamber.

The temperature rise time (heating up) then follows. This time is needed in order to heat the product to be sterilized up to the prescribed sterilization temperature.

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Afterwards, plateau time (sterilization) begins. This includes 3 phases; A, B and C.

A = Equalization time: interval between the sterilization temperature being reached in the chamber and the sterilization temperature being reached in all parts of the load.

B = Holding-elimination time: prescribed exposure time at the sterilization temperature for the elimination of those germs which are expected to occur most frequently.

C = Safety margin: the proportion of the exposure time accounted for by the different contaminants and equalization times of the products to be sterilized.

Drying takes place after the sterilization phase. The condensate present is vaporized and conducted away via the vacuum line. This is enhanced by the creation of a vacuum. Observance of the maximum load weight of 10 kg per sterilization unit facilitates the drying phase. The permissible residual moisture in accordance with DIN EN 285 amounts to 1% of net weight for textiles and 0.2% of net weight for metals.

After the drying phase, the sterilized product may be removed from the sterilization chamber. The product will still be hot and should not be placed directly onto cold placement surfaces or directly in a draft, because the sterile product will then cool down very quickly, which causes the residual moisture from sterilization to condense once more in the cold packaging. This would result in moist packaging. Moist sterile products may not be placed into storage.

i NOTE: *Long lumina, crevices and channels, as well as assembled instruments can only be sterilized without restrictions with the fractionated prevacuum method.*

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Chronological pressure/temperature profile of a fractionated prevacu-
um procedure

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Fig. 3: Chronological pressure/temperature profile of a fractionated prevacuum procedure

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6.9.6 Prevacuum

Some evacuations did not take place, which results in lower efficacy in crevices and lumina.

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Fig. 4: Chronological pressure/temperature profile of a prevacuum procedure

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6.9.7 Gravity displacement

Use of gravity, because heavy air is displaced downwards and forced out of the sterilization chamber by light steam.

CAUTION: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, crevices or stopcocks.

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Fig. 5: Chronological pressure/temperature profile of a gravity displacement procedure

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6.9.8 Flash sterilization (so called “ready to use” sterilization)

This takes place without a conditioning or drying phase.

CAUTION: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, crevices or stopcocks (e.g., scissors).

CAUTION: Flash sterilization is a method which, for hygiene reasons, is only used in exceptional cases. National regulations must be observed.

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Schematic representation of a flash sterilization procedure

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Fig. 6: Schematic representation of a flash sterilization procedure

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6.9.9 Steam quality

The necessary steam is produced from water (EN 285). Drinking water contains dissolved salts and minerals, which are partially vaporized along with the water and would lead to undesirable changes to or deposits on surfaces. This is why steam is created from completely demineralized water.

Pure steam: from demineralized water, chemically pure and pyrogen-free.

Saturated steam: water is vaporized depending on pressure and temperature. At 3.1 bar and 134°C, the pure saturated steam has zero liquid content.

In addition to the steam quality, the condition of the technical equipment is also decisive. The stainless steel steam generator prevents the steam created from absorbing corrosion particles or metal traces and bringing them into the sterilization chamber. The steam carrying lines must also be made from stainless steel in order to avoid contamination of the steam in transportation. Old lines or imperfect assembly may allow rust particles into the sterilization chamber and damage the product to be sterilized there.

i NOTE: The steam must satisfy the requirements of DIN EN 285 and ISO EN DIN 17665 in order to avoid discoloration and staining. Any cloths which are used to pack the instruments must not contain detergent residues, since this also causes discoloration.

i NOTE: In the case of sterilization with formaldehyde or H₂O₂ “plasma” sterilization, the instruments can be used immediately after the sterilization cycle. An additional airing time is not necessary.

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6.9.10 Low-temperature processes

Processes of this type are used when the medical devices concerned cannot be exposed to the high temperatures used in steam sterilization. Each of the processes named below uses specific packaging which is coordinated with the process.

- i** **NOTE:** Ensure that the sterilization agent can easily access the product to be sterilized during sterilization (contact points).
- i** **NOTE:** For greased and lubricated surfaces, low-temperature gas sterilization is not possible because the sterilization agent is only able to penetrate surfaces of this type to a limited extent.
- i** **NOTE:** Observe national laws and regulations for the gas sterilization methods.
- i** **NOTE:** Regarding alternative H_2O_2 sterilization methods, information for the approval of individual products should be requested directly from the manufacturer.

6.9.11 Hydrogen peroxide gas (H_2O_2)

These processes are suitable for medical devices whose material compatibility is assured. Limitations in respect of efficacy, e.g., due to the special geometry of the medical device to be sterilized (long, narrow lumen) or crevices, must be taken into account during use. The user must observe the specifications of the unit, instrument and packaging manufacturers before sterilization is carried out.

Low-temperature sterilization systems based on H_2O_2 are available from various manufacturers. These enable the sterilization of instruments which are sensitive to temperature, such as flexible endoscopes, amongst others.

- i** **NOTE:** The compatibility of the medical device to be sterilized with the respective procedure should be cleared with the manufacturer. The different limitations on the efficacy of the various H_2O_2 gas systems must also be considered.

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CAUTION: KARL STORZ recommends using only one cleaning agent and one sterilization method. Using several different procedures over a period of time can influence the external appearance and functionality of the medical device concerned. Since it is not possible to calculate the degree of material attack caused by the interactions between chemicals in cleaning, disinfection and sterilization processes, consistent use of only one method for cleaning, disinfection and sterilization is recommended.

NOTE: Sterilization with H₂O₂ may cause cosmetic changes to the instrument, however these have no influence on functionality or hygiene. Cosmetic changes to plastic parts and coated metal surfaces may also occur as a result of the upstream cleaning and disinfection processes.

CAUTION: When using disinfectants based on peracetic acid, a strong chemical interaction with hydrogen peroxide may cause serious irreversible material damage.

CAUTION: For flexible endoscopes, the pressure compensation cap (11025 E or Z10455) must be attached to the connection port for the leakage test before sterilization with H₂O₂. Do not use the yellow transport cap.

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- ASP STERRAD®
- Steris® AMSCO® V-PRO™ 1

WARNING: KARL STORZ advises that instruments which are sterilized in H₂O₂ must be disassembled prior to sterilization and put together afterwards in a sterile environment. The sterility of separable instruments cannot be assured if they are sterilized in a fully assembled state using H₂O₂ systems. This is particularly applicable if care procedures with instrument oil have been carried out.

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CAUTION: All instruments must be dried thoroughly before they are placed in the H₂O₂ sterilizer. Loads which contain moisture may cause the cycle to abort.

CAUTION: Only packaging approved by the respective manufacturer of the H₂O₂ device may be used.

CAUTION: KARL STORZ flexible gastrointestinal videoscopes are not compatible with the hydrogen peroxide gas sterilization process.

CAUTION: Each deviation from the recommended parameters for H₂O₂ sterilization must be validated by the user. Sterility is only assured if the manufacturer's specifications are observed.

NOTE: The classification of instruments by KARL STORZ as compatible with the STERRAD®, AMSCO® V-PRO™ 1 or other H₂O₂ sterilization processes was verified by at least 100 cycles in each case.

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6.9.13 ASP STERRAD®

 **NOTE:** KARL STORZ flexible endoscopes without lumina can be sterilized using the STERRAD® 50, 200 or 100S.

 **CAUTION:** Only use instrument trays which are expressly approved for use with STERRAD®. These trays are specially designed and validated for hydrogen peroxide sterilization.

 **CAUTION:** Sterilization is not possible on surfaces which have been greased and lubricated.

 **WARNING:** Flexible endoscopes with lumina may only be sterilized using the STERRAD® 50, 200 or 100S under certain conditions (use of a booster) whilst taking national regulations into account.

 **WARNING:** Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

 **NOTE:** To ensure that the medical device concerned can be sterilized using the different STERRAD® units, refer to the STERRAD® Sterility Guide (www.sterradsterilityguide.com).

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STERRAD® NX and 100NX systems (worldwide requirements)

Worldwide requirements	Inside diameter	Length	Special requirements
PE/PTFE tube	1 mm or greater	350 mm or shorter	NX only: Process with standard circulation
	1 mm or greater	1000 mm or shorter	For NX: Process with extended circulation
Flexible endoscopes with PE/PTFE lumina	1 mm or greater	850 mm or shorter	Flexible single channel endoscopes only (NX= extended circulation; 100NX only= FLEX circulation)
Stainless steel lumina	1 mm or greater	150 mm or shorter	NX only: Process with standard circulation
	2 mm or greater	400 mm or shorter	NX only: Process with standard circulation
	1 mm or greater	500 mm or shorter	NX only: Process with extended circulation
	0.7 mm or greater	500 mm or shorter	100NX only: Process with standard circulation

NX/100NX = max. 10 stainless steel lumina per load

Table 8b: Restrictions on effective sterilization with reference to instrument geometry

NOTE: These data were provided by ASP (2011)

CAUTION: For more detailed information, please consult the user handbook of the respective device.

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6.9.14 Steris® AMSCO® V-PRO™ 1

Various product categories were tested for material compatibility and validated in respect of efficacy of sterilization. They are suitable for sterilization without significant influences on the material or function. More detailed information regarding the current list of approved products is available from Steris®.

CAUTION: Sterilization is not possible on surfaces which have been greased and lubricated.

6.9.15 Ethylene oxide (EO)

WARNING: Instruments must be dismantled for sterilization with ethylene oxide (EO) and put together again afterwards in a sterile environment.

WARNING: If KARL STORZ instruments are sterilized with ethylene oxide in their fully assembled state, the sterility of the instruments cannot be guaranteed.

CAUTION: For flexible endoscopes, the pressure compensation cap (11025 E or Z10455) must be attached to the port for the leakage tester before sterilization with EO. Do not use the yellow transport cap.

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Ethylenoxid: ≤ 20 mg

Ethylenchlorohydrin: ≤ 12 mg

The procedure specified below is used by KARL STORZ for the validation of EO sterilization.

PARAMETERS (in the chamber)			
Temperature (setpoint value):	55°C (131°F)		
Humidity:	≥ 70% RH		
Vacuum (setpoint value):	1.3 psia 0.09 bar		
Adjusted hold time:	30 mins.		
STERILIZATION PARAMETERS			
Chemical:	100% ethylene oxide		
Temperature (setpoint value):	55°C (131°F)		
Humidity:	≥ 70% RH		
Adjusted hold time:	30-45 mins.		
Product load configuration:	①	②	③
EO gas concentration:	725 ± 30 mg/L	725 ± 30 mg/L	735 ± 30 mg/L
EO gas exposure time:	60 mins.	120 mins.	180 mins.
VENTILATION PARAMETERS			
Ventilation time:	12 hrs.		
Ventilation temperature:	51-59°C (120-130°F) setpoint value: 55°C		

Table 9: Process parameters for sterilization with 100% ethylene oxide. Maximum acceptable amount (ANSI/AAMI/ISO 10993-7) of residue following EO sterilization:

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Load configuration ①:

The 60-minute EO cycle is recommended for the sterilization of surfaces on instruments which have no lumina (without a sterilization tray). These are packed in two layers of single-ply soft packaging made from polypropylene or an equivalent material.

EO gas concentration:	725 ± 30 mg/L
EO gas exposure time:	60 mins.

Load configuration ②:

The 120-minute EO cycle is recommended for instruments (with and without lumina: placed in a sterilization tray and without a sterilization tray where instruments have lumina). These instruments are packed in two layers of single-ply soft packaging made from polypropylene or an equivalent material.

EO gas concentration:	725 ± 30 mg/L
EO gas exposure time:	120 mins.

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Load configuration ③:

The 180-minute EO cycle is recommended for the camera head product family. The camera head (with or without sterilization tray) is packed in two layers of single-ply soft packaging made from polypropylene or an equivalent material.

EO gas concentration:	735 ± 30 mg/L
EO gas exposure time:	180 mins.

Ethylene oxide:	25 ppm*/250 ppm**
Ethylene chlorohydrine:	25 ppm*/250 ppm**
Ethylene glycol:	250 ppm*/5000 ppm**

Table 12: Maximum acceptable levels (as per 21 CFR 812.100; US Federal Register, Vol. 43, No. 122, June 23 1978) of residues following EO sterilization

* for instruments which come into contact with blood.

** for instruments which come into contact with skin and mucous membrane.

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CAUTION: For gas sterilization with ethylene oxide, the airing times stipulated by the device manufacturer must be observed for instruments, due to gas absorption of the materials. The airing times depend crucially on the processing procedure of the EO sterilizer (concentration, process control).

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**6.9.16 Low-temperature steam and
formaldehyde sterilization (LTSF)**

Steam and formaldehyde (FO) serve as a sterilization agent in a fractionated vacuum. The fully automatic process control means that the release of FO can be avoided by working in a closed system. If an error occurs during the procedure, desorption and internal ventilation take place before the unit can be opened. The process temperature is between 50 and 60°C. The concentration of FO in the stock solution is 2-3%. The FO solution is automatically drawn from closed storage vessels. Residues are dissolved from the surface of the sterilized medical devices using a steam wash. The condensed mixture of water and FO is disposed of via a vacuum pump. In this process, the FO is only present in a strongly diluted form, so free, gaseous FO could only escape to a very slight extent.

i NOTE: Sterilization with FO is tested and approved by KARL STORZ in respect of material compatibility. However, validation in respect of efficacy was not carried out due to the large numbers of processes used in the global market. Therefore, the user must validate his specific processes individually in collaboration with the device manufacturer. KARL STORZ has had exclusively the aforementioned sterilization methods validated with regard to a reduction in microorganisms of the order of 10^{-6} SAL.

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6.9.17 Chemical low-temperature sterilization with peracetic acid

CAUTION: For flexible endoscopes, the pressure compensation cap (11025 E or Z10455) must be removed from the port for the leakage test before sterilization with peracetic acid.

Chemical sterilization (with the Steris® System 1/1E) is a special method which is particularly suitable for thermolabile medical devices, e.g., flexible endoscopes. Apart from telescopes with an eye-piece drive mechanism, all KARL STORZ products are suitable for sterilization with Steris®. During sterilization, color changes to the surfaces of the aluminum alloys may occur, which do not, however, affect the functioning of the instruments.

Refer to the instruction manual of the Steris® SYSTEM 1/1E for suitable sterilization parameters and use of the appropriate Quick Connect Kit (QCK). Endoscopes with channels can only be sterilized with the QCK which has been directly allocated to the respective model. This information should be requested from the device manufacturer.

Endoscopes without channels can be sterilized using the Steris® System 1/1E.

CAUTION: Sterility is only assured if the manufacturer's specifications for the Steris® System 1 are followed. This also applies to the use of the sterilization containers.

CAUTION: Chemical sterilization is not possible on surfaces which have been greased and lubricated.

NOTE: Parts with lumina must be prepared together using the Steris® QCK. All parts must be thoroughly cleaned before sterilization using the Steris® PROCESS™, just as it applies to the other sterilization procedures.

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KARL STORZ products which are compatible with and validated for the Steris® PROCESS™:

- Rigid telescopes without working channel, without eyepiece drive mechanism
- Semi-flexible telescopes without working channel
- Flexible hysteroscopes* and cystoscopes* for adults
- Flexible ureterscopes*, cystoscopes* for pediatrics and choledochoscopes* which were made after July 1, 1997
- Immersible video cameras
- Ureter light probes
- Light cables

* The flexible endoscopes are only compatible and validated if sterilization is carried out in the Steris® container for flexible endoscopes using the Steris® QCK.

KARL STORZ product categories which are compatible in terms of material, but were not validated in terms of sterility*:

- Rigid and semi-flexible telescopes with working channel
- Flexible bronchoscopes
- Flexible gastroscopes
- Working elements, proctoscopes and resectoscopes
- Trocars/sheaths
- Reusable tubing sets
- Insulated and non-insulated TAKE-APART® and CLICKLINE instruments
- Insulated and non-insulated forceps, scissors, punches and suture instruments
- Coagulation, irrigation and suction instruments

* The sterility of the parts which are described as “compatible in terms of material, but not validated” cannot be guaranteed if the Steris® PROCESS™ is used. These parts can however be gas-sterilized or steam sterilized.

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7 Maintenance

7.1 Troubleshooting/fault finding

i NOTE: Also refer to the section "Water and steam quality".

In cases of material damage, an analysis of the causes must be undertaken in collaboration with the device manufacturers, the manufacturers of the process chemicals and the internal technical service department. In particular, this involves the detailed analysis of each individual step in the process:

- Application in the OR (residues of pharmaceuticals, disinfectants, local anesthetics)
- (Safe) transportation (following drying/disposal)
- Impure side
- Preparation for cleaning
- Manual cleaning
- Machine cleaning (loading, connection, process chemicals)
- Disinfection
- Drying
- Inspection, maintenance, care
- Packaging (storage)
- Sterilization
- Handling (transportation, storage)

The results of the analysis then flow into the respective quality management system of the operator.

i NOTE: In case of changes to the material and surfaces, the Red Brochure "Proper maintenance of instruments" from the A-K-I (Working Group on Instrument Reprocessing) www.a-k-i.org provides indications of possible causes.

Removal of deposits on metallic surfaces: Deposits on metal surfaces can be removed with DURAGLIT polishing wool (Art. No. 27658) if necessary. The medical device must then be subjected to a complete preparation cycle once more.

Replacement parts: For this purpose, use the catalog "Replacement Parts".

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7.2 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

7.3 Important information

For reasons of hygiene and infection prevention, telescopes, instruments and devices must be cleaned, disinfected and, where necessary, sterilized before they are sent in for repair.

KARL STORZ reserves the right to return contaminated instruments/devices to the sender.

7.4 Repairs

Defective instruments must be repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

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7.5 Repair and replacement program

KARL STORZ maintains a repair and replacement warehouse which is normally adequate to ensure prompt replacement of damaged telescopes and instruments. Under the repair and replacement plan, you receive an identical as-new instrument and are only charged the repair costs for the defective medical device.

In Germany, you can refer the repairs directly to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

Service hotline: +49 7461/708 980
E-mail: technicalsupport@karlstorz.com

In other countries please contact your local KARL STORZ branch or authorized dealer.

8 Directive Compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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9 Appendix

9.1 Extract from the catalog “HYGIENE, Care, Sterilization and Storage Techniques”

Quality standards in the instrument cycle

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The objective is to create a reusable instrument, the materials, design and quality of which allow preparation using a suitable process that takes all applicable norms and laws into consideration, while maintaining the full functionality of the instrument.

Optimal handling throughout the instrument cycle goes further still and complements the instrument with individually coordinated accessories for optimized storage during the entire process – from the operating theater to transportation to the central sterile supply department, right down to the individual preparation steps.

Intelligent storage systems cut handling expenses. A reduced amount of repacking guarantees maximum protection for the instruments, optimizes operating theater procedures and ensures the desired results of the preparation.

KARL STORZ achieves this objective by continuously developing new storage systems and improving existing ones according to the demands of medical professionals.

i NOTE: The following list is merely an extract. A complete list of cleaning accessories can be found in the catalog “HYGIENE, Care, Sterilization and Storage Techniques”.

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





Tray holders for the cleaning, sterilization and storage of endoscopes

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstucio st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-it-iv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlyoftens 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gilwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

KARL STORZ Endoskop Magyarország Kft.
Toberek utca 2. fsz. 17/b
HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
E-Mail: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl
Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4
041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
E-Mail: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.*
Patriarhou Grigoriou E' 34
54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
E-Mail: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**
Gedik Is Merkezi B Blok
Kat 5, D 38-39, Bagdat Cad. No: 162
Maltepe Istanbul, Turkey
Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK
Derbenyevskaya nab. 7, building 4
115114 Moscow, Russia
Phone: +7 495 983 02 40, Fax: +7 495 983 02 41
E-Mail: info-ru@karlstorz.com

TOV KARL STORZ Ukraine
Obolonska naberezhna, 15
building 3, office 3
04210 Kiev, Ukraine
Phone: +380 44 42668-14, -15, -19
Fax: +380 44 42668-20
E-Mail: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG
Representation Office
Sabit Orudschow 1184, apt. 23
1025 Baku, Azerbaijan
Phone: +99 450 613 30 60
E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –
East Mediterranean and Gulf (Offshore) S.A.L.
Spark Tower 1st floor
Charles Helou St., Horch Tabet – Sin El Fil
Beirut, Lebanon
Phone: +961 1 501105, Fax: +961 1 501950
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoskopy Kasachstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & CO. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klontoy, Wattana,
10110 Bangkok, Thailand
Phone:
Fax:
E-Mail:

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2086
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyo-da-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongquiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang Intenational Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyo-da-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

 **KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen
Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com