Instruction Manual
Manuel d'instruction
Manual de instrucciones

Modell / Model / Modèle / Modelo

VIDEOSCOPE

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Instruction Manual

Manuel d'instruction

Manual de instrucciones

Modell / Model / Modèle / Modelo

VIDEOSCOPE

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1.1 Important note for users of KARL STORZ devices and instruments

Thank you for choosing a product by KARL STORZ, a name you can trust. You and your organization have decided in favor of a modern, high-quality item of equipment.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of this device from KARL STORZ. All essential details of the equipment and all actions required on your part are clearly presented and explained.

We thus ask that you read this manual carefully and that you keep it available for ready reference in a convenient and conspicuous location near the equipment.

KARL STORZ is continuously working on the enhancement of all its products. Please understand that changes to the scope of supply, design, equipment, and technology are possible for this reason. Therefore, no claims can be derived from the specifications, illustrations and descriptions contained in this information.

© All product pictures, photos and descriptions and texts are the intellectual property of KARL STORZ SE & Co. KG.

Utilization and reproduction by third parties require express authorization by KARL STORZ SE & Co. KG.

All rights reserved.

1.2 Explanation of symbols

Caution

Observe the instruction manual

Cleaning instructions

Repair instructions

Control of pollution caused by electronic information products (RoHS)

The device is subject to the requirements of the WEEE-Directive (see section 8.4 Disposal).

1 Généralités	1 Generalidades
1.1 Remarque importante	1.1 Indicaciones importantes
1.2 Signification des pictogrammes	1.2 Explicación de los pictogramas

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2.1 Warning symbols

Before using the device, please read this instruction manual thoroughly, particularly the chapter "Warnings and safety instructions", so as to avoid any risk to yourself or your staff.

WARNING: A WARNING indicates that the personal safety of the user or other persons may be endangered. Disregarding a WARNING could result in injury to the user or other persons.

CAUTION: A CAUTION indicates that particular service procedures or precautions must be followed to avoid possible damage to the unit.

NOTE: A NOTE contains additional information on the safe operation of the device.

2.2 General warnings and safety instructions

WARNING: Do not use this device for the inspection of bodily openings in humans or animals. This could result in the death of the person or animal. Keep the device out of the reach of children.

WARNING: Danger of an electric shock! On contact with live parts, the videoscope may be energized. Never touch live parts with the videoscope.

WARNING: Danger of an electric shock! Do not open the handle! Only have your equipment serviced by the manufacturer, KARL STORZ, or by manufacturer-authorized personnel. Any opening of the device components by unauthorized personnel will void the warranty.

WARNING: Do not insert the insertion part into an object to be inspected which is in operation or live. Otherwise, the distal end or sheath may be damaged if it gets stuck inside the inspected object, or similar. Also, an electric shock may occur if the videoscope comes into contact with the object.

2	Consignes de sécurité	2	Instrucciones de seguridad
2.1	Symboles d'avertissement	2.1	Símbolos de advertencia
2.2	Avertissements et consignes générales de sécurité	2.2	Advertencias e instrucciones de seguridad generales

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! **CAUTION:** When cleaning the device, always observe the instructions in the instruction manual.

! **CAUTION:** Never use the device without a lens installed or if a lens component is loose or damaged. Otherwise, the device or the object of examination may be damaged.

! **CAUTION:** Acclimatization time! If the device is stored/transported at temperatures below +10 °C before being put into operation, condensation may build up on or in the device. In this case, it is essential to allow for an acclimatization time of approx. 3 hours before use.

! **CAUTION:** Do not kink or squeeze the sheath under any circumstances and never place it under strain by forming overly tight coils which do not observe the minimum bending radius. It is best to suspend the videoscope by the handle for storage or keep it in the carrying case supplied, and to avoid pulling or carrying it by the sheath or connected cables.

i **NOTE:** The device must always be shut down before the lens is taken off or changed. This ensures that you will not be blinded by the illumination light. Furthermore, the videoscope must only be stored with the lens fitted to prevent the ingress of dirt and damage to the lens.

2.3 Intended use

The videoscope is designed for the endoscopic inspection of technical assemblies, such as turbines, motors, cavities in buildings and vehicle bodies, etc.

The videoscope may only be used with accessories and consumables which have been designated as suitable by KARL STORZ. Use of the product in fields other than those described above is not permitted for safety reasons.

Unauthorized conversions or modifications to the videoscope are not allowed for safety reasons and will immediately invalidate all warranty rights.

2.4 User qualification

The product must only be used by persons authorized by KARL STORZ who are suitably qualified and have received instruction in the use of the videoscope.

2.5 Safety measures at the location of use

- Workstations and rooms in which a videoscope with a LASER reference measurement system is used must be identified with a warning sign in compliance with the Safety of Laser Products standard IEC 60825-1.
- Prior to using the device, the user should confirm the safe function and proper condition of the device.

2 Consignes de sécurité	2 Instrucciones de seguridad
2.2 Avertissements et consignes générales de sécurité	2.2 Advertencias e instrucciones de seguridad generales
2.3 Utilisation conforme à l'emploi	2.3 Uso previsto
2.4 Qualification de l'utilisateur	2.4 Cualificación del usuario
2.5 Mesures de sécurité sur le lieu d'utilisation	2.5 Medidas de seguridad en el lugar de la aplicación

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- Any penetration of fluid MUST be avoided.
- Do not store fluids on or above equipment.

2.6 Safety measures during use

 **WARNING:** The videoscope must not be operated in potentially explosive atmospheres.

 **WARNING:** Risk of injury! Always wear suitable safety goggles when changing the interchangeable lenses on videoscopes with or without a laser measuring system.

 **WARNING:** DANGER OF BLINDING! Never look into the distal end of a videoscope or into the light beam of the LED light source which can be connected. Therefore, never look into the bright light outlet of the videoscope from a short distance.

 **WARNING:** RISK OF BURNS! The high-intensity light emitted by the light source may cause high temperatures at the tip of the videoscope.

 **WARNING:** RISK OF BURNS! Any contact with the tip of the light connector immediately after operation may cause burns. Do not place the connected videoscope on top of combustible materials, such as textiles, or in close proximity of textiles saturated with combustible fluids, since these may be ignited by the heat developing at the videoscope tip.

 **CAUTION:** The device must not be exposed to precipitation. Operation in condensing atmospheres (damp locations) is not permitted.

 **CAUTION:** Do not insert the sheath of the videoscope into an object to be inspected if the object is not within the permissible operating temperature range. Otherwise, this may result in damage to or impaired performance in the device.

 **CAUTION:** When removing the sheath of the videoscope from the cavity to be inspected, ensure that the exterior of the sheath is not damaged by burrs, edges, or excessive exertion of force.

 **CAUTION:** When removing the videoscope from the cavity to be inspected, ensure that the deflectable sheath end is set to a deflection of 0° (aligned straight). Otherwise, the videoscope may be damaged.

 **CAUTION:** When angling the distal tip, observe the hard stops and do not force the bending. Otherwise, the videoscope may be damaged.

 **CAUTION:** The device must not be operated at increased air pressure (diving stations) or low air pressure (at heights of over 3000 m above sea level).

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CAUTION: Handle the device in such a way that critical optical components, such as the light cable connection, distal end or display, are not damaged or contaminated.

NOTE: Protect the distal end in particular from impacts, blows and squashing, since it contains movable deflection joints.

2.7 Safety measures for storage and disposal

CAUTION: Do not store the device at ambient temperatures above +60 °C or below -10 °C.

CAUTION: Store the device in such a way that critical optical components, such as the light cable connection, distal end or display, are not damaged or contaminated.

CAUTION: Observe the following points when stowing the device away in the carrying case or removing it from the case:

- Ensure that the videoscope sheath is not stored in a twisted position and use the recess provided in the carrying case insert. Do not stow the videoscope away with an angled sheath and always release the brakes.
- Only stow the videoscope away once the distal end or light source has cooled down.
- When removing the videoscope sheath from the recess in the carrying case, do not apply any excessive force and do not remove the main unit by pulling on the insertion part.
- Do not actuate the rotary dial or brakes when the videoscope is stowed in the carrying case.

CAUTION: Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

NOTE: Avoid deflecting the distal end when the sheath is rolled up in the carrying case for storage. This may result in damage to the mechanical deflection.

2.8 General information on the use of lasers

WARNING: Danger of glare! Never look into the distal end of a videoscope. Always shut the device down before changing the lens.

The laser classification system is partly based on the 'blink reflex', a defensive reaction of the eye when suddenly exposed to strong light. This protective reflex does not occur if you deliberately look into a laser beam (or bright light source). For this reason, never look into a laser intentionally, even if it is classified as eye-safe. Exception: examination by an ophthalmologist, etc.

Intentionally looking into strong light sources and LEDs from a short distance can also endanger the eyes. Therefore, never look into the bright light outlet of the videoscope from a short distance.

Do not point a laser beam at people or animals.

2 Consignes de sécurité	2 Instrucciones de seguridad
2.6 Mesures de sécurité pour l'emploi	2.6 Medidas de seguridad durante la aplicación
2.7 Mesures de sécurité pour le stockage et l'élimination	2.7 Medidas de seguridad en la conservación y la eliminación
2.8 Informations d'ordre général sur l'utilisation des lasers	2.8 Generalidades sobre el empleo de láseres

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2.8.1. Changing lenses

The device must always be shut down before the lens is taken off or changed. This ensures that you will not be dazzled by either the laser or the illumination light.

Furthermore, the videoscope must only be stored with the lens fitted to prevent the ingress of dirt and damage to the lens.

2.8.2. Specifications

Reference measuring system
(without interchangeable lenses) and
Multipoint measuring system
(without interchangeable lenses)

Laser class in accordance with IEC EN 60825-1 standard:	laser class 3R
Laser power at the outlet:	< 5 mW
Wavelength:	639 nm
Modulation:	modulated 25/30 Hz, CW
Beam divergence:	approx. 2 mrad

2.8.3. Modifications

Dismantling or modifying the device may cause laser radiation of a higher class to be released. Therefore never carry out any modifications, repairs or maintenance on the device.

2.8.4. Laser warning plate

The information and warning plate is located on the front of the laser housing. The plate must not be removed.

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



2.8.5. Laser outlet apertures

The interchangeable lens fitted to the device divides the laser beam into 49 divergent measuring beams (multipoint reference measurement system), which then emerge from the rectangular glass window. The white illumination light also emerges here.

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With the 90° interchangeable lens, the arrangement of the optical components is different. Nevertheless, the two measuring beams still emerge at the rectangular glass window.

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2.9 Power supply

The laser of this videoscope is powered by the KARL STORZ TECHNO PACK® X / Xe / T LED. No other method of operation or power supply from other equipment is permitted. The supply connector must not be plugged into other devices (e.g., a PC) even if it would fit the socket.

2 Consignes de sécurité	2 Instrucciones de seguridad
2.8 Informations d'ordre général sur l'utilisation des lasers	2.8 Generalidades sobre el empleo de láseres
2.9 Alimentation électrique	2.9 Suministro de energía

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2.10 Instructions and extracts from the laser standard IEC 60825-1

10.7 Reflections

Great care should be taken to prevent unintended reflections of beams from laser equipment of Class 3R, Class 3B and Class 4. Mirrors, lenses and beam splitters should be firmly attached and should only be moved in a controlled manner while the laser is emitting.

Care should be taken to ensure that no unintended reflections of beams of laser equipment of Class 1M and Class 2M occur on surfaces that may focus the beam.

Reflective surfaces that appear diffuse may in fact reflect a significant proportion of the radiation, especially in the infrared spectral range. That may possibly be more dangerous at relatively great distances than might be expected for purely diffuse (Lambertian) reflections.

Potentially hazardous reflections occur on all surfaces of transparent optical components such as lenses, prisms, windows and beam splitters.

Potentially hazardous radiation may also pass through some reflective optical components such as mirrors (e.g. infrared radiation that passes through a reflector for visible radiation).

12.5 Laser equipment in laboratories and at the workplace

12.5.1 Laser equipment of class 1M, class 2, class 2M and class 3R

Precautionary measures are only necessary to prevent looking into the beam for a significant period; for Class 2 and Class 2M, momentary (0.25 s) radiation in the wavelength range of 400 nm to 700 nm, as may occur by accidentally looking into it, is not considered dangerous. However, the laser beam should not be intentionally directed at people. The use of optical aids (e.g. binoculars) together with laser equipment of Class 1M, Class 2M and Class 3R may increase the risk to the eyes. Additional precautionary measures for laser equipment of Class 1M, Class 2M and Class 3R are set down in 12.6.2.

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12.6.2 Laser equipment of class 1M, class 2M and class 3R for purposes of measurement, guidance and leveling

Only qualified and trained persons should be entrusted with setting up, adjusting and operating laser equipment.

Areas in which these lasers are used should be marked with an appropriate laser warning sign.

If it is possible in practice, mechanical or electronic means should be used for guiding the laser.

Preventative measures should be taken in order to ensure that people do not look directly into the beam. (Looking into the beam for lengthy periods may be dangerous). Even looking into the beam with optical aids (theodolites, etc.) may be dangerous, in particular with lasers of Class 1M and Class 2M that do not meet the conditions of Table 10, and should not be permitted, unless a laser protection officer expressly consents to it.

The laser beam should be closed off at the end of its course defined by its purpose. It should in any event be closed off if the hazardous course of the beam (safe distance) stretches beyond the laser monitoring area.

The course of the laser beam should be well above or below eye level, where this is possible in practice.

Precautionary measures should be taken to ensure that the laser beam does not fall unintentionally on surfaces that reflect like a mirror (in particular not on flat or concave reflective surfaces). Unused lasers should be stored in locations to which unauthorized persons have no access.

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3.1 Videoscope 8 mm / 6.3 mm / 6 mm / 5 mm / 4.2 mm / 4.4 mm / 3.8 mm

The videoscope is offered as a set.

The videoscope is equipped with a CCD or CMOS color camera chip which is integrated in the distal end of the flexible sheath. The 6.3 mm videoscope also has an opening (working channel) on the distal end, through which tools can be inserted for grabbing or magnetic grasping.

The flexible sheath contains the cable connection between the color camera chip and the handle, as well as the integrated fiberoptic light cable. The ergonomic handle has three integrated control keys with a triple function. For positioning the distal end, the device has 2 hand wheels that enable the deflection of the distal end in four directions.

The integral cable, which is fixed to the handle, branches at the T-piece into the fiberoptic light cable (with plug-in contact for the light outlet point on the cold light projector) and the video cable with connector for the camera control unit or mobile documentation unit.

The camera control unit or the mobile documentation unit contains the power supply and performs video processing of the camera signal. There, various standard video signals are output for further processing, e.g. for a monitor.

The interchangeable lenses are mounted on the tip of the videoscope with a safety locking system and allow the direction of view to be deflected.

The quick-action mount for the optional snap-on flat screen monitor enables fast and ergonomic attachment of the monitor to the handle itself.

The system case ensures the videoscope set's safe transport and protected storage.

3 Description de l'appareil	3 Descripción del equipo
3.1 Vidéoendoscope 8 mm / 6,3 mm / 6 mm / 5 mm / 4,2 mm / 4,4 mm / 3,8 mm	3.1 Videoendoscopio 8 mm / 6,3 mm / 6 mm / 5 mm / 4,2 mm / 4,4 mm / 3,8 mm

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
- ① CLICK4MOVE® deflection
- ② Integrated cable
- ③ Connection for monitor and universal adapter
- ④ Laser box
- ⑤ Light connector
- ⑥ Video cable
- ⑦ Control keys
- ⑧ Handle with bend protection
- ⑨ Sheath
- ⑩ "Laser" connector
- ⑪ Deflectable sheath end
- ⑫ Interchangeable lens
- ⑬ Protection cap for video connector
- ⑭ Video connector

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- ① Attachment knob for carrying strap
- ② Locking button for snap-on monitor
- ③ Quick-action mount for monitor
- ④ Integrated cable
- ⑤ T-piece
- ⑥ Video cable
- ⑦ Light connector
- ⑧ Protection cap for video connector
- ⑨ Video connector
- ⑩ Control key, black
- ⑪ Control key, silver
- ⑫ Control key, blue
- ⑬ Handle with bend protection
- ⑭ Sheath
- ⑮ Deflectable sheath end
- ⑯ Interchangeable lens
- ⑰ Locking lever
- ⑱ Large hand wheel
- ⑲ Rotary dial
- ⑳ Small hand wheel

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
4.1 Unpacking the device

The videoscope comes in a compact case which is easy to transport. Carefully unpack the videoscope and remove the device and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints immediately with the manufacturer or supplier. If possible, retain the original packing materials for later use; these may be useful if the unit has to be transported.

4.2 Scope of supply

The scope of supply is indicated on the shipping note.

***KARL STORZ continuously reviews the set configurations and will modify and adapt them according to current technical requirements.**

4.3 Controls

4.3.1 Videoscope 8 mm / 6.3 mm / 6 mm / 4.2 mm / 4.4 mm

- ① CLICK4MOVE® deflection of the movable tip to the right or the left
- ② CLICK4MOVE® deflection of the movable tip upward or downward
- ③ Control key, black
- ④ Control key, silver
- ⑤ Control key, blue

4 Mise en service*	4 Puesta en funcionamiento*
4.1 Déballage	4.1 Desembalaje
4.2 Équipement fourni	4.2 Suministro
4.3 Organes de commande	4.3 Elementos de control

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

- ① Large hand wheel for deflecting movable tip upward or downward
- ② Small hand wheel for deflecting movable tip to right or left
- ③ Rotary dial for locking movable tip in the respectively set right or left deflection position
- ④ Locking lever for locking the movable tip in the respectively set top or bottom deflection position
- ⑤ Release button for snap-on flat screen
- ⑥ Attachment knob for carrying system
- ⑦ Control key, black
- ⑧ Control key, silver
- ⑨ Control key, blue

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

- ① Large hand wheel for deflecting movable tip upward or downward
- ② Locking lever for locking movable tip in the respectively set up or down deflection position
- ⑤ Release button for snap-on flat screen
- ⑥ Attachment knob for carrying system
- ⑦ Control key, black
- ⑧ Control key, silver
- ⑨ Control key, blue

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
5.1 Preparing the device for use

5.1 Preparing the device for use

NOTE: The TECHNO PACK® X / Xe / T LED basic unit must be shut down when the camera connector and the laser box are inserted or removed.

5.1.1 Videoscope without measuring function with TECHNO PACK® X / Xe / T LED

- Remove the connector protection cap.

- Insert the light cable connector into the light outlet socket of the TECHNO PACK® X / Xe / T LED as far as it will go.
- Switch on the TECHNO PACK® X / Xe / T LED and light source.

- Insert the video connector in the camera socket of the TECHNO PACK® X / Xe / T LED.

5.1.2 Videoscope with measuring function with TECHNO PACK® X / Xe / T LED

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5 Operation

5.1 Preparing the device for use

- Insert the laser box and park support as far as they will go into the light extraction socket of the TECHNO PACK® X / Xe / T LED.

- Insert the connector into measuring box serial port of the TECHNO PACK® X / Xe / T LED.

- Insert the video connector in the camera socket of the TECHNO PACK® X / Xe / T LED.
- Switch on the TECHNO PACK® X / Xe / T LED and light source.

CAUTION: The unprotected connector contact of the video scope can be damaged if the connector falls down. It is therefore urgently recommended to protect the connector against damage with the protection cap fastened to the cable when it is unplugged.

5.1.3 Operation with 7" VGA flat screen for TECHNO PACK® X/Xe

- Insert the monitor holder into the quick-action mount of the video scope from top to bottom. To release, press down on the locking button and pull out the monitor holder.
- Insert light the cable connector into the light outlet socket of the TECHNO PACK® X / Xe as far as it will go.

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- Insert the connector into the serial port for the measuring box of the TECHNO PACK® X / Xe. In the case of a videoscope with measuring function, insert the connector into the laser box.
- Switching on the TECHNO PACK® X / Xe, light source and flat screen.

5.1.4 Operation with CLICK4MOVE® 7" flat screen for TECHNO PACK® X/Xe

- Push the locking switch on the monitor holder to the left and insert the monitor into the quick-action mount of the videoscope. To release the monitor, push the locking switch to the left and pull out the monitor holder.
- Insert the light cable connector into the light outlet socket of the TECHNO PACK® X / Xe as far as it will go.
- Insert the connector into the serial port for the measuring box of the TECHNO PACK® X / Xe. For a videoscope with measuring function, insert the connector into the laser box.
- Switch on the TECHNO PACK® X / Xe, light source and flat screen.

NOTE: The CLICK4MOVE® adaptor can also be purchased without flat screen.

5.2 White balance

The white balance adjusts the color reproduction of the system to match the color temperature of the cold light projector used. The white balance remains in memory after the base unit (TECHNO PACK® X/Xe only) is turned off. When the unit is switched on again a new white balance is only necessary if a different light source, a different lamp, a different light cable or a different videoscope is used.

- Point the videoscope towards a white surface while the camera control unit is turned on. Ensure that no objects with other colors can be seen in the image section.
- If the white balance is correct, the screen display is briefly inverted. If this does not happen, the white balance was not completed correctly. Either too much or too little light was received. In this case, increase or, as the case may be, decrease the distance between the videoscope tip and the white surface.
- Repeat the white balancing.

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5.3 Installation/removal of interchangeable lenses and protective cover

The Installation/Removal of the interchangeable lenses is done differently.

- With 8 mm / 6 mm / 4.2 mm / 4.4 mm videoscopes, the interchangeable lenses are fastened to the tip of the videoscope with a double safety screw thread.
- With 6.3 mm / 5 mm / 3.8 mm videoscopes, the interchangeable lenses are fastened to the tip of the videoscope with a safety bayonet lock.

CAUTION: In the case of incorrect operation and excessive use of force, the videoscope can be damaged.

5.3.1 Installation of an interchangeable lens with the videoscope 8 mm / 6 mm / 4.4 mm / 4.2 mm

The tip, or distal end, of the videoscope is equipped with two fine threads which are attached one in front of the other (main thread and safety thread). These prevent the interchangeable lens or protective cover from being lost in the space to be inspected.

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CAUTION: If the working steps described below cannot be carried out without using excessive force, the threads are probably not aligned correctly (are slanted). In this case, the protective cover or interchangeable lens must not be screwed onto the tip of the videoscope using excessive force or tools, as the thread may otherwise be damaged.

1. The lens head ① is located on the tip of the videoscope. This must be removed before the interchangeable lens is installed. To remove the lens head ① the tip is held at the base ② and unscrewed in the direction of the arrow.
2. Once you have passed the main thread ③ pass the safety thread ④ in the same way.
3. The interchangeable lens consists of the threaded sleeve ⑤ and the lens head ⑥. Two pins ⑦ are located at the tip of the 6 mm videoscope and one pin at the tip of the 4.2 mm and 4.4 mm videoscopes. The pins must engage in the notch in the lens head ⑥. The two marker pins ⑧ on the interchangeable lens and the tip of the videoscope are available as aids for this purpose and should be positioned in a line. The lens head can be moved back and forth within the threaded sleeve ⑤ to align the instrument guide.
4. The interchangeable lens is attached to the tip of the videoscope. The threaded sleeve ⑤ is pushed onto the safety thread ④ without slanting it and screwed on in the direction of the arrow. Once the safety thread ④ has been passed, the interchangeable lens on the threaded sleeve is pushed against the main thread ③ as far as it will go and, likewise, screwed on in the direction of the arrow. The interchangeable lens has been adequately secured if the lens head can no longer be moved within the threaded sleeve.

CAUTION: Do not tighten the threaded sleeve using excessive force or tools, as this may damage the thread. The interchangeable lens has been adequately secured if the lens head can no longer be moved within the threaded sleeve.

5.3.2 Removing an interchangeable lens from the videoscope 8 mm / 6 mm / 4.4 mm / 4.2 mm

An interchangeable lens is removed in the reverse order and using the opposite directions of rotation to the installation of an interchangeable lens described above. After the interchangeable lens has been removed, the protective cap for the 0° telescope or another interchangeable lens must be installed again.

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5.3.3 Validation tool V 0600 V for 6 mm videoscope

Checking with the validation tool V 0600 V ensures the measuring accuracy of the entire measuring system. The validation tool V 0600 V is used to check the configuration setting between the videoscope, interchangeable measuring lens and adjustment file (TECHNO PACK® X / Xe).

To check, please follow the instructions of the validation software.

- NOTE:** The entire measuring system (videoscope, interchangeable measuring lenses and TECHNO PACK® X / Xe) should be checked before each measuring procedure.
- NOTE:** Checking with the validation tool V 0600 V is only possible in connection with TECHNO PACK® Xe software version 1.18.4.4. and higher.
- NOTE:** No validation check is required for the TECHNO PACK® T LED, Software Version 29 or higher because this takes place as part of automatic level recognition.

5.3.4 Installation of an interchangeable lens with a videoscope 6.3 mm / 5 mm / 3.8 mm

1. The locking piece ① is located on the tip of the videoscope. It is spring-loaded towards the tip of the videoscope. For installing the interchangeable lens, the locking piece ① must be pressed all the way down against the spring pressure and kept in this position throughout the entire process.
2. With the other hand, the interchangeable lens ④ is aligned to the tip of the videoscope so that the base of the marking triangle ② is in line with the guide piece ③. The interchangeable lens ④ is then pushed onto the tip of the videoscope while overcoming a spring pressure.

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3. To lock the bayonet lock, the pushed-on interchangeable lens ④ is turned in the direction of the arrow until it noticeably engages. The locking piece ① can then be released.

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5.3.7 Removal of an interchangeable lens from a videoscope 6.3 mm / 5 mm / 3.8 mm

1. The locking piece ① is located on the tip of the videoscope. It is spring-loaded towards the tip of the videoscope. For installing the interchangeable lens, the locking piece ① must be pressed all the way down against the spring pressure and kept in this position throughout the entire process.
2. With the other hand, the interchangeable lens ④ is pressed back, overcoming a spring pressure and turned against the arrow direction of the marking triangle ② as far as it will go.
3. The interchangeable lens can now be removed from the tip. Afterwards, the locking piece ① can be released.

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


5.3.8 Installation of 5 mm interchangeable lenses with the interchangeable tool V 0500 AC

The interchangeable tool V 0500 AC is for better handling of the installation or removal of 5 mm interchangeable lenses on 5 mm videoscopes.

On the tip of the videoscope is the locking piece ①, which is permanently held in place towards the tip of the videoscope.

Place the interchangeable tool so that the lever handle can press down the locking piece ①.

When pressed, the interchangeable lens can be installed or removed as described previously.

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
5 Betrieb

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5 Operation

5.4 Controlling the viewing direction

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

5.4 Control of the viewing direction

The viewing direction of the videoscope can be modified by deflecting the moveable tip in four directions.

By turning the larger of the two hand wheels ① located on the handle, the movable tip of the videoscope is deflected downwards or upwards. The set direction of view downwards or upwards can be locked in place between the large hand wheel and the handle with the locking lever ②.

When the smaller of the two hand wheels ③ attached to the handle is turned, the movable tip of the videoscope is deflected to the left or to the right. The set viewing direction to the left or the right can be locked in place with the rotary dial ④ on the small hand wheel.

❗ **NOTE:** With the manual CLICK4MOVE® controls, the movable tip remains in its position after each deflection. Locking with a locking lever is not necessary.

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6.7 Control keys

5.5 Control keys

The assignment of the control keys on the handle of the videoscope depends on whether it is used together with a camera control unit or the KARL STORZ TECHNO PACK® X / Xe / T LED. All keys have two functions, which are distinguished by pressing the respective key(s) for a long or short time. A selected function is therefore not activated until the key(s) are released.

Control keys	Operation with the KARL STORZ camera control unit Techno Videoscope	Operation with the TECHNO PACK® X / Xe	Operation with the TECHNO PACK® T LED
Blue button (short pressure < 2 sec.)	activates the ACC1 function for controlling additional devices (e.g. a monitor etc.) via the camera control unit	activates the 'Still Image' function. Press it again and select either 'Save' (another short press on the blue button) or 'Do not save' (short press on the black button)	activates the 'Still Image' function. Actuating the control key again causes the photo to be saved
Sliver button (short pressure < 2 sec.)	activates the ACC2 function for controlling additional devices (e.g. a monitor etc.) via the camera control unit	activates the 'Carry out measurement' function (only with videoscopes with multipoint reference measurement system)	activates the 'Start measurement / End measurement' function (only with videoscopes with multipoint reference measurement system)
Black button (short pressure < 2 sec.)	Not assigned	activates the 'Laser on / off' function (only with measuring videoscopes; with other videoscopes, an error message may appear)	activates the 'Laser on / off' function (only with measuring videoscopes; with other videoscopes, an error message may appear)
Blue button (long pressure > 2 sec.)	activates the 'White balance' function	activates the 'Video recording' function. Actuating the control key again ends the video recording	activates the 'Video recording' function. Actuating the control key again ends the video recording
Sliver button (long pressure > 2 sec.)	activates the shutter function (increase in the shutter if the space being examined is inadequately lit) to manual 1/50 s or 1/60 s. This causes a slow-motion effect during the image reproduction	activates the 'Shutter and Integration Mode' function (reduces/increases the shutter if the space being examined is inadequately lit/overlit). This results in delayed image reproduction. Control is via the blue and black buttons. To exit, press the silver button for a short time	activates the 'Increase image brightness' function. Each press on the button increases the image brightness up to the maximum, after which it starts again at the minimum
Black button (long pressure > 2 sec.)	switches between center- and edge-weighted automatic exposure NOTE: In the case of pipe inspections, the space being examined in the center of the image is usually very dark. For this reason, overexposure can occur in the peripheral areas in the case of center-weighted exposure measurement. To avoid this effect, it is possible to switch over to edge-weighted exposure measurement. When edge-weighted automatic exposure is used, the upper and lower measurement windows are briefly displayed in succession with a greater brightness than the surroundings. When center-weighted automatic exposure is used, the center measurement window is briefly displayed	activates the 'Laser brightness' function (only with measuring videoscopes – with other videoscopes, an error message may appear). Control is via the blue and black buttons. To exit, press the silver button for a short time	activates the 'Increase laser intensity' function. Each press on the button increases the image brightness up to the maximum, after which it starts again at the minimum (only with measuring videoscopes – with other videoscopes, an error message may appear)

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6.1 Videoendoskop

6.1 Videoscope

Videoscopes	Ø 3.8 mm	Ø 4.2 mm	Ø 4.4 mm	Ø 5 mm	Ø 6 mm	Ø 6.3 mm*	Ø 8 mm
* with working channel diameter 1.43 +/- 0.026							
Color system	PAL and NTSC						
Color camera sensor technology	CCD						
Control keys	TECHNO PACK® setting: Still image, video recording function, manual shutter setting, measuring function, measuring laser settings or HighGain and zoom function						
Integrated cable	T-plug, integrated cable 2.3 m, video cable 0.7 m						
Deflection type	CLICK4MOVE®*			CLICK4MOVE®*			
Deflection [°] ± 10° * with 360° inspection in Generation II CLICK4MOVE®	Double deflection up/down 170°/90°	Quadruple deflection 4x 170°	Quadruple deflection 4x 170°	Quadruple deflection up/down 180°/120° right/left 150°/150°	Quadruple deflection up/down 180°/180° right/left 150°/150° For 4.5 m up/down 160°/160° right/left 130°/130° For 7.5 m 4x 90°	Quadruple deflection up/down 180°/120° right/left 120°/120° For 4.5 m up/down 160°/120° right/left 120°/120° For 7.5 m 4x 90°	Quadruple deflection up/down 120°/180° right/left 90°/90° For 7.5 m 4x 90°
Temperature sensor for MULTIPOINT videoscopes			 as of 80°	 as of 80°	 as of 80°	 as of 80°	 as of 80°
Lens lock	Bayonet system	Screw thread	Screw thread	Bayonet system	Screw thread	Bayonet system	Screw thread
Safety thread	Triple safety system	Double safety system	Double safety system	Triple safety system	Double safety system	Double safety system	Double safety system
Sheath Ø	Ø 3.8 mm ±0.1	Ø 4.0 mm ±0.2	Ø 4.0 mm ±0.2	Ø 4.9 mm ±0.05	Ø 6.0 mm +0.3	Ø 6.3 mm -0.3	Ø 8.0 mm ±0.2
Distal end, sheath	Ø 3.8 mm	Ø 4.2 mm -0.04	Ø 4.4 mm -0.04	Ø 4.9 mm ±0.05	Ø 6.0 mm +0.3	Ø 6.3 mm -0.3	Ø 8.0 mm ±0.2
Working lengths	1000 mm 2000 mm	2000 mm 3000 mm	2000 mm 3000 mm	2000 mm 3000 mm 4500 mm	2000 mm 3000 mm 4500 mm 7500 mm	3000 mm 7500 mm	2000 mm 3000 mm 4500 mm 7500 mm
Minimum bend radius							
Sheath	40 mm	60 mm	60 mm	35 mm	35 mm	35 mm	80 mm
Lens	10 mm	10 mm	10 mm	15 mm	15 mm	15 mm	25 mm
Rigid length visible distal end 0°	15 mm ±3	19.5 mm	18 mm	21 mm ±3	18 mm ±3	21 mm ±3	20 mm ±3
Rigid length visible distal end 90°	18 mm ±3	22 mm	22 mm	26 mm ±3	23 mm ±3	26 mm ±3	25 mm ±3
Temperature range	Sheath in water +10 °C – +30 °C Sheath in air -10 °C – +80 °C						

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Error description	Possible causes	Solution
Complete failure of device	AC power supply failed	Check AC power supply
	Defective fuse	Replace fuse
	Inadequate connection of AC power supply connector and device outlet	Press AC power supply connector firmly into device outlet
	Defective power cord	Check power cord and AC power supply connector
No image/faulty image	No connection between video output of the camera control unit and the monitor	Check the video cable connection and replace cable, if necessary
	Video cable is not plugged correctly into video outlet	Plug in video cable correctly
	Monitor set incorrectly	Check monitor settings
	Videoscope electronics defective	Send videoscope/camera control unit for repair
	Monitor defective	Check monitor
	Interchangeable lens has become loose	Screw interchangeable lens on tight or secure it in place
Fuzzy image, stripes, streaks, etc.	Videoscope telescope or interchangeable lens is soiled	Clean with cotton swab and alcohol solution
Incorrect measurement	Interchangeable lens has become loose	Refit and lock interchangeable lens
	Incorrect adjustment data selected	
The deflection of the distal end is stiff	The sheath is excessively coiled or twisted or stuck	Roll out the sheath and disconnect it from the voltage
	The distal end has insufficient room for maneuvering within the object of examination	

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8.1 Maintenance and care

Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

Care

To care for the videoscope, no special measures are necessary.

- Keep all the system components clean and treat them carefully.
- Wipe down regularly with a damp cloth.

Cleaning

The exterior surfaces of an endoscope are cleaned with a soft cloth soaked in a little isopropyl alcohol (100 %).

It is recommended to clean the optical surfaces (lens, eyepiece and light cable port) with a cotton swab soaked in isopropyl alcohol (100 %).

! **NOTE:** Cleaning immediately after use, especially in the case of contact with foreign materials, is urgently recommended. Virtually imperceptible soiling on the tiny lens of a videoscope will cause picture quality to deteriorate. In addition, removing caked-on dirt may easily damage the lens coating.

8.2 Repairs, the KARL STORZ repair program

On request, KARL STORZ Endoskope carries out repairs after first providing a cost estimate. When sending your unit for repair, please describe the fault and state the type and serial number of your unit. You will find this information on the identification plate of the device.

Preventive maintenance is not essential.

! **CAUTION:** The videoscope must not be opened. Defective units must only be repaired by KARL STORZ Endoskope or by persons authorized by KARL STORZ Endoskope. Only original parts from KARL STORZ Endoskope may be used for repairs.

8 Remise en état, entretien, élimination	8 Mantenimiento, conservación y eliminación de desechos
8.1 Maintenance et entretien	8.1 Mantenimiento y conservación
8.2 Réparations	8.2 Reparaciones

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KARL STORZ repair program

To prevent interruptions at your organization due to repairs, KARL STORZ Endoskope will provide units on loan. In Germany, you can obtain advice at the following address:

KARL STORZ SE & Co KG
 Industrial Group
 Dr.-Karl-Storz-Straße 34
 78532 Tuttlingen, Germany
 Phone: +49 7461 708-926
 Fax: +49 7461 78912
 E-mail: industrialgroup@karlstorz.com
 Internet: www.karlstorz.com

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8.3 Chemical resistance

In case of a short contact time, the insert can be immersed into the following substances:

Aviation fuels: kerosene, e.g. JetA-1 and TS-1, fuel for jet planes, e.g. JP4 / JP5 / JP8.

Motor oils and gear oils: mineral oil, e.g. hydraulic oil, partially synthetic oil, synthetic oil, e.g. Skydrol, Jet Oil 2

Transformer oil with inhibitors

Fuels: diesel fuel and bio-diesel, gasoline

Isopropyl alcohol

1.1.1 Trichlorethane

Water

mild saline solution (5 %)

The handle must not be inserted into the liquid but may come into contact with the chemicals mentioned above.

Clean the guide tube or the handle when they were in contact with the above-mentioned liquids or other liquids before using them again.

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8.4 Disposal

This device has been marked in accordance with the European Directive on waste of electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

 Control of pollution caused by electronic information products (China RoHS).

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US7486805; US7435218; US2008/9672A1
US2011/88498

9.1 Compliance with directives

Along with the TECHNO PACK® products, the electronic products described in this manual are provided with CE markings and meet the requirements of the following directives:

- LVD (Low Voltage Directive) 2006/95/EC
- EMC Directive 2004/108/EC (Electromagnetic Compatibility).

Laser Protection Class acc. IEC 60825-1:
Laser Protection Class 3R

9.2 Warranty

This device is warranted for a period of two years following initial operation. For further information please refer to our General Business Conditions. In the case of unauthorized opening, repairs, or modifications of the device by non-authorized persons, KARL STORZ Endoskope shall not be liable for the operational safety of the device. In case of any such action, all warranty claims shall be void. Any damage to the unit resulting from incorrect operation is not covered by the warranty.

9.3 Technical documentation

On request, the manufacturer, KARL STORZ, will provide available circuit diagrams, detailed parts lists, descriptions, adjustment instructions, and other documents useful for qualified personnel of the user when repairing device parts that have been classified as repairable by KARL STORZ. However, the provision of technical documents does not represent authorization by the manufacturer, KARL STORZ, to open and repair the device, even by technically trained personnel. Excluded from this are procedures described in the text of this manual. We reserve the right to make changes in the design, in particular in connection with the further development and improvement of our devices.

9.4 Property rights

The product is protected in the USA by at least one of the following US patents

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