


Scope Warmer for single use




Scope Warmer for Single Use

For warming the scope for fog-free work



Endoscope fogging represents a serious challenge in everyday clinical practice. The average room temperature in the OR is 20 °C while the patient's body temperature is approximately 37 °C. This differentiation between room and body temperature often leads to condensation of the scope when it is inserted inside the body, which considerably impedes the surgeon's view. To proceed with surgery, cleaning of the endoscope is necessary which is associated with prolonged operative duration and thus increased costs.



Product Features

- Warmer activated to heat endoscope via a metal plate
- Suitable for all rigid endoscopes
- Available in two different sizes (12 x 20 cm and 12 x 36 cm)
- No other heat source (electricity, warm water, etc.) required
- Product is sterile on delivery and is ready for immediate use
- Two scopes can be heated in the warmer at the same time

Using the Scope Warmer


1. Place the scope on the pad and close the warmer with the attached press buttons.



2. Bend the silver metal plate to activate the warmer.
For optimal warming results, it is recommended to activate the warmer approximately 5 minutes before placing the scope inside.



Available Sizes




- 100120-15 **Scope Warmer**, 12 x 20 cm, sterile, for single use, package of 15, for use with all standard rigid endoscopes, featuring optical systems with lengths that fit into the pad



- 100121-10 **Scope Warmer**, 12 x 36 cm, sterile, for single use, package of 10, for use with all standard rigid endoscopes, featuring optical systems with lengths that fit into the pad

Notes





75 Years

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com