

July 23, 2021

Security Notice CVE-2021-34527 update

Microsoft (CVE-2021-34527) Windows Print Spooler Remote Execution Vulnerability

Dear Sir or Madam,

KARL STORZ has evaluated our current, legacy and branded products. The table below shows what products are affected by the Windows Print Spooler Remote Execution Vulnerability and what support KARL STORZ is providing for these affected products.

Current Supported Products

Product number	Product Version(s)	Operating System	Solution	Patch Availability Date:
STREAMCONNECT	StreamConnect Content Management Solution Versions 4.0 and newer	Windows Server 2019 Windows Server 2012	Install Microsoft patch. KS verified patch	July 1 2021

Mitigation

- Apply Microsoft Security Patch for the PrintNightmare Vulnerability. KARL STORZ has tested and verified that the patch does not cause any adverse effects to the product once applied.

Impact

- None.

Product Risk Level (After mitigation/remediation applied)

- **Low**

Employ Good Network Hygiene Practices

- Ensure data has been backed up and stored according to your individual processes and disaster recovery procedures.
- Execute updates to malware protection, where available.

Customers that maintain Virtual Machine (VM) system patches independent of KARL STORZ patch release should ensure these actions are performed to maintain the correct security posture of the system(s).

For more information and instructions on how to apply the Microsoft security patch, please follow the link below:

<https://msrc.microsoft.com/update-guide/vulnerability/CVE-2021-34527>

If you have questions, please contact your KARL STORZ sales representative, or contact us at the following address: techsupport@karlstorz.com