



Set for Endoscopic Saphenous Vein Harvesting



Special Features:


- 45° angled telescope for optimal visualization of the operative field
- Reuseable
- Bipolar instruments and clip applicators enable dissection of side branches
- Open System: no carbon dioxide embolism
- Ergonomic handle
- Smoke evacuation
- The slender design allows greater freedom of movement



+ free

**1 x Endoscopic Vein Retractor,
FREIBURG model**

**2 x Vein Dissector
1 x Fiber Optic Light Cable**



EndoFax®
Endo @ Mail®

**KARL STORZ
GmbH & Co. KG**

Mittelstraße 8
78532 Tuttlingen/Germany
Postfach 230
78503 Tuttlingen/Germany

Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105

E-Mail: info@karlstorz.de
Web: www.karlstorz.com

CARDIO - 36 - E, 11/2011
1/2

Components of the KARL STORZ Bonus Offer

49205 FDZ **Endoscopic Vein Retractor**, FREIBURG model, for harvesting the saphenous vein, enhanced telescope guiding channel, distal width 25 mm, working length 27 cm, with integrated guide in handle for fiber optic light cable, **autoclavable**, for use with HOPKINS® II Telescope 49205 FA

including:

Cleaning Adaptor, for telescope channel


49205 FA **HOPKINS® II Forward-Oblique Telescope 45°**, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: black

free:

49201 VE **Vein Dissector**, blunt, distal angled to left, **autoclavable**, size 3 mm, working length 41 cm

49201 VD **Vein Dissector**, blunt, distal angled to right, **autoclavable**, size 3 mm, working length 41 cm

495 NA **Fiber Optic Light Cable**, with straight connector, diameter 3.5 mm, length 230 cm



EndoFax[®]
Endo @ Mail[®]

M-G36

To
KARL STORZ GmbH & Co. KG
Mittelstraße 8, D-78532 Tuttlingen/Germany

Via fax to
+49 (0)7461 708-105

I am interested in the KARL STORZ Bonus Offer

- Please send me a non-binding offer
- Please contact me regarding detailed information on this offer

Name, First Name _____ Tel. _____

Hospital / Office _____

Address _____

E-mail _____

**KARL STORZ
GmbH & Co. KG**

Mittelstraße 8
78532 Tuttlingen/Germany
Postfach 230
78503 Tuttlingen/Germany

Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105

E-Mail: info@karlstorz.de
Web: www.karlstorz.com

CARDIO - 36 - E, 11/2011
2/2