

LAPspay Systems for Dogs and Cats

Laparoscopic Ovariectomy / Ovariohysterectomy
Single and Multi-Port



LAPspay Multi-Port

Standard Laparoscope



62046 AA **HOPKINS® Straight Forward Telescope 0°**, enlarged view, diameter 5 mm, length 29 cm, **autoclavable**, fiber optic light transmission incorporated, color code: green

TERNAMIAN ENDOTIP Cannulas



60120 MS TERNAMIAN **ENDOTIP Cannula**, size 6 mm, color code: black-white including:
Cannula, with thread, working length 6.5 cm
Silicone Leaflet Valve

60120 MSS **Same**, with insufflation stopcock


Standard Trocar and Cannula, with multifunctional valve



62160 WZ **Trocar**, size 6 mm, color code: black-white including:
Trocar, only, with pyramidal tip
Cannula, without valve, with insufflation stopcock, working length 8.5 cm
Multifunctional Valve


LAPspay Single Port

Operating Laparoscope



26034 AAV **HOPKINS® Operating Telescope 0°**, with parallel eyepiece, diameter 10 mm, length 23 cm, **autoclavable**, fiber optic light transmission incorporated, with 6 mm working **channel**, including adaptor with Silicone Leaflet Valve 26034 V

TERNAMIAN ENDOTIP Cannula, with multifunctional valve



60103 TMR TERNAMIAN **ENDOTIP** Cannula, with thread and rotatable insufflation stopcock, size 11 mm, working length 10.5 cm, color code: green
including:
Cannula
Multifunctional Valve
Seal Set 60/10
Mandrel, for opening automatic valves

Standard Trocar and Cannula



62103 FZ **Trocar**, size 11 mm, color code: green
including:
Trocar, only, with pyramidal tip
Cannula, without valve, with insufflation stopcock, working length 8.5 cm
Automatic Valve



62141 DB **Reducer**, 11/5 mm

Whether you prefer the Multi-Port or Single-Port technique, minimally invasive ovariectomy and ovariectomy using a KARL STORZ LAPspay surgery system offers the following advantages:

- Unsurpassed optics for a clear, bright, enlarged view
- Less pain for the patient since tissues are cauterized and transected with precision, rather than torn
- Complications, such as wound dehiscence, ovarian remnants and ligated ureters, are reduced to a minimum
- Reduced chance of infection due to minimal exposure of the abdomen
- Faster recovery due to small incision size
- Satisfied clients who demand the best for their pets
- Practice expansion through state-of-the-art surgical services
- The opportunity to improve your laparoscopy skills, which can be applied to numerous other minimal access surgical techniques



Instruments, 5 mm



66322 MG

CLICKLINE Grasping Forceps, “tiger jaws”, rotating, dismantling, insulated, with connector pin for unipolar coagulation, with LUER-Lock irrigation connector for cleaning, single action jaws, 2 x 4 teeth, for grasping and removal of solid organs, size 5 mm, length 36 cm including:

Plastic Handle, with MANHES style ratchet

Metal Outer Sheath, insulated

Forceps Insert

33422 MG

Same, length 43 cm

Additional Instruments



62175 H

TANKERSLEY Ovariectomy Hook, with heavy handle, for percutaneous use, working length 8.5 cm

62175 HK

Same, with small hook, for small dogs and cats



62120 J

VERESS Pneumoperitoneum Needle, with spring-loaded blunt stylet, LUER-Lock, diameter 2.1 mm, length 10 cm

Camera System



- 69 2360 01 **Veterinary Video Camera III**, color system PAL, power supply 100-240 VAC, 50/60 Hz including:
Veterinary Video Camera III Camera Control Unit, color systems **PAL/NTSC**, with integrated digital fiberscope filters
Mains Cord
2x **ACC Connecting Cable**, for controlling KARL STORZ units, length 180 cm
BNC Video Cable, length 200 cm
S-Video (Y/C) Connecting Cable, length 180 cm
C-MOUNT Lens, f = 38 mm, soakable
VET C-MOUNT Camera Head, PAL

- 69 2361 01 **Same**, color system NTSC

Light Source



- 69 1315 01 **Cold Light Fountain XENON NOVA® 175**, power supply 100-240 VAC, 50/60 Hz including:
Mains Cord
- 20 1320 26 **XENON Spare Lamp**, 175 W, 15 V
- 69495 NA **Fiber Optic Light Cable**, with straight connector, diameter 3.5 mm, length 230 cm (not illustrated)

For other camera system and light source options, see catalog: *Veterinary Endoscopy, Small Animals*

Insufflator



62432501 CO₂ ENDO-ARTHROFLATOR-VET,
power supply 100-240 VAC, 50/60 Hz
including:
Silicon Tubing Set, sterilizable
Universal Wrench
Gas Filter, sterile, for single use,
package of 10

The pressure hose and CO2 bottles must be ordered separately





THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: vet@karlstorz.com
www.karlstorz.com