




EndoCAMeleon[®]


Change your Viewing Angle, not your Endoscope



STORZ
KARL STORZ — ENDOSKOPE

EndoCAMEleon® ARTHRO HOPKINS® Telescope

ENDOCAMELEON® is a multifunctional endoscope featuring HOPKINS® technology that has been expanded to include the versatility of a variable direction of view. ENDOCAMELEON® allows the user to select the desired direction of view without the changeover of telescopes. This combination guarantees very good image quality and excellent visualization.




ENDOCAMELEON® – Imitating Nature

The eye of the chameleon is the result of millions of years of evolution.

Perfectly adapted to the environment, the chameleon can look in almost any direction.

This “optical” natural wonder provided the inspiration for the ENDOCAMELEON® design.



Practical – Lightweight – Ergonomic

At only **90 g**, ENDOCAMELEON® is extremely lightweight and with a diameter of **4 mm** is optimally adapted to the requirements of **arthroscopy**.



ENDOCAMELEON® – Forward Vision

Almost as versatile as the eye of a chameleon, the ENDOCAMELEON® lens can be adjusted in five steps from 15° to 90° and thus allows a comprehensive view of intra-articular zones during surgical procedures.



It is recommended to check the suitability of the product for the intended procedure prior to use.

75 Years

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

