

SMART SCOPE



KARL STORZ SMART SCOPE

Do you need documentation possibilities that are quick and easy to use at any time?

The new smartphone adaptor from KARL STORZ allows straightforward documentation of endoscopic images on a smartphone. Combined with KARL STORZ endoscopes, the SMART SCOPE ensures remarkable image quality.



The data transfer function on the smartphone enables images and videos to be easily transferred to the doctor's office computer system. The SMART SCOPE can be adapted to all standard eyepieces and is compatible with various smartphone models. This ease of use provides the flexibility required in any setting – whether in the doctor's office or for consultation purposes.



Plug & Play



Can be used anywhere and at any time



Recording and transmission



Enables users to query colleagues and share information



Compatible with standard eyepieces and various smartphones

It is recommended to check the suitability of the product for the intended procedure prior to use.



...evolution continues

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: info@ksea.com



Health
made in Germany

www.health-made-in-germany.com