

# KARL STORZ GI SILVER SCOPE® SERIES



## For Flexible Videoendoscopy

KARL STORZ SILVER SCOPE series provides the surgeon with the best possible image to use intraoperatively in the OR for identifying disease or confirming a secure anastomosis. General, Bariatric, Thoracic, and Colorectal surgeons are adopting the use of GI scopes in the OR to ensure “function follows form” and to ensure the best outcome for their patients.


Today, KARL STORZ once again sets new standards in imaging for minimally invasive surgery with its IMAGE1 S™ technology. The modular design of IMAGE1 S™ allows customization to individual needs and integrates today’s video endoscopes from KARL STORZ to create a new camera platform via IMAGE1 S X-LINK.

# KARL STORZ GI SILVER SCOPE® SERIES

## Video Enhancements

The SILVER SCOPE series provides a clear crisp image and with IMAGE1 S™ allows the physicians to use the Picture in Picture (PiP) feature to view laparoscopic and endoscopic views on the same monitor. It enables all personnel in the OR to view the same image, improving patient safety.

Product Specifications:	13820NKSK Slim Video Gastroscope	13821NKSK Video Gastroscope	13920NKSK Video Sigmoidoscope	13925NKSK Video Colonoscope
Working Length	110 cm	110 cm	80 cm	160 cm
Outer Diameter	5.9 mm	9.3 mm	11.2 mm	12.9 mm
Working Channel	2.0 mm	2.8 mm	3.4 mm	3.8 mm
Deflection Up/Down	210°/100°	210°/100°	210°/80°	180°/160°
Deflection Right/Left	120°	120°	160°	120°



**Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.**

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)



KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100  
Phone toll free: 800 421-0837 (US only)  
Fax: +1 424 218-8525  
Fax toll free: 800 321-1304 (US only)  
E-Mail: [communications@karlstorz.com](mailto:communications@karlstorz.com)