

NEW FLEXIBLE VIDEO-NEURO ENDOSCOPE



Improved Visualization

KARL STORZ is the only company offering an FDA-cleared Flexible Video Neuro-Endoscope and comprehensive hydrocephalus portfolio to help you increase efficiency and surgeon satisfaction in the OR.

Patient Safety — The Flexible Video Neuro-Endoscope, with its smaller outer diameter, is particularly useful for treating patients with a narrow foramen of Monro (adult and pediatric patients).

Patient Safety and Outcomes — Studies have found that flexible endoscopes can enhance the ease of Endoscopic Third Ventriculostomy with Endoscopic Tumor Biopsy (ETV/ETB) procedures.

STORZ
KARL STORZ—ENDOSKOPE

NEW FLEXIBLE VIDEO-NEURO ENDOSCOPE

Improved Visualization

Hydrocephalus is the #1 reason for brain surgery in pediatrics and in the ventricles of the adult population. Our Flexible Video-Neuro Endoscope enhances visualization and access during diagnostic and therapeutic use during cranial procedures. This makes it ideal for Hydrocephalus treatment such as ETV, ETV+CPC, tumor biopsy and more.

Specifications — Flexible Video-Neuro Endoscope

- Shaft outer diameter of 2.90 mm with elliptical endtip that measures 2.42 mm X 3.22 mm
- Working length 35 cm and working channel 1.2 mm I.D.
- Peel-away sheath size: 11 Fr.
- Variable deflection up to 270°
- Integrated LED light source, so no light cable needed
- Works with manual arm fiberscope clamping jaw, 28272UF and 1 mm flexible instruments

Sterilization

- Ethylene oxide (EtO) sterilization
- STERRAD NX and 100NX sterilization systems (containers KSZ-39406AS or 39406AS)*
- V-PRO maX Sterilization System (container KSZ-39406AS) or V-PRO Tray (VP0004) and Mat (VP0033)**
- STERIS SYSTEM 1E Liquid Chemical Sterilant Processing System



Screen door effect



Image through Flexible Video Neuro-Endoscope

Ordering Information	
11161VK	CMOS Video Neuroscope 8.5 Fr. x 350 mm

* STERRAD, NX and 100NX are registered trademarks of Advanced Sterilization Products, a Johnson & Johnson company.
** STERIS, SYSTEM 1E and V-PRO are registered trademarks of STERIS Corp.

It is recommended to check the suitability of the product for the intended procedure prior to use.



KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com