

C-MAC® S IMAGER WITH C-MAC® SYSTEM INTERFACE



Opening Up New Possibilities

Our new generation of lightweight yet robust C-MAC® S Imager gives you entirely new possibilities for setup and use. Incorporated with the C-MAC® system interface, you gain the choice of connecting the unit to the C-MAC® PM or C-MAC® Monitor via the C-MAC® Cable.

- The future-proof C-MAC® system interface offers a new standard in image quality with its high-resolution chip.
- A specially coated shaft allows smooth insertion into C-MAC® S blades.
- One-button operation for documentation is provided by using the “BlueButton” function.
- Switch to a new blade in seconds using the open-to-intubate feature when connected with the new C-MAC® Pocket Monitor.

C-MAC® S IMAGER WITH C-MAC® SYSTEM INTERFACE

Opening Up New Possibilities

Expanding possibilities in the prehospital market, the C-MAC® S Imager offers one size of reusable adult imager (8403 XSI) that fits 3 different sizes of single-use adult blades (MAC #3, #4 and D-BLADE). The imager is reusable, waterproof and fully immersible for cleaning and disinfection.

Ordering Information	
8403 XSI	C-MAC® S Imager , with C-MAC® system interface
8403 XDK	C-MAC® Pocket Monitor Set , unit with LCD monitor and power supply for all C-MAC® video laryngoscopes
8403 XDL	Charging Unit , for one rechargeable battery 8403 XDA for C-MAC® Pocket Monitor 8403 XD
8403 ZXK	C-MAC® Monitor for CMOS Endoscopes
8403 X	C-MAC® Connecting Cable
051113-10	C-MAC® S Video Laryngoscope MAC #3 , with MACINTOSH laryngoscope blade, size 3
051114-10	C-MAC® S Video Laryngoscope MAC #4 , with MACINTOSH laryngoscope blade, size 4
051116-10	C-MAC® S Video Laryngoscope D-BLADE



It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

© 2017 KARL STORZ Endoscopy-America, Inc.