

PROTECTION FROM BLOCKED VIEWS C-MAC® BOEDEKER SUCTION TUBES



When the camera lens of a video laryngoscope blade is obstructed by blood or secretions, you can no longer see what you are doing. You must remove the blade for cleaning and attempt another intubation. Or, you may be forced to use a direct laryngoscope. KARL STORZ is the only company offering C-MAC® BOEDEKER Suction Tubes.

- Channeled C-MAC® blades hold the suction tubes securely
- Suction starts at the distal tip — blades can be inserted slowly while suctioning to protect against obstructed views
- Blade tip can be seen and used as a guide while suctioning
- Frees up “one hand,” as opposed to using a Yankauer
- Can be used with the oxygen port for passive oxygenation or to blow away secretions



PROTECTION FROM BLOCKED VIEWS

C-MAC® BOEDEKER SUCTION TUBES

An unobstructed airway is crucial to patient survival during intubation. The ability to rapidly clear the airway to improve ventilation or facilitate breathing is a critical aspect of airway management. C-MAC® BOEDEKER Suction Tubes help protect against a video laryngoscope black-out. The disposable suction tubes are especially useful in the ED, EMS and prehospital settings as these clinicians deal primarily with trauma patients and unexpected difficult airways.

In addition, because our blades use traditional Macintosh and Miller designs, your 'Plan B' is already in place. In the event a lens obstruction causes a sudden blackout, you can still revert to your direct laryngoscopy skills. You can use the C-MAC® blade as a regular direct laryngoscope — without losing critical time cleaning the blade or locating a direct laryngoscope replacement.

Ordering Information:

8401YST-20	C-MAC® BOEDEKER Suction tubes (box of 20) <i>Specially designed for C-MAC® blades.</i> <i>For use with:</i> 8403AX — Mac #3 w/channel 8403BX — Mac #4 w/channel 8403HX — D-BLADE w/channel
------------	---



It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com