


C-MAC® STYLET FOR ET TUBE PLACEMENT



051116-10 C-MAC S
Single-Use D-BLADE

8403 HX C-MAC®
Video Laryngoscope
D-BLADE

8403ZXK C-MAC® Monitor

8401DS C-MAC® Stylet

The Reusable Stylet for the C-MAC® D-BLADE Video Laryngoscope

The C-MAC® Stylet enables optimum placement of endotracheal tubes (ET tubes) with the C-MAC® D-BLADE Video Laryngoscopes. The angle of the stylet is adapted to the blade shape of the D-BLADE. Furthermore, gliding properties are enhanced by the device's special surface structure. The new stylet can be used if the ET tube requires added stability. The C-MAC® Stylet is suitable for use with both D-BLADE and C-MAC® MACINTOSH blades.

STORZ
KARL STORZ—ENDOSKOPE

C-MAC® STYLET FOR ET TUBE PLACEMENT

The C-MAC® D-BLADE, together with the C-MAC® reusable stylet, offers a simple, rapid intubation solution for unanticipated difficult airways, particularly in patients with problem anatomies. The D-BLADE is available in reusable and single-use versions.

Special features include:

- Reusable stainless-steel stylet with atraumatic tip
- Total length, 31 cm; each package contains 10 pieces
- Distal region adapted to the blade shape of the C-MAC® reusable and single-use D-BLADE Video Laryngoscopes
- For use with ET tubes from 6.5 mm and larger
- Integrated tube holder for fixation

Ordering Information:	
8401 DS	C-MAC® Guide, guide rod made of stainless steel with atraumatic tip, distal region is adapted to the blade shape of the C-MAC® video laryngoscope D-BLADE, fixation of endotracheal tube with the integrated tube holder possible, package of 10, for use with C-MAC® video laryngoscope D-BLADE 8401 HX and C-MAC® S single-use D-BLADE 051116-10
8403 HX	C-MAC® Video Laryngoscope D-BLADE, for difficult intubation, with catheter introduction sizes 16-18 Fr., for use with C-MAC® connection cable 8403 X for C-MAC® Monitor 8403 Z XK or C-HUB® II 20290301, as well as C-MAC® PM 8403 XD
051116-10	C-MAC® S Video Laryngoscope D-BLADE, with laryngoscope blade, for difficult intubation, for single use, package of 10, for use with KARL STORZ C-MAC® Monitor 8403 Z XK and C-MAC® S Imager 8403 XS or 8403 XSI



KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34
78532 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com