


## KARL STORZ Single-use Stone Baskets



## The ideal stone basket for each requirement

The numerous indications and methods of application in stone therapy require stone baskets with various properties – a challenge that KARL STORZ meets with its new product line of nitinol stone baskets. The stone baskets match the sizes of the working channels and can adapt to difficult anatomies thanks to the highly flexible nitinol wires.



## Design variations



Straight without tip



Straight with tip



Helical with tip

## Nitinol stone baskets for...

### ... rigid cystoscopy

Order no.	Size	Form	Wires	Working length
27023 TH	5 Fr.	with tip, straight	4	90 cm

### ... flexible cystoscopy and percutaneous nephrolithotomy

Order no.	Size	Form	Wires	Working length
27023 KF	3 Fr.	without tip, straight	4	70 cm

### ... semirigid uretero-rensoscopy

Order no.	Size	Form	Wires	Working length
27023 TD	2.5 Fr.	with tip, helical	4	120 cm
27023 LD	2.5 Fr.	without tip, straight	4	120 cm
27023 TF	3 Fr.	with tip, helical	4	120 cm
27023 LF	3 Fr.	without tip, straight	4	120 cm

### ... flexible uretero-rensoscopy

Order no.	Size	Form	Wires	Working length
27023 LB	1.9 Fr.	without tip, straight	4	120 cm

It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.



# 75 Years

*Shaping the Future  
of Endoscopy with you*



**STORZ**  
KARL STORZ—ENDOSKOPE  
THE DIAMOND STANDARD

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

