



KARL STORZ Single-use Stone Baskets

The right basket for every indication
in endourological stone therapy



The ideal stone basket for each requirement

The numerous indications and methods of application in stone therapy require stone baskets with various properties – a challenge that KARL STORZ meets with its new product line of nitinol stone baskets. The highly flexible nitinol wires allow baskets to adapt very easily to difficult anatomies and are also gentle on the working channel of your instruments.




Design variations



Straight without tip



Straight with tip



Helical with tip

Nitinol stone baskets recommended for...

... rigid cystoscopy

Order no.	Size	Form	Wires	Working length
27023 TH	5 Fr.	with tip, straight	4	90 cm

... flexible cystoscopy and percutaneous nephrolithotomy

Order no.	Size	Form	Wires	Working length
27023 KF	3 Fr.	without tip, straight	4	70 cm

... semirigid uretero-rensoscopy

Order no.	Size	Form	Wires	Working length
27023 TD	2.5 Fr.	with tip, helical	4	120 cm
27023 LD	2.5 Fr.	without tip, straight	4	120 cm
27023 TF	3 Fr.	with tip, helical	4	120 cm
27023 LF	3 Fr.	without tip, straight	4	120 cm

... flexible uretero-rensoscopy

Order no.	Size	Form	Wires	Working length
27023 LB	1.9 Fr.	without tip, straight	4	120 cm

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com

