

# The HyDOME™ System

The Solution for Hygiene in ERCP



**STORZ**  
KARL STORZ — ENDOSKOPE

# HyDOME

With the introduction of the HYDOME™ system, KARL STORZ has reinforced hygiene standards in ERCP procedures, retained the performance capabilities of a standard system and provided physicians with an option to work with brand new sterile modules every time.

This system, which is unique in the market, features an elevator wire channel that can be brushed, and ensures that neither the elevator nor the elevator wire are potential sources of cross-contamination.



STERILE



- Sterile elevator module
- Secured positioning with a single click
- Enables brushing of the elevator wire channel



It is recommended to check the suitability of the product for the intended procedure prior to use.



...evolution continues

**STORZ**  
**KARL STORZ — ENDOSKOPE**

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)