

# LARYNGOBLOC

## SINGLE-USE LARYNGOSCOPES



### LED Laryngoscope for single use

LARYNGOBLOC is a single-use intubating laryngoscope with a one-piece-design to improve patient safety and eliminate reprocessing and sterilization.

- Extremely bright, 80k Lux LED white light optimally illuminates the oral cavity
- Reusable lithium battery insert reduces disposal of environmentally critical material
- Individually clean packed, single-use blades reduce risk of cross-contamination
- One piece design increases tensile strength, making LARYNGOBLOC more robust than conventional single-use laryngoscopes
- Polypropylene material eliminates potential for deformation and sensitivity to cold temperatures
- Recycling program reduces impact of medical waste on the environment

# LARYNGOBLOC

## SINGLE-USE LARYNGOSCOPES

LED Laryngoscope for single use

Ordering Information	
8512 AF	<b>MACINTOSH (Low Profile Blade)</b> , size 4, standard handle
8512 BF	<b>MACINTOSH (Low Profile Blade)</b> , size 3, standard handle
8512 C	<b>MACINTOSH</b> , size 2, standard handle
8512 ASF	<b>MACINTOSH (Low Profile Blade)</b> , size 4, short handle
8512 BSF	<b>MACINTOSH (Low Profile Blade)</b> , size 3, short handle
8512 CS	<b>MACINTOSH</b> , size 2, short handle
8530 B	<b>MILLER</b> , size 3, standard handle
8530 C	<b>MILLER</b> , size 2, standard handle
8530 DN	<b>MILLER</b> , size 1, slimline
8530 EN	<b>MILLER</b> , size 0, slimline
8533 F	<b>OXFORD</b> , size 0, standard handle
8533 FS	<b>OXFORD</b> , size 0, short handle
8547 HL	<b>LED Battery Insert</b> , 2 batteries, long
8547 HS	<b>LED Battery Insert</b> , 1 battery, short
8547 HL	<b>LED Battery Insert</b> , 2 batteries, slimline
121306 KS	<b>AA LR06 1.5v</b> , Qty. 2
8400 A	<b>Bag for Intubation Set</b> , splash-protected
KS-BLBR	<b>Recycling Service</b> , set of 10 recycling boxes, 12x12x20



**It is recommended to check the suitability of the product for the intended procedure prior to use.**

KARL STORZ SE & Co. KG  
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
 Postbox 230, 78503 Tuttlingen/Germany  
 Phone: +49 (0)7461 708-0  
 Fax: +49 (0)7461 708-105  
 E-Mail: info@karlstorz.com  
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.  
 2151 East Grand Avenue  
 El Segundo, CA 90245-5017, USA  
 Phone: +1 424 218-8100  
 Phone toll free: 800 421-0837 (US only)  
 Fax: +1 424 218-8525  
 Fax toll free: 800 321-1304 (US only)  
 E-Mail: communications@karlstorz.com