


# NEW! AUTOCLAVABLE SHUNTSCOPE



ShuntScope™

## Endoscopic Placement of Ventricular Catheters During Shunt Surgery

KARL STORZ is the only company offering a reusable SHUNTSCOPE, and a comprehensive portfolio for endoscopic treatment of hydrocephalus.

# NEW! AUTOCLAVABLE SHUNTSCOPE

## We Can Help You:

**Enhance Clinical Outcomes** — More than 30,000 shunt procedures are performed in the US each year.<sup>1</sup> A recent study found 33% of shunt failures within 90 days were due to potentially preventable causes — the most common being inaccurate proximal catheter placement.<sup>2</sup>

**Promote Cost Effectiveness** — Our reusable SHUNTSCOPE offers an economical solution to support endoscopic-assisted placement of a ventricular catheter, Endoscopic Third Ventriculostomy (ETV) and Endoscopic Biopsy.

**Use Clinical Advantages** — The SHUNTSCOPE and operating sheath, with its smaller outer diameter (2.6 mm), can be useful when treating patients with a narrow foramen of Monro.

Ordering Information	
28164SSA	SHUNTSCOPE 0° Semi-Rigid Telescope, 1 mm x 16 cm
11577KA-NSM	Examination Sheath, 1.3 mm x 16 cm
11577KE-NSM	Operating Sheath with 1.15 mm working channel, 1.3 mm x 2.6 mm oval, length 16 cm
28160TV	Flexible Ventriculostomy Forceps, double action outer serrated jaws, 1 mm x 30 cm
28160TJ	Flexible Serrated Grasping Forceps with double action jaws, 1 mm x 30 cm
28160ZJ	Flexible Biopsy Forceps with double action jaws, 1 mm x 30 cm
28160KS/4	Flexible Monopolar Spatula Electrode, 1 mm x 45 cm, single use, sterile, 4/pkg
28160KA/4	Flexible Monopolar Electrode, 1 mm x 53 cm single use, sterile, 4/pkg
26006M-NSM	Monopolar High Frequency Cord, for Valleylab Electrosurgical units
495NTA	Fiber Optic Light Cable, diameter 2.5 mm, length 230 cm



- PATWARDHAN RV, NANDA A. Implanted ventricular shunts in the United States: The billion-dollar-a-year cost of hydrocephalus treatment. *Neurosurgery*. 2005;56:139–4
- VENABLE GT, ROSSI NB, JONES GN, et al. The Preventable Shunt Revision Rate: a potential quality metric for pediatric shunt surgery. *Journal of Neurosurgery*. July 2016;18 (1):7–15.

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

KARL STORZ SE & Co. KG  
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
 Postbox 230, 78503 Tuttlingen/Germany  
 Phone: +49 (0)7461 708-0  
 Fax: +49 (0)7461 708-105  
 E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)

[www.karlstorz.com](http://www.karlstorz.com)



KARL STORZ Endoscopy-America, Inc.  
 2151 East Grand Avenue  
 El Segundo, CA 90245-5017, USA  
 Phone: +1 424 218-8100  
 Phone toll free: 800 421-0837 (US only)  
 Fax: +1 424 218-8525  
 Fax toll free: 800 321-1304 (US only)  
 E-Mail: [communications@karlstorz.com](mailto:communications@karlstorz.com)