

LARYNGOFORCE® II CLIP FORCEPS



45° Angled Jaws Enable Precise Titanium Clip Application

The KARL STORZ LARYNGOFORCE® II Clip Forceps, 8665 LL and 8665 RR, are modified versions of the LARYNGOFORCE® II Clip Forceps, 8665 L and 8665 R. These are used in laryngoscopy and tracheoscopy for the application of Ethicon LT200 Medium Titanium Clips

- The 45° angled jaws of the LARYNGOFORCE® II Clip Forceps can be visualized more clearly under the microscope than with the standard model
- Improved visualization of the jaws allows for easier, more precise application of the titanium clip

LARYNGOFORCE® II CLIP FORCEPS

Enables Precise Titanium Clip Application

Ordering Information	
8665 LL	LARYNGOFORCE® II Clip Forceps , jaws curved 45° to left, with cleaning connector, working length 22 cm, for use with titanium clip Ethicon LT200 Medium Titanium Clips
8665 RR	LARYNGOFORCE® II Clip Forceps , jaws curved 45° to right, with cleaning connector, working length 22 cm, for use with titanium clip Ethicon LT200 Medium Titanium Clips

**Contact your local KARL STORZ representative for more information
about these or other KARL STORZ products, or call 800-421-0837.**

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com