

NEW PYLOROMYOTOMY KNIFE



Introducing The New Pyloric Knife for Pediatric Surgery







The KARL STORZ instrument set for laparoscopic pyloromyotomy has been expanded to include a new pyloric knife that facilitates pediatric surgery. The knife includes an integrated ratchet mechanism that allows the knife insert to be locked in 8 positions and extended to the desired length in a controlled manner.

- Sterile single-use knife blades (package of 6 per box) provide optimal cutting
- Ergonomic and reusable handle
- Laser-etched millimeter markings on the knife blade ensure proper lengths
- Easy button release for safe assembly and disassembly of the knife blade

STORZ
KARL STORZ—ENDOSKOPE

NEW PYLOROMYOTOMY KNIFE

Introducing The New Pyloric Knife for Pediatric Surgery

New Pyloric Knife			
	26167 KS	Handle with outer sheath, size 3 mm, length 10 cm, for use with knife Insert	
	26167 KRS-S	Sterile knife insert, for single use, package of 6, for use with handle with outer sheath 26167 KS	
	26167 DRS-KT	Pyloric knife kit consisting of: <ul style="list-style-type: none"> • Handle with outer sheath • Sterile knife insert, for single-use, package of 6 	
3 mm diameter, 20 cm length, (2-Piece System)			
	30310 PS	CLICKline outer sheath with working insert, GEIGER pyloric grasper, with LUER-Lock connector for cleaning, double action jaws	
	30361 PS	CLICKline GEIGER pyloric grasper, rotating, for percutaneous use, consisting of: <ul style="list-style-type: none"> • CLICKline outer sheath with working insert, GEIGER pyloric grasper • Metal handle 	
	3 mm diameter 20 cm length (3-Piece System)	3.5 mm diameter 20 cm length (3-Piece System)	
	31110 PS		CLICKline forceps insert, percutaneous pyloric spreader
	31100 SM	31100 M	CLICKline metal outer sheath, with LUER-Lock connector for cleaning*
	31161 PSSM	31161 PSM	Complete instrument with CLICKline handle 33161**
	31110 VR		CLICKline Forceps Insert, for TAN pyloric spreader and grasping forceps, double action jaws, jaws with rounded tips, in- and outside serrated
	31100 SM	31100 M	CLICKline metal outer sheath, with LUER-Lock connector for cleaning*
	31161 VRSM	31161 VRM	Complete instrument with CLICKline handle 33161**
	31110 VT		CLICKline Forceps Insert, for TAN pyloric spreader and grasping forceps, double action jaws, in- and outside serrated dissecting and grasping forceps, jaws with rounded tips
	31100 SM	31100 M	CLICKline metal outer sheath, with LUER-Lock connector for cleaning*
	31161 VTSM	31161 VTS	Complete instrument with CLICKline handle 33161**

* Also suitable with KARL STORZ CLICKline outer sheath, insulated: 31100 S (3 mm diameter) and 31110 (3.5 mm diameter).

** Also suitable to all KARL STORZ CLICKline handles of your choice. See main LAPAROSCOPY catalogue for the large variety of CLICKline handles.

It is recommended to check the suitability of the product for the intended procedure prior to use.



KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 (0)7461 708-0
 Fax: +49 (0)7461 708-105
 E-Mail: info@karlstorz.com

www.karlstorz.com

KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com