

C-MAC® VIP POLE STAND VIDEO INTUBATION PLATFORM



Keep Vital Airway Equipment Organized and Accessible

An unanticipated difficult airway can strike at any time ... keeping you as prepared as possible for your “Plan B” is our mission. That is why we designed the Video Intubation Platform carts to keep your essential C-MAC® devices close at hand. The VIP pole stand provides immediate access to video laryngoscopes and flexible scopes on a slim and smooth rolling cart. Secure organization of your video laryngoscopes and flexible scopes is provided while also ensuring immediate access during emergent events.

STORZ
KARL STORZ – ENDOSKOPE


C-MAC® VIP POLE STAND VIDEO INTUBATION PLATFORM

Keep Vital Airway Equipment Organized and Accessible

- **Articulating Arm** fixed height, cable management, bar tilt up or down, pan left and right for customized positioning
- **Two Scope Holders** help ensure secure storage and protection of flexible scopes
- **Two Storage Bins** made of polycarbonate include cable management hooks on the back
- **Curved Handle** is covered in a medical-grade Santoprene* material, making it smooth to the touch for function and added comfort



Articulating Arm



Two Scope Holders



Power Supply Storage



Two Storage Bins



**C-MAC® Video
Intubation Platform**

Part Number: 9700VIP	
Dimensions	Cart: 23 in. length / 19 in. width / 55 in. height Bins: 12 in. length / 6.5 in. width / 5 in. height
Weight	24 lbs.

* Santoprene is a registered trademark of Exxon Mobil Corp

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2018 KARL STORZ Endoscopy-America, Inc.