

KARL STORZ C-MAC® BLADE CONTAINER SYSTEM



Protect Your C-MAC® Blade Investment

Our KSZ-39301CMAC Container System is specially designed to protect two C-MAC® blades during transport, reprocessing and storage.

- Houses two C-MAC® video laryngoscopes
- Save time and money over traditional blue-wrap trays
- Anodized to prevent pitting and corrosion

STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ C-MAC® BLADE CONTAINER SYSTEM

Protect Your C-MAC® Blade Investment

The KSZ-39301CMAC blade sterilization container system is designed to house your KARL STORZ C-MAC® video laryngoscopes through transport, sterilization and storage. It offers hassle-free, cost-effective sterilization for your valuable C-MAC® blade investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They help eliminate the need for blue sterilization wrap, while saving time and money.

Validated for STERRAD 50, 100S, 200, NX and 100NX systems, as well as for STERIS V-PRO 1 Plus, V-PRO 60 and V-PRO MAX systems.*

Accessories For Steam Sterilization	
KSZ-SCF01	Disposable paper filters for steam sterilization, 1000 per case. 100% cellulose, and water repellent, 7.5" round.
KSZ-SCL01	Steam Sterilization Dual Process Indicator ID Cards, 1000 per case.
KSZ-SCI001	Steam Sterilization Dual Process Indicator Strips, 1000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure with up to two 4" long indicator strips.
KSZ-SCS01B	Steam Sterilization Tamper-Evident Seals 1000 per pack.

STERRAD, NX and 100NX are registered trademarks of Advanced Sterilization Products, Div. of Ethicon US, LLC.

STERIS and V-PRO are registered trademarks of STERIS, Inc.



Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 (0)7461 708-0
 Fax: +49 (0)7461 708-105
 E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com