

KARL STORZ PCNL CONTAINER SYSTEM



Protect Your PCNL System Investment

Our KSZ-39312J Container System is specially designed to protect your Percutaneous Nephrolithotomy (PCNL) system during transport, reprocessing and storage.

- Validated and approved for steam sterilization
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion

STORZ
KARL STORZ—ENDOSKOPE


KARL STORZ PCNL CONTAINER SYSTEM

Protect Your PCNL System Investment

The KSZ-39312J Sterilization Container System is designed to house your KARL STORZ PCNL equipment through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable PCNL equipment investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They eliminate the need for blue sterilization wrap, saving time and money.

Ordering Information

KSZ-SCF01	Disposable paper filters for steam sterilization, 1000 per case. 100% cellulose, and water repellent, 7.5" round.
KSZ-SCL01	Steam Sterilization Dual Process Indicator ID Cards, 1000 per case.
KSZ-SCI001	Steam Sterilization Dual Process Indicator Strips, 1000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure with up to two 4" long indicator strips.
KSZ-SCS01B	Steam Sterilization Tamper-Evident Seals, 1000 per pack.



**Contact your local KARL STORZ representative for more information
about these or other KARL STORZ products, or call 800-421-0837.**

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com