

MINIATURE OPTICAL FORCEPS FOR PREMATURE & NEONATE PATIENTS











Because Your Smallest Patients Deserve The Best

Complete your foreign body removal set for your smallest patients with the miniature optical forceps. The miniature optical forceps can be used with the smallest pediatric bronchoscopes, size 2.5 & 3.

MINIATURE OPTICAL FORCEPS FOR PREMATURE & NEONATE PATIENTS

Continuously improving our products so you can better care for your patients

Ordering Information		
	10030 AA	Miniature semi-rigid telescope , 1.3 mm X 30.6 cm
	10374 HF	Miniature optical forceps with alligator jaws , for the controlled grasping of hard foreign bodies, with force-limited handle, incl. cleaning adaptor, for use with miniature straight forward telescope 10030 AA
	10374 L	Miniature optical forceps, spoon-shaped , round, for biopsy, incl. cleaning adaptor, for use with miniature straight forward telescope 10030 AA

For Use With		
	10339 F	Bronchoscope tube , size 2.5, O.D. 4.2 mm, I.D. 3.5 mm, length 18.5 cm
	10339 E	Bronchoscope tube , size 3, O.D. 5 mm, I.D. 4.3 mm, length 18.5 cm
	10339 D	Bronchoscope tube , size 3.5, O.D. 5.7 mm, I.D. 5 mm, length 18.5 cm
	10339 EE	Bronchoscope tube , size 3, O.D. 5 mm, I.D. 4.3 mm, length 26 cm
	10339 DD	Bronchoscope tube , size 3.5, O.D. 5.7 mm, I.D. 5 mm, length 26 cm

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 (0)7461 708-0
 Fax: +49 (0)7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com