

HD-VIEW DIGITAL CYSTOSCOPE



HD-VIEW

See the Difference

IMAGE1 S™ High-Definition Image Quality with Built-In LED Light Source – See more color, contrast, tissue differentiation and depth cues without the need for an additional light source

Atraumatic Tip – Contoured shape for smooth access and increased patient comfort


Optional Suction and Larger Working Channel – Enhanced flow for rapid bladder draining and a clear surgical field of view

Designed for Durability – Working channel and distal articulation cover are fortified to resist scrapes and punctures

HD-VIEW DIGITAL CYSTOSCOPE

Customize your HD-VIEW with detachable working channel port options.

Product Specifications	
Outside Diameter	16 Fr.
Distal Tip	10.9 Fr.
Working Channel	7 Fr.
Deflection	210°/140°
Direction/Angle of View	0°/100°
Light Source	Integrated LED
Working Length	37 cm

Ordering Information	
11272VHUK	HD-VIEW Video Cystoscope (US Deflection)
11272VHK	HD-VIEW Video Cystoscope (Reverse Deflection)
11014LU	Double LUER Port
11014L	LUER Port with Stopcocks
11301P	Suction Port Seal
11301CE1/20	Suction Valve, Sterile (Pack of 20)
	NEW IMAGE1 S™ Office/Clinic Video Solution (IMAGE1S-X-KT) and Sony 21" HD Monitor with Table Top Stand (VTG0011)



HD-VIEW:
With Suction



HD-VIEW:
Without Suction

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE
THE DIAMOND STANDARD

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2019 KARL STORZ Endoscopy-America, Inc.