



ENDOCAMELEON® HOPKINS® TELESCOPE






Change your Viewing Angle, not your Endoscope

The KARL STORZ ENDOCAMELEON® HOPKINS® telescope with rod-lens system and a rigid 10 mm sheath offers a variable direction of view, as well as enhanced illumination. The ENDOCAMELEON® telescope offers a number of key capabilities and functions.

- Variable direction of view, 0° to 90°
- Adjusting knob for selecting the desired direction of view
- Paired with our new 4K camera allows for optimal 4K viewing
- Suitable for use in anatomically narrow working spaces

ENDOCAMELEON® HOPKINS® TELESCOPE

Change your Viewing Angle, not your Endoscope

Ordering Information		
	26003 EC	ENDOCAMELEON® HOPKINS® Telescope , diameter 10 mm, length 31 cm, autoclavable, variable direction of view 0° – 90°, with adjusting knob with fin for selecting the direction of view, fiber optic light transmission incorporated, color code: gold
	495 NCSC	Fiber Optic Light Cable , with straight connector, extremely heat-resistant, enhanced light transmission, with safety lock, diameter 4.8 mm, length 250 cm
	KSZ-39301C1S	No-Wrap Sterilization Container , holds one 34 cm rigid telescope up to 10 mm in diameter plus light cable, external dimensions (w x d x h) 18.5" x 6" x 3"

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 (0)7461 708-0
Fax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com