

# KARL STORZ FLEXIBLE GI SILVER SCOPE®



## Think Flexible — Be Flexible

The KARL STORZ Flexible GI SILVER SCOPE® represents a remarkably adaptable addition to your IMAGE1 S™ architecture.



- Optimize the performance of your GI Silver scope by ensuring availability of accessory components.
- Choose accessories for your GI scopes, Xenon cold light sources and reprocessing

**STORZ**  
KARL STORZ — ENDOSKOPE

# KARL STORZ FLEXIBLE GI SILVER SCOPE®

## Think Flexible – Be Flexible

Ancillary support equipment for new generation GI scopes and proprietary light source.  
The following accessories are necessary for appropriate use of the GI Videoscope System:

Accessories for Scopes	
 13991DA Rubber Lip Valve	 13991SSV Air/Water Valve
 13991SAV Suction Valve	 11025XE ETO Sealing Cap

Xenon Cold Light Source 20132601-1	
 13992BS Water Bottle	 20100031 Water Bottle Holder
 20132526 Spare Light Bulb	 13991SW Rinsing Connector

Reprocessing	
 13991CD Channel Separator	 13991CR Air/Water Channel Cleaning Adaptor
 13991CS Suction Channel Cleaning Adaptor	 13991SS Working Channel Cleaning Adaptor
 13991SWJ Water Jet Channel Cleaning Adaptor	 9990/50 Short Cleaning Brush
 110910-50 Long Cleaning Brush	 110917-50 Long Cleaning Brush for 13920NKSK & 13925NKSK Only
 27652/3 Flat Cleaning Brush	 13991SRV Cleaning Valve
	 13242LX Leakage Tester

GI Scopes	Product #	Distal Tip Outer Diameter (mm)	Instrument Channel Diameter (mm)	Working Length (cm)	Up (degrees)	Down (degrees)	Left (degrees)	Right (degrees)
Gastroscope	13820NKSK	5.9	2.0	110	210°	100°	120°	120°
Gastroscope	13821NKSK	9.3	2.8	110	210°	100°	120°	120°
Sigmoidoscope	13920NKSK	11.2	3.4	80	210°	180°	160°	160°
Colonoscope	13925NKSK	12.9	3.8	160	180°	180°	160°	160°

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105  
E-Mail: info@karlstorz.com  
  
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100  
Phone toll free: 800 421-0837 (US only)  
Fax: +1 424 218-8525  
Fax toll free: 800 321-1304 (US only)  
E-Mail: communications@karlstorz.com