

# The End of the Stone Age – CALCULASE® III




**STORZ**  
KARL STORZ – ENDOSKOPE

# CALCULASE® III – The End of the Stone Age

## 35 W Holmium Laser

The new 35 Watt desktop laser CALCULASE® III from KARL STORZ offers the user an enhanced performance and optimal user comfort with an excellent price-performance ratio\*.



### Fibers

- RFID fiber identification
- Backward compatibility of fibers to CALCULASE® II → delivers full performance capacity with CALCULASE® II, even with the thinnest fibers
- Single-use and reusable fibers (230 µm, 365 µm, 600 µm)



## Power

- Rating: Max. 35 W
- Energy: Max. 4 Joules
- Frequency: Max. 30 Hertz



## Fields of application

CALCULASE® III is versatile in its use. This laser system can be used for the lithotripsy of calculi and the treatment of tissue.



## Other Features

- Optimized pulse modes
  - Short: For stone fragmentation
  - Long: For dusting
  - Burst: For reducing the repulsion of stone fragments
- 10" HD touch screen
- Predefined procedures
- Compact desktop design

\* Comparison with predecessor  
CALCULASE® II

It is recommended to check the suitability of the product for the intended procedure prior to use. Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

The logo features the number '75' in a large, bold, blue sans-serif font. The '5' has a white circle in its center. Below the '75' is the word 'Years' in a blue, cursive script font.

*Shaping the Future  
of Endoscopy with you*

**STORZ**  
**KARL STORZ—ENDOSKOPE**

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

