

Silver is the New Black

LIVE – the new **LED** Integrated **V**ideo **E**ndoscopes
for bronchoscopy



LIVE – LED Integrated Video Endoscopes

Take a look into the future with the right technology:

A modern way of working with CMOS and HD without the need for a separate light source



Easy-to-use, flat suction valve, safely secured with bayonet coupling

Camera buttons individually programmable

Compatible with IMAGE1 S™ and TELE PACK+

one step deflection control



Ergonomic handle in elegant silver design

Easier maneuverability in the bronchial system due to a soft sheath

6.5 mm

5.5 mm

It is recommended to check the suitability of the product for the intended procedure prior to use.
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

More than
75
Years

*Shaping the Future
of Endoscopy with you*