

# Arthroscopic Vaporization Electrodes

Tissue removal through cutting or coagulating during arthroscopic procedures



# Arthroscopic Vaporization Electrodes

Fitting your requirements in arthroscopic procedures



- + Different electrodes with and without suction
- + 90°, 70° and 45° angled tips
- + Sterile electrodes for single use in packaging of 10 pcs.
- + Reusable handle with bipolar cable
- + To be used with KARL STORZ bipolar high-end generator AUTOCON® III 300
- + One extra long electrode 180 mm
- + Two electrodes with different hook sizes
- + Foot switch operation

Standard length:



Bipolar arthroscopic electrode (shaft Ø 5 mm) with suction, 90°, length 150 mm – single use, sterile.  
Cutting and coagulation mode, for tissue removal in the shoulder and knee.

---

Bipolar arthroscopic electrode (shaft Ø 5 mm) with suction, 45°, length 150 mm – single use, sterile.  
Cutting and coagulation mode, for tissue removal, ideal in the difficult-to-reach areas of the knee.

---

Bipolar arthroscopic electrode (shaft Ø 5 mm, double curved) with suction, 45°, length 150 mm – single use, sterile.  
Cutting and coagulation mode, for precise soft tissue removal in small and difficult-to-reach areas of the hip and knee.

---

Bipolar arthroscopic electrode (shaft Ø 5 mm, double curved) with suction, 45°, length 180 mm – single use, sterile.  
Cutting and coagulation mode, for precise soft tissue removal in difficult-to-reach areas of the hip.

---

Bipolar arthroscopic electrode (shaft Ø 4 mm), 70°, length 150 mm – single use, sterile.  
Cutting and coagulation mode, for the knee and shoulder.

---

Bipolar arthroscopic hook electrode (shaft Ø 3.3 mm), 90°, length 150 mm – single use, sterile.  
For tissue cutting only, with a short tip.


---

Bipolar arthroscopic hook electrode (shaft Ø 4 mm), 90°, length 150 mm – single use, sterile.  
For tissue cutting only, reinforced, with a long tip.

---

Bipolar arthroscopic electrode (shaft Ø 5 mm), 90°, length 150 mm – single use, sterile.  
Cutting and coagulation mode in the shoulder and knee.

---



It is recommended to check the suitability of the product for the intended procedure prior to use.  
Please note that the described products in this medium may not be available yet in all countries due to different regulatory requirements.

More than  
**75**  
Years

*Shaping the Future  
of Endoscopy with you*

**STORZ**  
KARL STORZ — ENDOSKOPE

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)

