

KARL STORZ RIGID TELESCOPE CONTAINER SYSTEMS



Protect Your Rigid Telescope Investment

Our Rigid Telescope Container Systems are specially designed to protect your rigid scope inventory during transport, reprocessing, and storage.

- Validated and approved for Steam, STERRAD® 100S®, NX® and 100NX®, as well as STERIS® AMSCO® V-PRO® sterilization systems*
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion

* STERRAD®, NX® and 100NX® are registered trademarks of Advanced Sterilization Products.
STERIS®, AMSCO® and V-PRO® are registered trademarks of STERIS Corp.

KARL STORZ RIGID TELESCOPE CONTAINER SYSTEMS

Protect Your Rigid Telescope Investment

The KARL STORZ Rigid Telescope Sterilization Container Systems are designed to house your KARL STORZ rigid telescopes through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable rigid telescope investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, while saving time and money.

HOPKINS® Rigid Telescope Sterilization Container System		
	KSZ-39301AS	4 mm Telescope Set , 2 scopes up to 4 mm diameter, 18 cm working length. With holder for light post adaptor and silicone scope holders. Includes basket, sealed container base and lid. 12" x 6" x 2" (30.5 x 15.24 x 5.08 cm).
	KSZ-39311AS	4 mm Telescope and Light Cable Set , 4 scopes up to 4 mm diameter, and 18 cm working length. With holder for light post adaptor, light cable and silicone scope holders. Includes basket, sealed container base and lid. 18" x 11.16" x 4.32" (45 x 28.35 x 40.64 cm).
	KSZ-39301BS	5 mm Telescope Set , 2 scopes up to 5 mm in diameter, 31 cm working length. With holder for light post adaptor and silicone scope holders. Includes basket, sealed container base and lid. 18.5" x 6" x 3" (47 x 15.24 x 7.62 cm).
	KSZ-SCS01B	10 mm Telescope Set , 2 scopes up to 10 mm diameter, 31 cm working length. With holder for light post adaptor and silicone scope holders. Includes basket, sealed container base and lid. 18.5" x 6" x 3" (47 x 15.24 x 7.62 cm).
	KSZ-39301C1S	10 mm Telescope and Light Cable Set , 1 scope up to 10 mm diameter, 31 cm working length. With holder for light cable and silicone scope holder. Includes basket, sealed container base and lid. 18.5" x 6" x 3" (48.26 x 15.88 x 9.02 cm).
	KSZ-39314FS	Four Telescope and Light Cable Set , 4 scopes up to 10 mm diameter, 31 cm working length. With additional light cable tray. Includes basket, sealed container base and lid. 23.5" x 11" x 4" (59.7 x 28 x 10 cm).

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 (0)7461 708-0
 Fax: +49 (0)7461 708-105
 E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2019 KARL STORZ Endoscopy-America, Inc.