


KARL STORZ FLEXIBLE SCOPE CONTAINER SYSTEMS



Protect Your Flexible Scope Investment

Our Flexible Scope Container Systems are specially designed to protect your flexible scopes during transport, reprocessing and storage.

- Validated and approved for STERRAD® NX® “Advanced” Cycle and 100NX® “Flex” Cycle, as well as STERIS® AMSCO® V-PRO® maX “Flexible” Cycle sterilization systems*
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion




* STERRAD®, NX® and 100NX® are registered trademarks of Advanced Sterilization Products.
STERIS®, AMSCO® and V-PRO® are registered trademarks of STERIS Corp.

KARL STORZ FLEXIBLE SCOPE CONTAINER SYSTEMS

Protect Your Flexible Scope Investment

The Flexible Scope Sterilization Container Systems are designed to house your KARL STORZ flexible scopes through transport, sterilization and storage. It offers cost-effective, hassle-free sterilization for your valuable flexible scope investment. KARL STORZ sterilization systems are made from aircraft-grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, saving time and money.

HOPKINS® Rigid Telescope Sterilization Container System	
KSZ-39402AS	Flexible Fiberscope Container , Includes base, basket and lid. 23.5" x 11" x 4" (59.7 x 28 x 10 cm)
KSZ-39403AS	Flexible CCD Video Endoscope Container , Includes base, basket and lid. 23.5" x 11" x 4" (59.7 x 28 x 10 cm)
KSZ-39406AS	Flexible CMOS Video Endoscope Container , Includes base, basket and lid. 23.5" x 11" x 4" (59.7 x 28 x 10 cm)

Accessories for Vaporized H ₂ O ₂ /STERRAD® Sterilization		
	KSZ-SCFO2	Disposable Polypropylene Filters , round, SteriTite® filter for vaporized H ₂ O ₂ /STERRAD® sterilization. 1000 per package. 7.5" diameter
	KSZ-SCLH2023	Large Load Cards , process indicator cards for containers, low temperature, SteriTite® cards for vaporized H ₂ O ₂ /STERRAD® sterilization. 1000 per package.
	KSZ-SCS01W	Tamper-Evident Seals , low temperature, for vaporized H ₂ O ₂ /STERRAD® sterilization. 1000 per package.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com