


# POLYMER RIGID SCOPE BASKETS



## Hassle-Free Telescope Protection

The KARL STORZ Polymer Rigid Scope Baskets are designed to minimize telescope breakage by decreasing the amount of manual handling during reprocessing. Combining the unique advantages of conventional wire trays with those of high-performance, lightweight polymer. The result is an innovative solution that can reduce the cost of telescope ownership.





- Reduce overall rigid scope repairs
- Optimized for machine processing
- Validated for STERRAD® sterilization methods\*

\* STERRAD® is a registered trademarks of Advanced Sterilization Products.

# POLYMER RIGID SCOPE BASKETS

## Telescope Protection For Cleaning, Sterilization & Storage

To protect the investment in rigid telescopes, the use of scope baskets allows telescopes to be processed through automatic washers, limiting manual exposure – a significant cause of damage. The Polymer Rigid Scope Baskets also protect delicate rigid telescopes from impact and vibration during transportation and storage.

HOPKINS® Rigid Telescope Sterilization Container System		
	39305 C1S	<b>Polymer Rigid Scope Basket</b> , for cleaning, sterilization, and storage of one rigid telescope up to 5 mm diameter and 34 cm working length. With holder for light post adaptor, silicone telescope holders and lid. 8.5" x 9" x 2" (47 x 10.5 x 5 cm)
	39305 C2S	<b>Polymer Rigid Scope Basket</b> , for cleaning, sterilization, and storage of two rigid telescopes up to 5 mm diameter and 34 cm working length. With holder for light post adaptor, silicone telescope holders and lid. 8.5" x 9" x 2" (47 x 10.5 x 5 cm)
	39305 L1S	<b>Polymer Rigid Scope Basket</b> , for cleaning, sterilization, and storage of one rigid telescope up to 10 mm diameter and 34 cm working length. With holder for light post adaptor, silicone telescope holders and lid. 8.5" x 9" x 2" (47 x 10.5 x 5 cm)
	39305 L2S	<b>Polymer Rigid Scope Basket</b> , for cleaning, sterilization, and storage of two rigid telescopes up to 10mm diameter and 34cm working length. With holder for light post adaptor, silicone telescope holders and lid. 8.5" x 9" x 2" (47 x 10.5 x 5 cm)

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG  
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
 Postbox 230, 78503 Tuttlingen/Germany  
 Phone: +49 7461 708-0  
 Fax: +49 7461 708-105  
 E-Mail: info@karlstorz.com  
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.  
 2151 East Grand Avenue  
 El Segundo, CA 90245-5017, USA  
 Phone: +1 424 218-8100  
 Phone toll free: 800 421-0837 (US only)  
 Fax: +1 424 218-8525  
 Fax toll free: 800 321-1304 (US only)  
 E-Mail: communications@karlstorz.com