

ENDOMAT® SELECT HYSTEROSCOPIC FLUID MANAGEMENT PUMP



The Choice is Clear

The ENDOMAT® SELECT is the latest innovation in fluid pumps, providing a smooth and consistent performance. Designed for optimal irrigation of fluids during hysteroscopic procedures, the compact single-roller pump can be conveniently placed in the doctor's office or clinic. Features include:

- User-friendly interface – easy to follow touch-screen operation
- Efficient procedural flow – fast and accurate pressure control
- Streamlined workflow – suitable for use in the doctor's office or clinic
- Optimal visualization – provides clear visualization during a procedure

ENDOMAT® SELECT HYSTEROSCOPIC FLUID MANAGEMENT PUMP

The Choice is Clear

Designed for use in hospitals, clinics, and physician practices, the ENDOMAT® SELECT from KARL STORZ is the optimal pump for all hysteroscopic needs.

Ordering Information	
UP210-GYN	Endomat® Select Pump
Software	
UP602	Hysteroscopy Software, license (allows selection of "HYS" procedure)
UP610*	Advanced Software, license (extends functions of installed software packages)
Accessories	
031523-10-GYN	Tubing Set, irrigation, pressure-controlled, sterile, for single use, package of 10
Optional Accessories	
030247-10	Outflow Y-Tubing, sterile, for single use, package of 10
26163 FD	Kimberly Clark Collection Drape, sterile, for single use, 1 pk
26018 FM-5	Fluid Management Pack, package of 5



* UP610 Advanced Software

- Adjust flow: 200/400/600 ml/min
- Adjust pump height according to patient position
- Customize key tone volume
- Modify brightness of touchscreen
- SCB integration

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com