

Quality You Trust, Available In Seconds

When securing the airway, time is critical



STORZ
KARL STORZ — ENDOSKOPE

C-MAC®

Your premium system
for single use



C-MAC® S – video laryngoscopes for single use

- Indirect and direct laryngoscopy – standard MACINTOSH blades in sizes 3 and 4
- Pediatric intubation – standard MILLER blades in sizes 0 and 1
- For the particularly difficult airway – hyperangulated D-BLADE
- Robust design – withstands even strong forces
- Blade and handle of C-MAC® S single-use video laryngoscopes form one continuous piece

It is recommended to check the suitability of the product for the intended procedure prior to use.

C-MAC® S – streamlined and improved

The current scientific recommendation “Effective use of video laryngoscopy in emergencies¹” recommends the use of video laryngoscopes with MACINTOSH blades since this blade shape also permits direct laryngoscopy if necessary.

Due to considerable material reduction, the single-use video laryngoscope blades of the C-MAC® system now feature similar dimensions as standard MACINTOSH blades and the proven reusable C-MAC® video laryngoscope blades.



The slimmer design improves the view during video laryngoscopy and direct laryngoscopy.

The video laryngoscopes are also suitable for use in training and continued education since the instructor can watch the intubation on the monitor.



¹Hossfeld B et al. Effective use of video laryngoscopy in emergencies. Notarzt 2017; 33: 242-249

For indirect & direct laryngoscopy



OTI® – Open To Intubate

Flipping the OTI® display automatically switches the device on and off. The more compact design facilitates stowing of the C-MAC® S video laryngoscope blade with the attached and folded C-MAC® Pocket Monitor.

OTI Open To Intubate®

The logo features the number '75' in a large, bold, blue sans-serif font. The '5' has a white circle in its center, resembling an endoscopic lens. Below the '75' is the word 'Years' in a blue, cursive script font.

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ — ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com