

KARL STORZ SEMI-RIGID URETEROSCOPE CONTAINER SYSTEMS



Protect your Semi-Rigid Ureteroscope Investment

Our KSZ-39301SR and KSZ-39301SP Container Systems are specially designed to protect your semi-rigid ureteroscopes during transport, reprocessing, and storage.



- Designed for one pediatric semi-rigid ureteroscope and small accessories
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion





STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ SEMI-RIGID URETEROSCOPE CONTAINER SYSTEMS

Cost-effective, hassle free sterilization

The Semi-Rigid Ureteroscope sterilization system are designed to house your KARL STORZ semi-rigid ureteroscopes through transport, sterilization, and storage. They offer cost-effective sterilization for your valuable semi-rigid ureteroscope investment. KARL STORZ sterilization systems are made from aircraft grade anodized aluminum to prevent pitting and corrosion. They eliminate the need for blue sterilization wrap, saving time and money.

Part Number		Description
KSZ-39301SR		Container System, Standard Semi-Rigids (27" x 11" x 6")
KSZ-39301SP		Container System, Pediatric Semi-Rigids (18.1" x 11" x 5")

Accessories for steam sterilization		Description
KSZ-SCF01		Disposable Paper Filters for Steam Sterilization , 1000 per case. 100% Cellulose and water repellent, 7.5" round.
KSZ-SCL01		Steam Sterilization Dual Process Indicator ID Cards , 1000 per case.
KSZ-SCI001		Steam Sterilization Dual Process Indicator Strips , 1000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure up to two 4" long indicator strips.
KSZ-SCS01B		Steam Sterilization Tamper Evident Seals , 1000 per pack.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com