

KARL STORZ SIALENDOSCOPY INSTRUMENT SET TRAY



Convenient Transportation, Sterilization and Storage


KARL STORZ now offers a sterilization tray system for the packaging and protection of the Sialendoscopy Instrument Set for sterilization, transportation, storage and aseptic presentation. Used with steam sterilization systems.

- Transparent lid facilitates quick identification and aseptic presentation of the procedural set, without the need for excessive handling.
- Insert tray provides additional level for storage protection for up to 3 telescopes
- Perforated base with silicone inserts ensures superior protection for individual probes, dilators, and other delicate components

KARL STORZ SIALENDOSCOPY INSTRUMENT SET TRAY

Convenient Transportation, Sterilization and Storage

Providing a secure, protective environment for the reprocessing of delicate instrumentation is critical to safeguarding your technology investments. The KARL STORZ Sialendoscopy Set Tray facilitates reliable sterilization, safe transportation, and secure storage of your Sialendoscopy telescopes, instruments and accessories.

Ordering Information		
	39312A-SIA	Sialendoscopy Instrument Set Tray , for sterilization and storage, perforated with transparent lid, external dimensions: 21.55" x 10.18" x 4.37" (548 x 259 x 110 mm), for the use with 3 sialendoscopes, instruments, and accessories.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com