

TELE PACK+

Compact Endoscopy

DISPLAY +

LIGHT +

CAMERA +

DOCU +

NETWORK +

STORZ
KARL STORZ ENDOSKOP

tele pack +
TP 101

STORZ
KARL STORZ — ENDOSKOPE

TELE PACK+

Compact Endoscopy

TELE PACK+ is an ALL-IN-ONE system that combines a monitor, LED light source, FULL HD camera control unit and documentation with integrated network function in a portable and compact unit. A very strong feature is its compatibility with a wide range of rigid, flexible and single-use endoscopes. TELE PACK+ is suitable for universal use in doctors' offices, day clinics, emergency rooms, intensive care units and outpatient settings.

DISPLAY +

- 18.5" touch screen monitor +
- On-screen keyboard +
- FULL HD +
- Display with no edges – ideal for cleaning +
- Straightforward handling thanks to intuitive menu navigation and self-explanatory user interface +

LIGHT +

- Integrated LED light source +
- Bright and uniform illumination +
- Automatic light source control +





CAMERA +

- + Can be connected to rigid as well as flexible and single-use endoscopes
- + Compatible with X-LINE and C-LINE endoscopes

DOCU +

- + Archiving of patient data on internal memory, on USB memory devices or via the network integration function
- + Encrypted patient data transfer
- + Patient data management
- + 50 GB internal memory

NETWORK +

- + Network function in File Share mode or Worklist mode
- + Export and import of patient data in/out of HIS/PACS

More than
75
Years

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ—ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

