

ENDOMAT® SELECT FLUID MANAGEMENT



The Choice is Clear

ENDOMAT® SELECT is the latest innovation in fluid management, providing smooth, reliable, and consistent performance. Designed for optimal irrigation or suction during urology procedures, this compact pump can be conveniently placed in the doctor's office or clinic. Features include:

- **Easy Set Up** with simple system interface and tubing set recognition
- Designed for exceptional **Clarity in Clinical Visualization**
- **Upgradable** software designed for multi-department use, upgradable via USB drive

STORZ
KARL STORZ – ENDOSKOPE

ENDOMAT® SELECT FLUID MANAGEMENT

Product Specifications	
Product Dimensions	305 mm x 110 mm x 260 mm (12.01" x 4.33" x 10.24")
Inflow/Outflow	Both
Flow Rate	Suction: 100 - 1000 ml/min Irrigation: 200 - 400 ml/min (up to 600 ml/min with Advanced Software)
Pressure	Max 150 mmHg



Ordering Information	
UP210	Endomat Select
UP604	Urology Software
030647-10	Direct Suction Tubing Set, sterile, pack of 10
031523-10	Irrigation Tubing Set, Pressure Controlled, sterile, pack of 10

Optional Advanced Accessories	
UP610	<p>Advanced Software:</p> <ul style="list-style-type: none"> • Adjustment of inflow to 600 ml/min • Modify duration of boost of pressure • Adjustment of software according to pump height • Customize loudness of the button tones • Fine-tune brightness of touchscreen • SCB integration
20014130	Single Pedal Footswitch

Urologic Procedures

Suction RES	Irrigation URS
Irrigation CYST	Irrigation PCN

Tubing Sets

<p>Direct Suction 030647-10</p> 	Designed for thorough removal of tissue chips during resection procedures.
<p>Irrigation 031523-10</p> 	Ideal for clearing the field of view for stone management and cystoscopy.

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com