

NEXT-LEVEL CLEAN – SILCUT® PRO AND SILGRASP® PRO HAND INSTRUMENTS



Instruments Designed With Cleaning In Mind

Hospital readmissions due to orthopedic infections are a real concern for surgeons and administrators alike. The financial costs associated with surgical-site infections can be significant—but the cost to a hospital's reputation can be even more profound. Proper cleaning and disinfection should always be the first line of defense.

If you're ready to take your instrument cleaning to the next level, KARL STORZ has got you covered! Our comprehensive line of arthroscopy instruments are specifically designed to provide the highest degree of simplicity—and certainty—in cleaning, all while maximizing surgical performance.

Our SILCUT® PRO punchers and SILGRASP® PRO graspers are a prime example, designed with an easy-to-clean inner shaft.

www.karlstorz.com

STORZ
KARL STORZ – ENDOSKOPE

NEXT-LEVEL CLEAN – SILCUT® PRO AND SILGRASP® PRO HAND INSTRUMENTS

Instruments Designed With Cleaning In Mind

The SILCUT® PRO and SILGRASP® PRO are the only hand instruments on the market designed with CLEANFIT® technology. A luer port near the handle makes it easy to flush out the inner channel: simply uncap the port, attach a syringe or irrigation tube, and send a stream of cleaning solution through the instrument. This helps prevent debris and residual biomaterial from remaining in the shaft and jaws. It's hard to believe that no one else offers this critical, common-sense feature!

The SILCUT® PRO and SILGRASP® PRO set new standards in arthroscopic resection and grasping. A wide variety of jaw configurations are available to suit the full range of surgical needs. Features include:

- Precise, sharp cutting
- Optimal balance between cutting force and cutting sensation
- Extremely flat instrument jaws
- Different edges for different tissue types
- Confident cleaning



Ordering Information	
28571AB	SILCUT® PRO Punch w/cleaning port, straight jaws
28571BB	SILCUT® PRO Punch w/cleaning port, jaws curved 15° upwards
28571BBL	SILCUT® PRO Upbiter w/cleaning port, jaws curved 15° upwards, sheath curved 15° left and 10° upwards
28571BBR	SILCUT® PRO Upbiter w/cleaning port, jaws curved 15° upwards, sheath curved 15° right and 10° upwards
28571AG	SILCUT® PRO Alligator Grasping Forceps w/cleaning port, serrated, straight jaws

For more information on the SILCUT® PRO, SILGRASP® PRO, and related instruments, please contact us at **1-800-421-0837**.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com