

TELE PACK +

Compact Endoscopy

DISPLAY +

LIGHT +

CAMERA +

DOCU +

NETWORK +

STORZ

tele pack +
TF 101

STORZ
KARL STORZ — ENDOSKOPE

TELE PACK +

Compact Endoscopy

TELE PACK+ is an ALL-IN-ONE system that combines a monitor, LED light source, FULL HD camera control unit and documentation with integrated network function in a portable and compact unit. A very strong feature is its compatibility with a wide range of rigid, flexible and single-use endoscopes. TELE PACK+ is suitable for universal use in doctors' offices, day clinics, emergency rooms, intensive care units and outpatient settings.

DISPLAY +

- 18.5" touch screen monitor +
- On-screen keyboard +
- FULL HD +
- Display with no edges – ideal for cleaning +
- Straightforward handling thanks to intuitive menu navigation and self-explanatory user interface +

LIGHT +

- Integrated LED light source +
- Bright and uniform illumination +
- Automatic light source control +





CAMERA +

- + Can be connected to rigid as well as flexible and single-use endoscopes
- + Compatible with X-LINE and C-LINE endoscopes

DOCU +

- + Archiving of patient data on internal memory, on USB memory devices or via the network integration function
- + Encrypted patient data transfer
- + Patient data management
- + 50 GB internal memory

NETWORK + *

- + Network function in File Share mode or Worklist mode
- + Export and import of patient data in/out of HIS/PACS

* Networking functionalities in the U.S. will require a future software release and the StreamConnect® media management application.

It is recommended to check the suitability of the product for the intended procedure prior to use.



75 Years

*Shaping the Future
of Endoscopy with you*

STORZ
KARL STORZ — ENDOSKOPE

THE DIAMOND STANDARD

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Telefon: +49 7461 708-0
Telefax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com