

ENDOFLATOR® 50 AND S-PILOT® STABLE CAVITY, SMOKE-FREE IMAGE



It's Time We Cleared the Air...

The ENDOFLATOR® 50 insufflator and S-PILOT® regulation for smoke-evacuation form an essential piece of the imaging chain, in order to provide excellent visibility during endoscopic interventions. This combination of units was created to have a stable cavity while removing dangerous toxins and unpleasant odors when needed.

STORZ
KARL STORZ – ENDOSKOPE

ENDOFLATOR® 50 AND S-PILOT®

STABLE CAVITY, SMOKE-FREE IMAGE

It's Time We Cleared the Air...

- With a maximum flow of 50 l/min, the ENDOFLATOR® 50 rapidly creates a stable working cavity
- The heating element warms the CO₂ gas to a body temperature of 37° C
- The SECUVENT® safety system ensures that overpressure is immediately detected and excess gas is released
- Maximum flow and pressure can be preset to suit the needs of specific patients and procedures
- The S-PILOT® has a simple pinch valve and footswitch design that allows smooth integration with existing systems or use as a stand-alone unit

Ordering Information	
UI500	ENDOFLATOR® 50, with integrated SCB module
UP501S1	S-PILOT®, including footswitch
031210-10	Heated Insufflation Tubing Set, with gas filter, sterile, for single use, package of 10
031200-10	Insufflation Tubing Set, with gas filter, sterile, for single use, package of 10
031447-10	Tubing Set for Smoke, Gas and Fluid Suction, sterile, for single use, package of 10
031111-10	Smoke-Evacuation Filter, unsterile, for single use, package of 10
UP004	S-PILOT® Connecting Cable, diameter 3.5 mm, length 300 cm
UP005	S-PILOT® Connecting Cable, diameter 2.5 mm, length 300 cm



ENDOFLATOR® 50



S-PILOT®

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com