

THE ONLY RETRACTABLE DURA KNIFE FOR ENDOSCOPIC ENDONASAL SURGERY



Enhance Your Surgical Workflow

The KARL STORZ Dura Knife is optimally designed for endoscopic endonasal skull base surgery. The knife includes a retractable blade that allows for safe insertion through the nasal pathways and precise dural incisions.



- Enable safe introduction of blade during endoscopic endonasal surgery
- Retractable blade minimizes the risk of trauma or damage to surrounding tissue
- Sterile, single-use blade provides sharp, precise cut every procedure
- Long, thin shaft to optimize utilization alongside an endoscope

STORZ
KARL STORZ – ENDOSKOPE

THE ONLY RETRACTABLE DURA KNIFE FOR ENDOSCOPIC ENDONASAL SURGERY

Enhance Your Surgical Workflow

The KARL STORZ Dura Knife includes a reusable handle and sheath that accommodates two styles of single-use blade inserts. See below for product information.

Dura Knife Sets, with Retractable Blade Inserts		
	28164KK-KT	de DIVITIIS-CAPPABIANCA Scalpel , with sickle-shaped retractable blade
	28164M-KT	de DIVITIIS-CAPPABIANCA Scalpel , with pointed retractable blade

Dura Knife, Blade Inserts		
	28164MC-S	Micro Knife , pointed, sterile, for single use, package of 3
	28164MS-S	Micro Knife , sickle-shaped, sterile, for single use, package of 3

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2020 KARL STORZ Endoscopy-America, Inc.