

ENDOMAT® SELECT CLEARVISION® LENS CLEANING SYSTEM



The Choice is Clear

The ENDOMAT® SELECT is the latest innovation in lens cleaning that provides smooth, reliable, and consistent performance. Designed to maintain a clear view, the CLEARVISION® module must be used in combination with irrigation sheaths. Features include:







- Easy and Simple Setup
- User-friendly Touchscreen Interface
- Tubing Cassette Recognition
- Multi-specialty Fluid Management System in One Console

STORZ
KARL STORZ – ENDOSKOPE

ENDOMAT® SELECT CLEARVISION® LENS CLEANING SYSTEM

The Choice is Clear

The CLEARVISION® module flushes the distal lens with fluid before the continuous oscillation of the pump head draws the drops of liquid away from the lens.

Ordering Information	
UP210-HN	ENDOMAT® SELECT Pump
Software	
UP606	ENT/NEURO Software , license (allows selection of CLEARVISION® procedure)
UP610*	Advanced Software , license (extends functions of installed software packages)
Accessories	
	031529-10-HN Tubing Set, irrigation , CLEARVISION®, sterile, for single use, package of 10
	20014230-HN One-Pedal Footswitch , digital, two-stage
	7230AS Irrigation and Suction Sheath 0° , outer diameter 4.8 x 6 mm, working length 14 cm, for use with HOPKINS® Telescopes 7230AA, 28132AA, TIPCAM®1 S 3D 7240AA 3D and 28164AA 3D as well as ENDOMAT® SELECT UP210. Suction and irrigation sheath also available in 30°, 45° and 70°
	7241AS Irrigation Sheath 0° , outer diameter 4.8 mm, working length 14 cm, connection facing upwards, for use with TIPCAM®1 S 3D 7240AA 3D or HOPKINS® Telescope 7230AA and ENDOMAT® SELECT UP210. Suction and irrigation sheath also available in 30°, 45° and 70°
	7230AES Irrigation and Suction Sheath 15°-90° , outer diameter 4.8 mm, working length 14 cm ENDOCAMELEON® ENT HOPKINS® Telescope 7230AE, ENDOCAMELEON® NEURO HOPKINS® Telescope 28132AE and ENDOMAT® SELECT UP210
	28164CBA-ENT Irrigation Sheath 0° , outer diameter 5 mm, working length 14 cm, for use with HOPKINS® Telescopes 7230AA, 28132AA, TIPCAM®1 S 3D 7240AA3D and 28164AA 3D. Suction and irrigation sheath also available in 30°, 45° and 70°.

*Optional UP610 Advanced Software

- Adjust pump height according to patient position
- Modify brightness of touchscreen
- Customize key tone volume
- SCB integration

**Contact your local KARL STORZ representative for more information
about these or other KARL STORZ products, or call 800-421-0837.**

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com