



WUIS9019 COLLABORATOR PROCESSOR PRO

PROCESSOR UPGRADE PACKAGE (SOFTWARE V4.7.5 OR LATER)

HIGHLIGHTS

- Includes additional processing power to display up to (8) HD video sources simultaneously on the CollaboratOR screen via (2) optionally attached WUIS9021 Quad Hardline Input Modules
- Replaces the standard CollaboratOR Processor



HARDWARE

FEATURE	DETAIL	ADDITIONAL NOTES
Connections	(4) Thunderbolt 3 (2) USB 3.0 (1) 3.5 mm Audio Out (2) HDMI Out (8) Optional HDMI In	Max Resolution (3840 x 2160) Landscape or Portrait Modes Requires 2x WUIS9021 Quad Hardline Input Modules
Network	(2) 10G RJ45 Ethernet Optional: Wi-Fi 802.11ac or 802.11 a/b/g/n	
Power	Voltage: 100 VAC - 240 VAC Frequency: 50 Hz - 60 Hz Power Consumption: 1280 W (12A - 6A) BTU: 4367 BTU Max Operating Temperature: 50°F - 95°F (10°C - 35°C)	
Physical	Dimensions (W x H x D): 19" x 8.7" x 21.2" (482 x 221 x 539 mm) Weight: 38.8 lbs. (17.6 Kg)	5 RU Height


Please note that product details are subject to change without prior notice.
For the most up-to-date product information, refer to the KARL STORZ website (www.karlstorz.com).

WUIS9019 COLLABORATOR PROCESSOR PRO

PROCESSOR UPGRADE PACKAGE (SOFTWARE V4.7.5 OR LATER)

TYPICAL SYSTEM

- Diagram below shows the typical CollaboratOR hardware with the WUIS9019 Processor Pro upgrade package



HARDWARE

FEATURE	DETAIL	ADDITIONAL NOTES
Power	Voltage: 100 VAC - 240 VAC Frequency: 50 Hz - 60 Hz Combined Power Consumption: 1324 W (13A - 6A) Combined BTU: 4673 BTU Max Operating Temperature: 50°F - 95°F (10°C - 35°C)	
Physical	Combined Dimensions (W x H x D): 19" x 19.25" x 21.2" (482 x 489 x 539 mm) Combined Weight: 88.8 lbs. (40.3 Kg)	11 RU Height Combined

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com