


WUIS9020 AND WUIS9021

HARDLINE INPUT MODULES (COLLABORATOR SOFTWARE V4.7.5 OR LATER)

HIGHLIGHTS

- Ingest (2) or (4) 1080p HDMI or DVI video feeds for simultaneous display on the CollaboratOR
- Optional TouchControl of up to (2) or (4) devices from the CollaboratOR display via USB
- Gangable to allow for higher counts of hardline input video feeds (max 8 per CollaboratOR)



HARDWARE

FEATURE	DETAIL	ADDITIONAL NOTES
Connections	(2) Thunderbolt (1) Gigabit Ethernet (2) or (4) HDMI Inputs Optional: (1) to (4) USB-B	WUIS9020 Dual Hardline Input: (2) HDMI Inputs WUIS9021 Quad Hardline Input: (4) HDMI Inputs
Power	Voltage: 100 VAC - 240 VAC Frequency: 50 Hz - 60 Hz Power Consumption: WUIS9020 - 16 W WUIS9021 - 22 W BTU: WUIS9020 - 105 BTU Max WUIS9021 - 153 BTU Max	WUIS9020 Dual Hardline Input: (2) HDMI Inputs WUIS9021 Quad Hardline Input: (4) HDMI Inputs
Physical	Dimensions (W x H x D): 19" x 5.25" x 14" / 483 x 134 x 356 mm / 3 RU Weight: 25 lbs. / 11.3 Kg	

Please note that product details are subject to change without prior notice.
For the most up-to-date product information, refer to the KARL STORZ website (www.karlstorz.com).

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com

www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com