

# BOOST EFFICIENCY WITH THE KARL STORZ FLUID-MANAGEMENT PACK



## One Convenient Package, Ready for Quick Clinical Setup!

Office procedures are growing in popularity, thanks to the many advantages it offers for the patient, physician, and overall healthcare system. The KARL STORZ high-quality customized all-in-one fluid-management pack combines ease and efficiency, meeting all procedural requirements in one convenient package.

- Reduce preparation time
- Fast and hassle-free post-procedural clean-up
- Increase procedural volume

KARL STORZ: Your partner for **QUALITY, EFFICIENCY, and SAFETY.**

**STORZ**  
KARL STORZ – ENDOSKOPE

# BOOST EFFICIENCY WITH THE KARL STORZ FLUID-MANAGEMENT PACK

## KARL STORZ FLUID-MANAGEMENT PACK – A customized solution for increased efficiency.

Imagine the time and costs that can be saved by including all components in a single, convenient and cost-effective package. At the end of the procedure, there should be NO FLUID on the floor or table – this allows for rapid room turnover between procedures!



### UNDER-BUTTOCKS DRAPE WITH GRADED MARKINGS

- Adhesive tape on the back secures it to the bed
- Comfortable paper – not cold plastic – for the patient to lie on
- Hand pockets to facilitate drape placement
- Metal malleable wire keeps the fluid collection pouch open to collect all the fluid
- Graduated collection pouch helps estimate the amount of fluid collected



### INFLOW IRRIGATION SET WITH LUER LOCK

- Cysto tubing diameter
- Luer-lock connection keeps the tubing firmly attached to the hysteroscope /sheath during the procedure
- Firm rubber distal end prevents kinking of the tubing during scope manipulation



### OUTFLOW DRAINAGE COLLECTION BAG

- Luer-lock connection
- Outflow tubing into a 2000 ml collection bag enables both fluid collection and monitoring of inflow and outflow to reveal potential fluid deficits

### Part Numbers

260814FM-5	Fluid Management Pack - 5 Packs/BOX
------------	-------------------------------------

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 7461 708-0  
Fax: +49 7461 708-105  
E-Mail: info@karlstorz.com  
  
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100  
Phone toll free: 800 421-0837 (US only)  
Fax: +1 424 218-8525  
Fax toll free: 800 321-1304 (US only)  
E-Mail: communications@karlstorz.com