

KARL STORZ CAMERA HEAD CONTAINER SYSTEM



Protect Your Camera Head Investment




- Validated and approved for Steam, STERRAD® 100S®, NX® and 100NX®, as well as STERIS® AMSCO® V-PRO® sterilization systems*
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion





* STERRAD® 100S®, NX® AND 100NX® are registered trademarks of Advanced Sterilization Products.
STERIS®, AMSCO® and V-PRO® are registered trademarks of STERIS Corp.

KARL STORZ CAMERA HEAD CONTAINER SYSTEM

Protect Your Camera Head Investment

The KSZ-39301HCTS Sterilization Container System is designed to house your KARL STORZ camera heads** through transport, sterilization, and storage. It offers cost-effective, hassle-free sterilization for your valuable camera head investment. KARL STORZ sterilization systems are made from aircraft grade anodized aluminum to prevent pitting and corrosion. They also help eliminate the need for blue sterilization wrap, saving time and money.

Accessories for H ₂ O ₂ /STERRAD® Sterilization		
	KSZ-SCF02	Polypro Single Use Filter for vaporized H ₂ O ₂ /STERRAD® sterilization (1,000 per case), hydrophobic non-woven polypropylene
	KSZ-SCLH2023	Load Cards (1,000 per case) Large process indicator cards for vaporized H ₂ O ₂ /STERRAD® sterilization
	KSZ-SCS01W	Tamper Evident Seals for vaporized H ₂ O ₂ /STERRAD® sterilization (1,000 per pack) W=White for Hydrogen Peroxide H ₂ O ₂ /STERRAD® & V-PRO®

Accessories for Steam Sterilization (for Autoclavable Camera Heads Only)		
	KSZ-SCF01	Paper Single Use Filter (1,000 per case) 100% cellulose and water repellant, 7.5" round
	KSZ-SCL01	Dual Process Indicator Cards (1,000 per case) ID card with a dual-indicator dot
	KSZ-SCI001	Dual Process Indicator (1,000 per case) Dual purpose indicator strip with room for process identification, 8" long with perforation to reconfigure up to two 4" long indicator strips
	KSZ-SCS01B	Tamper Evident Seals (1,000 per pack) B=Blue for steam, KSZ-SCS01R (red pictured)

**Not recommended for sterilization processing of Pendulum camera heads.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postbox 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-Mail: info@karlstorz.com
www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100
Phone toll free: 800 421-0837 (US only)
Fax: +1 424 218-8525
Fax toll free: 800 321-1304 (US only)
E-Mail: communications@karlstorz.com

Caution: Federal (USA) law restricts this device to sale by or on the order of physician.

© 2020 KARL STORZ Endoscopy-America, Inc.