

KARL STORZ VIDEO MEDIASTINOSCOPY SET CONTAINER SYSTEM



Protect your Video Mediastinoscopy Instrument Set Investment

Our KSZ-39312VMS Container System is specially designed to protect your video mediastinoscopy instrument sets during transport, sterilization, and storage.





- Designed for video mediastinoscopy instruments
- Save time and money over blue wrap
- Anodized to prevent pitting and corrosion

STORZ
KARL STORZ – ENDOSKOPE

KARL STORZ VIDEO MEDIASTINOSCOPY SET CONTAINER SYSTEM

Cost-effective, hassle free sterilization

The KSZ-39312VMS Video Mediastinoscopy Set sterilization system is designed to house your KARL STORZ video mediastinoscopy instruments through transport, sterilization, and storage. It offers cost-effective sterilization for your valuable video mediastinoscopy instrument set investment. KARL STORZ sterilization systems are made from aircraft grade anodized aluminum to prevent pitting and corrosion. They eliminate the need for blue sterilization wrap, saving time and money.

Accessories for steam sterilization		Description
	KSZ-SCF01	Disposable paper filters for steam sterilization, 1,000 per case. 100% Cellulose and water repellent, 7.5" round.
	KSZ-SCL01	Steam Sterilization Dual Process Indicator ID Cards, 1,000 per case.
	KSZ-SCI001	Steam Sterilization Dual Process Indicator Strips, 1,000 per case. A dual purpose indicator strip with room for process identification. 8" long with perforation to reconfigure up to two 4" long indicator strips.
	KSZ-SCS01B	Steam Sterilization Tamper Evident Seals 1,000 per pack.

Contact your local KARL STORZ representative for more information about these or other KARL STORZ products, or call 800-421-0837.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ SE & Co. KG
 Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
 Postbox 230, 78503 Tuttlingen/Germany
 Phone: +49 7461 708-0
 Fax: +49 7461 708-105
 E-Mail: info@karlstorz.com
 www.karlstorz.com



KARL STORZ Endoscopy-America, Inc.
 2151 East Grand Avenue
 El Segundo, CA 90245-5017, USA
 Phone: +1 424 218-8100
 Phone toll free: 800 421-0837 (US only)
 Fax: +1 424 218-8525
 Fax toll free: 800 321-1304 (US only)
 E-Mail: communications@karlstorz.com